



Report No.: 10489 Date: 03.06.2003

Contact: HP

Task:

Lab Unlir

Application field: Mineralogy / Metallurgy

Material: Lignite and brown coal

Feed size: 0-50 mm

Feed quantity: 1250 g (per batch)

Material medium-hard, moist

specification(s):

Customer < 3 mm

requirement(s):

Subsequent Determination of humidity and calorific value

Solution:

analysis:

Selected SM 100 Cutting Mill

instrument(s):

Configuration(s): Bottom sieves of stainless steel, square holes 6 mm;

Ring-type filter with trapezoid hole body for collecting

receptacle 5 l; Standard hopper

Parameter(s): Revolution speed 1450 rpm

Time: 5 min. (per batch)

Achieved predominantly < 3 mm

result(s):

Remark(s):

Recommendation: For sample preparation of lignite and brown coal, the

Cutting Mill SM 100 is suitable under the above mentioned

conditions.



Report No.: 10489 Date: 03.06.2003

Contact: HP

Pictures of the sample





Fig. 1: Lignite with 60 % rest humidity Fig. 2: Brown coal, 45 % humidity

Tallaght Business Park Whitestown, Dublin 24, Ireland D24 RFK3

Tel: (01) 4523432 Fax: (01) 4523967

Web: www.labunlimited.com Web: www.labunlimited.co.uk

Quatro House, Frimley Road, Camberley, United Kingdom **GU16 7ER**

Tel: 08452 30 40 30 Fax: 08452 30 50 30

E-mail: info@labunlimited.com E-mail: info@labunlimited.co.uk