





**Report No.:** 11409 Date: 01.06.2006

Contact: DE

Task:

**Application field:** Building Materials

Material: 1. Cement clinker; 2. Slag; 3. Mud sediment

Feed size: 3-7 mm

**Feed quantity**: 30 g (+ 3 tablets "Spectromelt", for each sample)

Material hard brittle, dry

specification(s):

**Customer** < 100 μm; for X-ray fluorescence analysis

requirement(s):

**Subsequent** X-ray Fluorescence Analysis

Solution:

analysis:

Selected Planetary Ball Mill PM 100

instrument(s):

Configuration(s): Grinding jar "comfort" tungsten carbide 125 ml;

7 x grinding ball tungsten carbide ø 20 mm

Parameter(s): Revolution speed 450 rpm

Time: 5 min. (sample 1 and 2); 2 min. (sample 3)

Achieved 1. Cement clinker:  $100 \% < 100 \mu m$  result(s): 2. Slag:  $99 \% < 100 \mu m$ 

3. Mud-sediment:  $93 \% < 100 \mu m$ 

**Remark(s):** To avoid the agglomeration of the powders, 3 tablets

"Spectromelt" were added before grinding.

Due to its two large initial particle size the mud sediment

(sample 3) had to be pre-crushed manually with a

hammer.

**Recommendation**: For the grinding of hard brittle clinker and slag samples the

Planetary Ball Mill PM 100 is suitable according to the above

mentioned conditions.



Report No.: 11409 Date: 01.06.2006

Contact: DE

## Pictures of the sample




Fig. 1: Cement clinker sample before and after grinding in PM 100 for 5 min.



Fig. 2: Slag sample before and after grinding in PM 100 for 5 min.



Fig. 3: Mud sediment sample before and after manual pre-crushing with a hammer



Fig. 4: Mud sediment sample after grinding in PM 100 for 5 min.



Tallaght Business Park Whitestown, Dublin 24, Ireland D24 RFK3

Tel: (01) 4523432 Fax: (01) 4523967

Web: www.labunlimited.com Web: www.labunlimited.co.uk

Quatro House, Frimley Road, Camberley, United Kingdom **GU16 7ER** 

Tel: 08452 30 40 30 Fax: 08452 30 50 30

E-mail: info@labunlimited.com E-mail: info@labunlimited.co.uk





**Report No.: 17047** Date: 14.03.2016

Contact: MLU

Task:

Application field: Environment / Recycling

Material: Soil

Feed size: 2 mm

Feed quantity: 30 g

Material

specification(s):

Customer

212 µm

dry

requirements(s):

. ( )

Subsequent

analysis:

XRF - X-ray Fluorescence Analysis ICP-MS - Mass Spectroscopy

Heavy metals

Solution

Selected Planetary Ball Mill PM 100

Instrument(s):
Configuration(s)

Configuration(s) 1 x PM 100, 230 V, 50/60 Hz

Item nos.:

1 x Grinding jar "comfort", tungsten carbide, 80 ml

5 x Grinding ball, tungsten carbide, 20 mm ø

Please note: Other electrical versions of the instrument(s) are

available with different item numbers.

Parameter(s): Revolution speed: 200 rpm

Time: 1 min

Achieved result(s): d90: < 212 µm

**Remark(s):** The sample was pre-crushed with a Jaw Crusher BB 200.

Recommendation: The Planetary Ball Mill PM 100 is suitable to grind the sample material

under the above mentioned conditions.





**Report No.:** 17047 Date: 14.03.2016

Contact: MLU

## Pictures of the sample:



Picture 1: Pre-crushed sample



Picture 2: Sample after grinding process with **PM 100** 



Tallaght Business Park Whitestown, Dublin 24, Ireland D24 RFK3


Camberley, United Kingdom GU16 7ER

Tel: (01) 4523432 Fax: (01) 4523967 Tel: 08452 30 40 30 Fax: 08452 30 50 30

Web: www.labunlimited.com Web: www.labunlimited.co.uk

E-mail: info@labunlimited.com E-mail: info@labunlimited.co.uk

Quatro House, Frimley Road,







**Report No.:** 10411
Date: 07.03.2003

Contact: HP

Task:

**Application field:** Agriculture

Material: Soil samples, with high content of clay and minerals

Feed size: 0-10 mm (different to the samples)

Feed quantity: 120 q

**Material** dry

specification(s):

Customer  $100 - 150 \text{ ml appx.} < 50 \text{ }\mu\text{m}$ 

requirement(s):

Subsequent analysis:

**ient** ICP Inductive Coupling Plasma

**Solution:** 

Selected PM 100 CM Planetary Ball Mill

instrument(s):

**Configuration(s):** Grinding jar 250 ml, Sintered alumina (Al<sub>2</sub>O<sub>3</sub>);

5 grinding balls, Sintered alumina, 30 mm Ø

**Parameter(s):** Revolution speed 450 rpm

**Time:** 10 min. (per sample)

**Achieved** 98 % < 50  $\mu$ m result(s): see samples

**Remark(s):** Moistly samples must be dried before grinding.

To avoid sticking effects of the material during grinding, it is necessary to add 5 - 8 drops of methanol into the jar.

**Recommendation:** For fine grinding of dry soil samples, the Planetary Ball Mill

PM 100 CM is suitable under the above mentioned

conditions



Report No.: 10411 Date: 07.03.2003 Contact: ΗP

## Pictures of the sample



Fig. 1: Soil with clay



Fig. 2: Soil with minerals

Tallaght Business Park Whitestown, Dublin 24, Ireland D24 RFK3

Tel: (01) 4523432 Fax: (01) 4523967 Quatro House, Frimley Road, Camberley, United Kingdom **GU16 7ER** 

Tel: 08452 30 40 30 Fax: 08452 30 50 30 E-mail: info@labunlimited.com E-mail: info@labunlimited.co.uk Web: www.labunlimited.com Web: www.labunlimited.co.uk