# Data sheet

# 2X HotBegan™ Red-Taq Master Mix

Cat. No: P0320 Cat. No: P0321

#### Introduction

2X HotBegan<sup>™</sup> Red-Taq Master Mix is an optimized ready-to-use solution containing HotBegan Taq DNA Polymerase (hot start performance), dNTPs, MgCl<sub>2</sub> and stabilizers. It is inactive at room temperature and only requires addition of template, primers, and water.

HotBegan Taq DNA polymerase is a Taq DNA polymerase bound to a proprietary antibody that blocks polymerase activity until denaturation step occurs. The heat labile antibodies are rapidly inactivated by raising the temperature (4 minutes at 95-97°C). This prevents or minimizes primer-dimer and nonspecific products.

The mix also contains an agarose loading buffer including a red dye for visual tracking of DNA migration and a dense compound to facilitate the drop-down of the samples into the well agarose gels. Like the Tag polymerase, the enzyme has  $5' \rightarrow 3'$ polymerase activity and a weak  $5' \rightarrow 3'$  exonuclease activity but no 3′→5′ exonuclease (proofreading). Before enzyme activation none of enzyme activities are detectable.

#### **Features**

- Inactive at room temperature.
- Adds extra nucleotides (preferentially adenine) without template at 3'ends leaving 3'overhangs PCR fragments. This fact allows the popular TAcloning or GC cloning.
- Both save times in the PCR process and in agarose loading samples.
- Amplifies from a femtograms of DNA targets.

#### **Applications**

- PCR fragments amplification for TA or GC cloning
- Design for high throughput applications.
- Amplification from a limited DNA template or low copy number genes.

| Components                         | P0320      | P0321        |
|------------------------------------|------------|--------------|
| 2X HotBegan Red-Taq<br>Master Mix* | 5x 1.25 mL | 10 x 1.25 mL |
| 50mM MgCl <sub>2</sub> Solution**  | 1.5 mL     | 1.5 mL       |

- \*2X HotBegan Red-Taq Master Mix include HotBegan Taq DNA polymerase, 2X Red buffer, 0.4 mM of each dNTP, 5 mM Mg2<sup>+</sup> and 5% Glycerol.
- \*\*Separate tube 50 mM MgCl<sub>2</sub> solution is provided for further optimisation. In some cases, we recommend to optimize Mg2<sup>+</sup> concentration.

### Assay conditions

25mM Tris-HCl pH9.0 at 25°C, 50mM KCl, 2mM MgCl<sub>2</sub>, 0.1mg/mL gelatine, 200 µM of dATP, dGTP, dTTP,  $100\mu M[\alpha 32-P]dCTP$  (0.05 $\mu Ci/nmol$ ) and 12.5  $\mu g$  activated salmon sperm DNA.

Unit definition: One unit is defined as the amount of enzyme required to catalyse the incorporation of 10 nanomoles of dNTPs into acid-insoluble material in 30 minutes at 74°C.

# **Quality Certifications**

- Functionally tested in PCR.
- Not detectable activity of nucleases (endo-, exo, and ribo-).

**Storage:** Upon receipt, store the entire kit at -20 °C

(Continued on reverse side)



Tallaght Business Park Whitestown, Dublin 24, D24 RFK3

Tel: (01) 4523432 Fax: (01) 4523967 Quatro House, Frimley Road, Camberley, United Kingdom GU16 7ER

Tel: 08452 30 40 30 Fax: 08452 30 50 30 E-mail: info@labunlimited.com E-mail: info@labunlimited.co.uk Web: www.labunlimited.com Web: www.labunlimited.co.uk





# Recomended PCR assay

|                                 | 20 µl assay                          | 50 µl assay                          |
|---------------------------------|--------------------------------------|--------------------------------------|
| 2X HotBegan Red- Taq Master Mix | 10 μl (1X)                           | 25 μl (1X)                           |
| Forward Primer (15μM)           | 0.75 pmol/μL                         | 0.75 pmol/μL                         |
| Reverse Primer (15μM)           | 0.75 pmol/μL                         | 0.75 pmol/μL                         |
| Template DNA                    | plasmid: 30-75ng;<br>gDNA: 100-500ng | plasmid: 30-75ng;<br>gDNA: 100-500ng |
| Nuclease-free water             | up to 20 μL                          | up to 50 μL                          |

## **Cycling instructions:**

1x 94°C 10:00; 30x (94°C 0:35, Tm 0:35, 72°C 1'/kb); 1x 72°C 7:00; 1x 4°C  $\infty$ 

This procedure is intended for use as a guide only and may need optimization. The optimal conditions will vary from reaction to reaction and are dependent on the template/primers used.

| Red dye Agarose Mobility*        |                               |                     |  |
|----------------------------------|-------------------------------|---------------------|--|
| Agarose Gel<br>Concentration (%) | Effective separation of: (bp) | Migration Rate (bp) |  |
| 0,7                              | 800-12000                     | 3000                |  |
| 1,0                              | 400-8000                      | 1500                |  |
| 1,5                              | 200-3000                      | 900                 |  |
| 2,0                              | 100-2000                      | 300                 |  |
| 3,0                              | 25-1000                       | > 100               |  |

<sup>\*</sup> in TAE Buffer

#### PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively for research purposes and in vitro use only. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals. Please refer to www.canvaxbiotech.com for Material Safety Data Sheet of the product.



Tallaght Business Park Whitestown, Dublin 24, Ireland D24 RFK3

Tel: (01) 4523432

Web: www.labunlimited.com

Quatro House, Frimley Road, Camberley, United Kingdom GU16 7ER

Tel: 08452 30 40 30 Fax: 08452 30 50 30 E-mail: info@labunlimited.com E-mail: info@labunlimited.co.uk Web: www.labunlimited.co.uk

