GE Healthcare

Whatman

Filtration Devices

Mini-UniPrep Syringeless Filter

Prepare HPLC Samples Three Times Faster.

Whatman Mini-UniPrep™ Syringeless Filters provide a faster, easier way to remove particulates from samples being prepared for High Performance Liquid Chromatography (HPLC) analysis. Add up the time savings, plus the money saved from cutting multiple consumables out of the sample preparation process, and you'll see huge benefits for your lab.

Features and Benefits

- **Ease of use** Filtering and storage vial allows you to process sample loads in one-third the time.
- Wide range of membranes Suits a variety of sample types.
- Available in 0.2 and 0.45 µm pore sizes Meets specific sample requirements.
- Compatible with all major autosamplers Allows you to speed processes.
- Fewer consumables required Lowers costs by up to 40 percent.

Applications

- Routine analysis
- Composite assays
- Content uniformity
- Protein precipitation
- Solubility testing
- Dissolution testing
- Sample filtration





The Mini-UniPrep Syringeless Filter on the left is shown with fluid in the chamber. On the right, the filter plunger is shown compressed with the sample ready for analysis.

Protein Precipitation 'Spin' Method vs. Mini-UniPrep Method



A Variety Of Special Mini-UniPrep Filters To Meet Your Specific Needs

In a process of continuous improvement and innovation, Whatman has listened to customers and created a whole family of Mini-UniPrep filters to meet specific needs. For customers using robotics to maximize throughput, Whatman offers Slit Septa Mini-UniPrep. For customers who need to protect samples from light damage, there's Amber Mini-UniPrep.

Slit Septa Mini-UniPrep Syringeless Filter

For High Throughput Automation.



Features and Benefits

- Slit septa cap Enables Mini-UniPrep use with curent robotics on HPLC instruments for high throughput automation.
- Increases needle longevity Eliminates coring problems associated with repeated sampling.
- Available in 0.2 and 0.45 µm pore sizes Meets specific sample requirements.
- **Durable yet flexible** Slit septa cap has been specially designed for instruments with sensitive sampling needs.

Applications

• Use with standard robotics on HPLC instruments with sensitive needles, allowing for higher throughput.

Amber Mini-UniPrep Syringeless Filter

Protects Samples From Light Damage.

Features and Benefits

- Amber colorant Prevents photo degradation of light-sensitive samples. Same colorant used in pharmaceutical containers designed to meet USP specifications for light resistance.
- Translucent amber chamber and plunger Easy visual inspection.

Applications

• Use with any compound that requires protection from light, such as catecholamine or vitamins.

Six-Position Compressor Makes Your Job Even Easier



The optional compressor allows you to process up to six Mini-UniPrep filters at the same time, further speeding your workflow and reducing the risk of hand stress.



Mini-UniPrep Makes It Faster And Easier To Process HPLC Samples.



Simple.

Place unfiltered sample in chamber.

Innovative.



Compress filter plunger into sample chamber. Clean filtrate fills reservoir from bottom up.

Convenient.

The Mini-UniPrep vial shape fits easily in autosamplers.

The filter membrane is housed at one end of a plunger, with a cap/ septum on the other. By pressing the plunger through the sample placed into the chamber, positive pressure forces the filtrate into the reservoir of the plunger. Air escapes through vent holes until the evaporation seal is engaged, providing an airtight seal. Within seconds the Mini-UniPrep can be placed into your auto sampler for injection into your instrument.

Choose the Right Mini-UniPrep Filtering Media

Sample Type	Suitable Mini-UniPrep Media	
Particulate laden liquids	Glass Microfiber (GMF)	
Aqueous/Organic samples in 3 to 10 pH range	Nylon (NYL)	
General filtration media/Solvent based samples	Polypropylene (PP)	
Chemically aggressive solutions	Polytetrafluroethylene (PTFE)	
Biological samples requiring low protein binding media	Polyethersulfone (PES) or Regenerated Cellulose (RC)	
Aqueous/Organic solvents – low non-specific protein binding media	Polyvinylidenefluoride (PVDF)	
Aqueous/Organic solvents - high flow and loading capacity	Polypropylene Depth (dp PP)	

Technical Specifications

Dimensions	Equivalent in size to 12 mm x 32 mm vials
Materials of Construction	Housing and Cap: Polypropylene Filter Media: As specified Septa: PTFE Membrane and Silicone Rubber
Filtering Capacity	0.4 mL
Nominal Force Needed to Compress	Approximately 8 psi (0.6 Bar)
Max. Operating Temperature	120 °F (50 °C)

Ordering Information - Mini-UniPrep Syringeless Filters

Standard Cap - Translucent Housing					
Catalog Number	Pore Size (µm)	Pack Size			
Nylon Filtration Media					
UN203NPENYL	0.2	100			
UN503NPENYL	0.2	1000			
UN203NPUNYL	0.45	100			
UN503NPUNYL	0.45	1000			
PES Filtration Media					
UN203NPEPES	0.2	100			
UN503NPEPES	0.2	1000			
UN203NPUPES	0.45	100			
UN503NPUPES	0.45	1000			
PVDF Filtration Media					
UN203NPEAQU	0.2	100			
UN503NPEAQU	0.2	1000			
UN203NPUAQU	0.45	100			
UN503NPUAQU	0.45	1000			
RC Filtration Media					
UN203NPERC	0.2	100			
UN503NPERC	0.2	1000			
UN203NPURC	0.45	100			
UN503NPURC	0.45	1000			
PTFE Filtration Media					
UN203NPEORG	0.2	100			
UN503NPEORG	0.2	1000			
UN203NPUORG	0.45	100			
UN503NPUORG	0.45	1000			
PP Filtration Media					
UN203NPEPP	0.2	100			
UN503NPEPP	0.2	1000			
UN203NPUPP	0.45	100			
UN503NPUPP	0.45	1000			
DEPTH PP Filtration Media					
UN203NPUDPP	0.45	100			
UN503NPUDPP	0.45	1000			
GMF Filtration Media					
UN203NPUGMF	0.45	100			
UN503NPUGMF	0.45	1000			

Ordering Information - Mini-UniPrep Syringeless Filters (Continued)

Slit Septum Cap - Translucent Housing						
Catalog Number	Pore Size (µm)	Pack Size				
Nylon Filtration Media						
US203NPENYL	0.2	100				
US503NPENYL	0.2	1000				
US203NPUNYL	0.45	100				
US503NPUNYL	0.45	1000				
PES Filtration Media						
US203NPEPES	0.2	100				
US503NPEPES	0.2	1000				
US203NPUPES	0.45	100				
US503NPUPES	0.45	1000				
PVDF Filtration Media						
US203NPEAQU	0.2	100				
US503NPEAQU	0.2	1000				
US203NPUAQU	0.45	100				
US503NPUAQU	0.45	1000				
PTFE Filtration Media						
US203NPEORG	0.2	100				
US503NPEORG	0.2	1000				
US203NPUORG	0.45	100				
US503NPUORG	0.45	1000				
PP Filtration Media						
US203NPEPP	0.2	100				
US503NPEPP	0.2	1000				
US203NPUPP	0.45	100				
US503NPUPP	0.45	1000				
DEPTH PP Filtration Media						
US203NPUDPP	0.45	100				
US503NPUDPP	0.45	1000				
GMF Filtration Media						
US203NPUGMF	0.45	100				
US503NPUGMF	0.45	1000				

Amber Housing (for Light Sensitive Samples) - Standard Cap								
Catalog Number		Pore Size (µm)		Pack Size				
Nylon Filtration Media								
UN203APENYL		0.2		100				
UN203APUNYL		0.45		100				
PES Filtration Media								
UN203APEPES		0.2		100				
UN203APUPES		0.45		100				
PVDF Filtration Media								
UN203APEAQU		0.2		100				
UN203APUAQU		0.45		100				
PTFE Filtration Media								
UN203APEORG		0.2			100			
UN203APUORG		0.45			100			
PP Filtration Media								
UN203APUDPP		0.2		100				
UN203APUPP	UN203APUPP			100				
DEPTH PP Filtration Media								
UN203APUDPP		0.45		100				
GMF Filtration Media								
UN203APUGMF		0.45	0.45 100		100			
Amber Housing (for Light Sensitive Samples) - Silt Septum Cap								
Catalog Number		Description	Pore	Size (µm)	Pack Size			
Nylon Filtration Membrane								
UN203APUNYL	Mini Un	UniPrep Syringeless Filter		0.45	100			
Accessories								
Catalog Number		Description	Pa		nck Size			
Six Position Compressor								
CR0000006		Mini UniPrep Compressor		1				
Adapter								
UN2031545PP		Adapter Sleeve for A Sampler Trays Made 15 x 45 mm Vials	to Fit	to 150 p Fit				

Distributed by:



ADVANCED APPLIED TECHNOLOGIES

Contact Us: Irl Ph: 01 4523432 UK Ph: 08452 30 40 30 Web: www.carlstuart.com Email: info@carlstuart.com



imagination at work

GE, imagination at work and GE monogram are trademarks of General Electric Company.

Whatman and Mini-UniPrep are trademarks of GE Healthcare companies.

All third party trademarks are the property of their respective owners.

© 2008 General Electric Company – All rights reserved.

First published August 2008.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

Whatman Inc., a General Electric Company, going to market as GE Healthcare.