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**Pressureless plastic mixed bed ion exchange cartridge consisting of 01.1753 and 04.1609**

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**Pressureless plastic mixed bed ion exchange cartridge Art.-Nr.: 01.1705 consisting of 01.1750 and 04.1612**

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**Read these operating instructions carefully before installing and starting-up the system!**

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## ECC-Declaration of conformity

**According to EG Directive 89/336/EEC**  
**- Electromagnetic Compatibility -**

We hereby declare that the machine specified below conforms, in its design and construction, and in the versions we have introduced into the market, to the fundamental safety and health requirements of the EEC Directive 89/336/EEC.

This declaration loses its validity when changes which were not agreed to by us are made to the machine.

Designation of the machine: **Ion Exchange System**  
Machine type: **TKA DI 750**  
Article number: **01.1705 / 04.1609**  
Applicable EEC Directives: **EEC Directive Low voltage (73/23/EEC)**  
**EEC Directive Electromagnetic Compatibility (89/336/EEC)**  
Applied Standards: **EN 55014**

Niederelbert, 10. March 2003

  
Manufacturer / authorized representative

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Unit for direct connection to a water supply, with free flow through and out of the open outlet. Supplies deionized water for laboratories, doctors and dentists practices.

The benefits:

- Supplies purified water at a very favourable cost
■ Easy to use
■ Extremely effective
■ Lightproof blue cartridge avoids growth of algae
■ Wall mounting (option) allows the purified water supply to be most practically positioned at the workplace

Table with 2 columns: Specification and Value. Title: Maße und Leistungen TKA DI 750 (Art.-Nr. 01.1753). Includes rows for flow rate, capacity, water quality, temperature, power supply, dimensions, weight, and article number. Also includes a section for Leitfähigkeitsmessgerät (04.1609) and delivery scope (Seite 4).

Table with 2 columns: Specification and Value. Title: Performance + Specifications TKA DI 750 (01.1753). Includes rows for flow rate, capacity, water quality, temperature, power supply, dimensions, weight, and article number. Also includes a section for Conductivity meter (04.1609) and equipment supplied (page 4).

Table with 2 columns: Specification and Value. Title: Maße und Leistungen TKA DI 750 (Art.-Nr. 01.1705). Includes rows for flow rate, capacity, water quality, temperature, power supply, dimensions, weight, and article number. Also includes a section for Ersatzpatrone (01.1750) and delivery scope (Seite 6).

Table with 2 columns: Specification and Value. Title: Performance + Specifications TKA DI 750 (01.1705). Includes rows for flow rate, capacity, water quality, temperature, power supply, dimensions, weight, and article number. Also includes a section for Replacement cartridge (01.1750) and equipment supplied (page 6).

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
**Pressureless plastic  
mixed bed ion exchange  
cartridge,  
consisting of  
01.1753 and 04.1609**

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### Mounting Instructions

01. Screw the conductivity meter onto the cartridge.
02. Use the screws and plugs provided to mount the bracket near the water tap.
03. Connect the feedwater hose to the water tap and to the unit. Take care that the hose is not buckled.
04. Screw the outlet hose onto the 3/4" connection positioned on the neck of the conductivity meter. Connect the unit to 230V outlet.
05. The unit is ready for use: Simply turn on the water tap for deionized water to flow out of the outlet.

**CAUTION: Do not block the free flow of deionized water out of the outlet hose in any way, as the plastic cartridge could burst under build-up pressure!**

06. The unit can be used until the pointer of the conductivity meter leaves the green section or indicates 20 µS/cm.


### Change of Cartridge

The cartridge is exhausted, and must be replaced, when the pointer of the conductivity meter leaves the green section or indicates 20 µS/cm.

01. Close the water tap. Unplug the unit from the electrical supply. Disconnect the inlet hose from the cartridge. Allow the water to run out of the hose until it is empty.
02. Remove the cartridge from the bracket. Let the water run out of the unit.
03. Unscrew the conductivity meter from the exhausted cartridge and screw the 3/4" cap on. Send the exhausted cartridge in for regeneration.
04. Screw the conductivity meter onto the new cartridge. Attach the feedwater hose to the cartridge.
05. Fit the cartridge, with the feedwater hose attached, in the bracket.
06. Plug it into the 230V supply. The unit is again ready for use.

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**Pressureless plastic  
mixed bed ion exchange  
cartridge  
Art.-Nr.: 01.1705  
consisting of  
01.1750 and 04.1612**

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## Safety Precautions:

01. Using the cartridge in any other way than that described in these operating instructions invalidates the guarantee.
02. **Deionized water must on no account be consumed in any way. It is not fit for drinking!**

## Mounting Instructions

01. Use the screws and plugs provided to mount the bracket near the water tap.
02. Connect the feedwater hose to the water tap and to the unit. Take care that the hose is not buckled.
03. Ease the outlet hose onto the outlet hose nipple on the neck of the conductivity meter (warm the hose if necessary). Connect the unit to 230V outlet.
04. The unit is ready for use: Simply turn on the water tap for deionized water to flow out of the outlet.

**CAUTION: Do not block the free flow of deionized water out of the outlet hose in any way, as the plastic cartridge could burst under build-up pressure!**

05. The unit can be used until the pointer of the conductivity meter leaves the green section or indicates 20  $\mu\text{S}/\text{cm}$ .

## Change of Cartridge

The cartridge is exhausted, and must be replaced, when the pointer of the conductivity meter leaves the green section or indicates 20  $\mu\text{S}/\text{cm}$ .

01. Close the water tap. Unplug the unit from the electrical supply. Disconnect the inlet hose from the cartridge. Allow the water to run out of the hose until it is empty.
02. Remove the cartridge from the bracket. Let the water run out of the unit.
03. Unscrew conductivity meter and draw the tube out of the resin.
04. Before now connecting the new cartridge, turn it upside down to loosen the resin in it. Press the tube down into the resin and screw the conductivity meter on. Attach the feedwater hose to the new cartridge.
05. Fit the cartridge, with the feedwater hose attached, in the bracket.
06. Plug it into the 230V supply. The unit is again ready for use.

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**Carl Stuart Limited**
**TKA Wasseraufbereitungssysteme GmbH**