

# Data sheet

## pASSEMBLE™ Lentiviral Packaging System

Cat. No: ME0044

Cat. No: ME0044-Plus

50 transfection reactions in wells of 6-well plates.

### Description

To produce virions, once cloned your gene of interest into a viral vector (known as transfer vector), you need to use of packaging viral vector(s). Transfer vectors have been modified to reduce the hazardous level and highest generation of transfer vectors are lacking of *gag*, *pol*, *env* viral genes. Packaging vectors provide all the viral proteins required for transcription and packaging of your expression cassette into recombinant viral particles.

ASSEMBLE Lentiviral Packaging System includes an optimized mix of three vectors with sequences of *gag*, *pol* and *rev* genes from human immunodeficiency virus (HIV-1) and the envelope gene from vesicular stomatitis virus (VSV-G). VSV-G envelope confers a wide range of tropism because of this glycoprotein binds to phospholipid receptor universally expressed in mammalian cells. This packaging system requires lentiviral transfer vectors of 3rd generation or higher.

### Biosafety Features

- ✓ Biosafety Level 2 (BL-2).
- ✓ Mix of three packaging vectors. One of them contains *gag/pol* genes; other, *env* gene and the third, *rev* gene.
- ✓ Rev regulation sequence is located on *gag/pol* expression vector adding an additional biosafety element: in the absence of *rev*, *gag/pol* are not produced.
- ✓ Viruses from SIN vectors result in the transcriptional inactivation of the provirus in the infected cell.
- ✓ Disposable Packaging vectors are 3'LTR SIN and packaging signal deleted ( $\Delta\Psi$ -GAG).

### Kit Components

Components	ME0044	ME0044-Plus
pASSEMBLE™ Lentiviral Packaging Mix (0.5 µg/µL)	150 µL	150 µL
CANFAST Transfection Reagent	-	0.5 mL
eGFP Lentiviral Transfer Control Vector (0.5 µg/µL)	-	30 µL

The Kit is enough for 50 transfection reactions in wells of 6-well plates.

### Storage

If stored properly, shelf life is 1 year from the date of shipment.

The pASSEMBLE™ Lentiviral Packaging Mix and eGFP Lentiviral Transfer Control Vector should be stored at -20 °C. Avoid multiple freeze thaw cycles.

**Important:** Store CANFAST Transfection Reagent at 4°C. **DO NOT FREEZE.**

(Continued on reverse side)

Distributed by:

**Lab Unlimited**  
CARL STUART GROUP

Tallaght Business Park  
Whitestown, Dublin 24,  
Ireland  
D24 RFK3

Tel: (01) 4523432  
Fax: (01) 4523967  
E-mail: [info@labunlimited.com](mailto:info@labunlimited.com)  
Web: [www.labunlimited.com](http://www.labunlimited.com)

Quatro House, Frimley Road,  
Camberley,  
United Kingdom  
GU16 7ER

Tel: 08452 30 40 30  
Fax: 08452 30 50 30  
E-mail: [info@labunlimited.co.uk](mailto:info@labunlimited.co.uk)  
Web: [www.labunlimited.co.uk](http://www.labunlimited.co.uk)

 **canvax**

## Protocol

### Culturing Packaging Cell Lines

We recommended use HEK-293 cell line to produce viral supernatant.

### Virus Production (Protocol for 6-well plate)

This protocol is for 6-well plates, this protocol can be scaled up or down if desired.

1. For transfection of lentiviral vector, seed your packaging cell line at 60-80% confluence the day before transfection. The number of cells recommended for 6-well plates are  $2.5-8.0 \times 10^5$  cells/well.
2. On the day of transfection, prepare Transfection Mix as follows:

For each well of 6-well plates prepare:

- A. **DNA solution:** Dilute 0.5  $\mu\text{g}$  DNA plasmids and 1.5  $\mu\text{g}$  Lentivirus Packaging Mix in 100  $\mu\text{L}$  Medium without serum and mix gently.
- B. **CANFAST solution:** Dilute 6  $\mu\text{L}$  of CANFAST Transfection Reagent in 100  $\mu\text{L}$  Medium without serum and mix gently.
- C. Prepare the **Transfection Mix** adding **CANFAST solution** drop to drop into **DNA solution** which is gently stirring at vortex.

Tissue Culture Vessel	DNA Solution		CANFAST solution		Transfection Mix
	DNA	Medium without serum	CANFAST Reagent	Medium without serum	
1 well of a 6-well plate	2 $\mu\text{g}$	100 $\mu\text{L}$	6 $\mu\text{L}$	100 $\mu\text{L}$	200 $\mu\text{L}$

3. Incubate **Transfection Mix** 15-20 minutes at room temperature. Add **Transfection Mix** into each well by leaking and shake **gentle** the plate.
4. Incubate for 48 – 72 hours to obtain high viral titer in the supernatant (*usually titers of  $10^5 - 10^6$  cfu/mL*).
5. Collect the viral supernatant and centrifuge by low speed (about 500 x g) for 5 min at 4°C or pass it through a 0.45  $\mu\text{m}$  cellulose acetate or PES filter for removing cell debris. (Do not use a nitrocellulose filter because the virus is destroyed). Aspirate the supernatant and storage it at -80°C or infect target cells.

### Viral Infection of Adherent Cells (Protocol for 6-well plate)

1. **Plate the cells**, the day before infection, at a cell density of 60-70% confluence. The number of cells recommended for 6 well plates are  $2.5-5.0 \times 10^5$  cells / well.
2. On the day of infection, removed from the culture medium of the well and **add viral supernatant**. Use 2ml viral supernatant/well or if you know the viral titer, use the volume of viral supernatant which you need and complete with complete medium up to 2ml/well. If you need improve yours results, you can do multiple rounds of infection. Allow cells to rest for 12-24 hr between each infection.
3. **Add polybrene** to a final concentration of 8  $\mu\text{g}/\text{mL}$ . Polybrene range is 2 – 12  $\mu\text{g}/\text{mL}$ . Could be necessary optimized this parameter.
4. **Replace medium** with fresh medium after 12-24h of infection. You can replace medium after 4 – 6hr after infection if your target cells require it.
5. After 48- 72hr infection you can analyse the results of infection process.

### Viral Infection of Suspension Cells (Protocol for 6-well plate)

1. **Collect the cells** by centrifugation (250g for 5 min) and **resuspend**  $1 \times 10^6$  cells in 1 mL of the filtered viral supernatant or if you know the viral titer, use the volume of viral supernatant which you need and complete with complete medium up to 1ml.
2. **Add polybrene** to a final concentration of 8  $\mu\text{g}/\text{mL}$  and dispense the mix into a well of the 6-well plate. Polybrene range is 2 – 12  $\mu\text{g}/\text{mL}$ . Could be necessary optimized this parameter.
3. **Centrifuge the plate** at 1000g, 1 hr 32°C (32°C increases viral half-life).
4. **Incubate** the plate for 12-24 hr.
5. **Collecting cells and replace medium** with 2mL of fresh medium by centrifugation (250g for 5min) and **incubate** the cells at 37°C, 5%  $\text{CO}_2$ .
6. After 48- 72hr infection you can analyse the results of infection process.

**Restricted Use:** For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Distributed by:

**Lab Unlimited**  
CARL STUART GROUP

Tallaght Business Park  
Whitestown, Dublin 24,  
Ireland  
D24 RFK3

Tel: (01) 4523432  
Fax: (01) 4523967  
E-mail: info@labunlimited.com  
Web: www.labunlimited.com

Quatro House, Frimley Road,  
Camberley,  
United Kingdom  
GU16 7ER

Tel: 08452 30 40 30  
Fax: 08452 30 50 30  
E-mail: info@labunlimited.co.uk  
Web: www.labunlimited.co.uk