
## **OPERATING MANUAL**

ba77094e01 06/2015



# photoLab® 7600 UV-VIS

**SPECTROPHOTOMETER** 

Distributed by:



Tallaght Business Park Whitestown, Dublin 24, Ireland D24 RFK3

Tel: (01) 4523432 Fax: (01) 4523967 Quatro House, Frimley Road, Camberley, United Kingdom GU16 7ER

Tel: 08452 30 40 30 Fax: 08452 30 50 30 E-mail: info@labunlimited.com E-mail: info@labunlimited.co.uk Web: www.labunlimited.com Web: www.labunlimited.co.uk



a xylem brand



The most current version of the present operating manual is available on the Internet at <u>www.WTW.com.</u>

Copyright© Weilheim 2015, WTW GmbH<br/>Reproduction in whole - or even in part - is prohibited without the<br/>express written permission of WTW GmbH, Weilheim.<br/>Printed in Germany.

## Contents

1	Ove	Overview				
	1.1	Overvi	ew of the instrument	7		
	1.2	Keypa	d the second s	8		
	1.3		/			
2	Safe	ety		11		
	2.1	Safety	information	11		
		2.1.1	Safety information in the operating manual 1	11		
		2.1.2	Safety signs on the product 1			
		2.1.3	Further documents providing safety information 1	11		
	2.2	Safe o	peration	12		
		2.2.1	Authorized use	12		
		2.2.2	Requirements for safe operation 1	12		
		2.2.3	Unauthorized use 1	12		
	2.3	User q	ualification	12		
	2.4	Handli	ng of hazardous substances 1	13		
	2.1		5			
3			ັ oning			
3		nmissio	oning	14		
3	<b>Con</b> 3.1	n <b>missic</b> Scope	oning	<b> 4</b> 14		
3	<b>Con</b> 3.1 3.2	n <b>missic</b> Scope Genera	oning	1 <b>4</b> 14 14		
3	<b>Con</b> 3.1	n <b>missic</b> Scope Genera Initial c	oning	1 <b>4</b> 14 14 15		
3	<b>Con</b> 3.1 3.2	n <b>missic</b> Scope Genera Initial c 3.3.1	oning       1         of delivery       1         al notes on handling       1         commissioning       1         Inserting the buffer batteries       1	14 14 15 15		
3	<b>Con</b> 3.1 3.2	n <b>missic</b> Scope Genera Initial c	oning       1         of delivery       1         al notes on handling       1         commissioning       1         Inserting the buffer batteries       1         Connecting the power supply       1	14 14 15 15		
3	<b>Con</b> 3.1 3.2	nmissic Scope Genera Initial c 3.3.1 3.3.2	oning       1         of delivery       1         al notes on handling       1         commissioning       1         Inserting the buffer batteries       1         Connecting the power supply       1         Switching on the photometer for the first time       1	1 <b>4</b> 14 15 15 16		
3	<b>Con</b> 3.1 3.2	Scope Genera Initial c 3.3.1 3.3.2 3.3.3	oning       1         of delivery       1         al notes on handling       1         commissioning       1         Inserting the buffer batteries       1         Connecting the power supply       1	1 <b>4</b> 14 15 15 16 17		
3	<b>Con</b> 3.1 3.2	nmissic Scope Genera Initial c 3.3.1 3.3.2 3.3.3 3.3.4 3.3.5	oning       1         of delivery       1         al notes on handling       1         commissioning       1         inserting the buffer batteries       1         Connecting the power supply       1         Switching on the photometer for the first time       1         Setting the language       1         Setting the date and time       1	14 14 15 15 16 17 17		
3	<b>Con</b> 3.1 3.2 3.3	nmissic Scope Genera Initial c 3.3.1 3.3.2 3.3.3 3.3.4 3.3.5	oning       1         of delivery       1         al notes on handling       1         commissioning       1         Inserting the buffer batteries       1         Connecting the power supply       1         Switching on the photometer for the first time       1         Setting the language       1	14 14 15 15 16 17 18		
3	<b>Con</b> 3.1 3.2 3.3	nmissic Scope Genera 1nitial c 3.3.1 3.3.2 3.3.3 3.3.4 3.3.5 Conne	oning       1         of delivery       1         al notes on handling       1         commissioning       1         inserting the buffer batteries       1         Inserting the buffer batteries       1         Connecting the power supply       1         Switching on the photometer for the first time       1         Setting the language       1         Setting the date and time       1         cting optional accessories       1	14 14 15 15 16 17 18 18		
3	<b>Con</b> 3.1 3.2 3.3	Conne 3.3.1 3.3.2 3.3.3 3.3.4 3.3.5 Conne 3.4.1	oning       1         of delivery       1         al notes on handling       1         commissioning       1         inserting the buffer batteries       1         Connecting the power supply       1         Switching on the photometer for the first time       1         Setting the language       1         Setting the date and time       1         communication interfaces       1	14 14 15 15 16 17 18 18 18 20		
3	<b>Con</b> 3.1 3.2 3.3	nmissic Scope Genera Initial c 3.3.1 3.3.2 3.3.3 3.3.4 3.3.5 Conne 3.4.1 3.4.2 3.4.3 3.4.4	oning       1         of delivery       1         al notes on handling       1         commissioning       1         inserting the buffer batteries       1         Connecting the power supply       1         Switching on the photometer for the first time       1         Setting the language       1         Setting the date and time       1         cting optional accessories       1         Communication interfaces       1         PC/printer       2         USB memory device       2         PC keyboard       2	<b>14</b> 14 15 15 16 17 18 18 20 21		
3	<b>Con</b> 3.1 3.2 3.3	nmissic Scope Genera 1nitial c 3.3.1 3.3.2 3.3.3 3.3.4 3.3.5 Conne 3.4.1 3.4.2 3.4.3	oning       1         of delivery       1         al notes on handling       1         commissioning       1         inserting the buffer batteries       1         Connecting the power supply       1         Switching on the photometer for the first time       1         Setting the language       1         Setting the date and time       1         cting optional accessories       1         Communication interfaces       1         PC/printer       2         USB memory device       2	<b>14</b> 14 15 15 16 17 18 18 20 21 21		

4	Оре	eration.		23	
	4.1	Switching on or off the photometer			
	4.2	Genera	l operating principles	25	
		4.2.1	Navigating with function keys and menus		
		4.2.2	Display of navigation paths in short form	26	
		4.2.3	Entry of numerals, letters and characters	27	
		4.2.4	Detailed operating example: Changing the language	29	
	4.3	Photom	eter settings and system administration	30	
		4.3.1	Language		
		4.3.2	Date/Time	30	
		4.3.3	Display settings	31	
	4.4	Zero ad	ljustment	32	
	4.5	Measur	ing in <i>Concentration</i> mode	37	
	-	4.5.1	Measuring cell tests with barcode		
		4.5.2	Measuring reagent tests with AutoSelector		
		4.5.3	Measuring reagent-free tests and user-defined methods	39	
		4.5.4	Exceeding the upper or lower limits of the measuring range		
		4.5.5	Selecting a method manually		
		4.5.6	Settings for <i>Concentration</i> mode	44	
		4.5.7	Measuring diluted samples	45	
		4.5.8	Sample blank value		
		4.5.9	Reagent blank value		
		4.5.10	User calibration (standard adjustment)		
			Automatic Turbidity correction		
			Programming / modifying user-defined methods		
			The <i>IQ-LabLink</i> procedure		
	4.6		ing the Absorbance / % Transmission		
		4.6.1	General information		
		4.6.2	Measuring the absorbance or transmission		
		4.6.3	Measuring against the Reference absorbance		
	4.7		/ Multi wavelengths methods		
		4.7.1	Basic information on Special / Multi wavelengths measurements		
		4.7.2	Programming / modifying the Special / Multi wavelengths methods		
		4.7.3	Selecting a Special / Multi wavelengths method		
		4.7.4	Carrying out Special / Multi wavelengths measurements		
	4.8		measurement		
		4.8.1	General information		
		4.8.2	Selecting OptRF measurement		
		4.8.3	Carrying out measurements		
		4.8.4	Setting the dilution		
		4.8.5	User calibration		
	4.9		im		
		4.9.1	General information		
		4.9.2	Recording a Spectrum		
		4.9.3	Loading/editing a spectrum		
		4.9.4	Saving / exporting a spectrum	104	

4.10	Kinetics	
	4.10.1 Creating/editing profiles for Kinetics recordings	
	4.10.2 Loading a profile for kinetic recording	
	4.10.3 Recording the Kinetics	
	<ul><li>4.10.4 Saving / exporting a Kinetics record</li></ul>	
	4.10.6 Editing a kinetic record	
4 11	Timer	
7.11	4.11.1 User defined timer	
	4.11.2 Analysis timer	
4.12	Memory	
	4.12.1 Overview	
	4.12.2 Instructions on using USB memory devices	122
	4.12.3 Measurement datasets	
	4.12.4 Saving measurement datasets manually	
	4.12.5 Saving measurement datasets automatically	
	4.12.6 Displaying measurement data memory	
	<ul><li>4.12.7 Filtering measurement datasets</li></ul>	
	<ul><li>4.12.8 Inverting filters</li></ul>	
	4.12.10 Saving kinetic recordings, spectra and AQA files	
	4.12.11 Saving data as a pdf file	
4.13	Saving / exporting files	
	4.13.1 Copying all measurement data files to a USB memory device	
	4.13.2 Copying user-defined methods / profiles to a USB memory device	
	4.13.3 Copying files to a PC	133
4.14	Importing files	
	4.14.1 Importing spectra or kinetic recordings from a USB memory device	
	4.14.2 Importing methods / profiles from a USB memory device	
	4.14.3 Importing files from a PC	
4.15	Printing the data (USB)	
	4.15.1 Printer and terminal programs	
	<ul><li>4.15.2 Settings for data transmission</li></ul>	
	4.15.3 Printing measurement datasets	
	4.15.5 Printing Kinetics records	
4 16	Analytical quality assurance (AQA)	
	4.16.1 General information	
	4.16.2 Photometer monitoring (AQA1)	
	4.16.3 Total system monitoring (AQA2)	
	4.16.4 AQA3/MatrixCheck	148
4.17	User management	
	4.17.1 User levels and user rights	
	4.17.2 Activating or deactivating the User management function	
	4.17.3 Creating, changing or deleting a user account	
	<ul><li>4.17.4 Login with active user management</li></ul>	

	4.18	Reset	161			
	4.19	Photometer information ([Info])				
	4.20	Lamp counter	162			
	4.21	Software and methods update				
		4.21.1 Software- and method update using a USB memory device				
		4.21.2 Software- and method update via PC				
5	Mair	ntenance and cleaning				
•	5.1	Exchanging the buffer battery				
	5.2	Cleaning				
	0.2	5.2.1 Cleaning the enclosure				
		5.2.2 Cleaning the cell shaft				
		5.2.3 Cleaning the detector lens	168			
6	Wha	at to do if				
•	6.1	Actions in the case of a broken cell				
	6.2	Error causes and remedies				
7	Tool	hnical data	174			
1	7.1	Measurement characteristics				
	7.1 7.2	Measured value documentation and quality assurance				
	7.2 7.3	General meter data				
	7.3					
8	Acce	essories and options	180			
	8.1	Accessories	180			
	8.2	Test equipment	181			
	8.3	Optional equipment	181			
	8.4	Connection cable:	181			
9	App	endix				
-	9.1	Glossary				
	9.2	List of trademarks				
	9.3					
	-					

## 1 Overview

**1.1** Overview of the instrument






Fig. 1-1 Front of the instrument with operating elements



#### Keypad 1.2

**Overview** 



- Function keys F1 to F4 (function menu-depending)
   Alphanumeric keypad
   Keys with dedicated function

Fig. 1-3 Keypad


**Key functions** The keys on the right side of the keypad have the following functions:

Key	Designation	Functions
<mark>ل</mark>	<on off=""></on>	<ul> <li>Switches the photometer on and off</li> </ul>
НОМЕ	<home></home>	<ul> <li>Switches to the main menu from any operating situation. Actions that are not completed are can- celed.</li> </ul>
PRINT	<pre>PRINT&gt;</pre>	<ul> <li>Outputs the displayed measured value to an interface if the <i>Printer</i> symbol is displayed in the status line.</li> </ul>
STORE	<store></store>	<ul> <li>Saves a displayed measured value or spectrum if the Save symbol is displayed in the status line.</li> </ul>
ZERO • BLANK	<zero·blank></zero·blank>	<ul> <li>Starts one of the following measurements, depending on the operating situation:</li> <li>Zero adjustment</li> <li>Blank value measurement</li> <li>Baseline measurement</li> <li>User calibration</li> </ul>

Key	Designation	Functions
TIMER	<timer></timer>	- Opens the menu, <i>Timer</i> .
ESC	<esc></esc>	<ul> <li>Cancels the running action.</li> <li>Entries that have not yet been accepted are discarded.</li> </ul>
		<ul> <li>Switches to the next higher menu level.</li> </ul>
START · ENTER	<start.enter></start.enter>	<ul> <li>Starts an action (e.g. measure- ment)</li> </ul>
		<ul> <li>Opens a selected menu</li> </ul>
		<ul> <li>Confirms a selection or entry</li> </ul>
	< <b>▲&gt;</b> or < <b>▼&gt;</b>	<ul> <li>Moves the selection in menus and lists one position up or down</li> </ul>
	< <b>∢</b> >	<ul> <li>When entering characters, deletes the character left of the cursor</li> </ul>
(arrow keys)		<ul> <li>Moves the cursor to the left in a spectrum or kinetic diagram</li> </ul>
	<▶>	<ul> <li>Moves the cursor to the right in a spectrum or kinetic diagram</li> </ul>

**Function keys** The function keys F1 to F4 have different functions depending on the operating situation. The current functions are displayed in the function key menu at the bottom edge of the display (see section 4.2.1).

#### 1.3 Display



- 1 Status line (current state, date and time)
- 2 Display range for menus and measurement results3 Function keys menu

Fig. 1-4 Display

Symbols in the	Symbol	Designation	Function
status line		Save	The <b><store></store></b> key is active. You can store the displayed data with <b><store></store></b> (see section 4.12).
	ð	Printer	The <b><print></print></b> key is active. You can output to an interface the displayed data with <b><print></print></b> (see section 4.15).

## 2 Safety

#### 2.1 Safety information

#### 2.1.1 Safety information in the operating manual

This operating manual provides important information on the safe operation of the product. Read this operating manual thoroughly and make yourself familiar with the product before putting it into operation or working with it. The operating manual must be kept in the vicinity of the product so you can always find the information you need.

Important safety instructions are highlighted in this operating manual. They are indicated by the warning symbol (triangle) in the left column. The signal word (e.g. "CAUTION") indicates the level of danger:



#### WARNING

indicates a possibly dangerous situation that can lead to serious (irreversible) injury or death if the safety instruction is not followed.

### CAUTION

indicates a possibly dangerous situation that can lead to slight (reversible) injury if the safety instruction is not followed.

#### NOTE

indicates a situation where goods might be damaged if the actions mentioned are not taken.

#### 2.1.2 Safety signs on the product

Note all labels, information signs and safety symbols on the product. A warning symbol (triangle) without text refers to safety information in this operating manual.

#### 2.1.3 Further documents providing safety information

The following documents provide additional information, which you should observe for your safety when working with the measuring system:

- Operating manuals of other components of the photoLab<sup>®</sup> 7600 UV-VIS (accessories)
- Safety datasheets for chemicals.

#### 2.2 Safe operation

#### 2.2.1 Authorized use

The authorized use of the photometer consists exclusively of the carrying out of photometric measurements according to this operating manual. Follow the technical specifications of the cells in chapter 7 TECHNICAL DATA. Any other use is considered **unauthorized**.

#### 2.2.2 Requirements for safe operation

Note the following points for safe operation:

- The product may only be operated according to the authorized use specified above.
- The product may only be supplied with power by the energy sources mentioned in this operating manual.
- The product may only be operated under the environmental conditions mentioned in this operating manual.
- The product may not be opened.

#### 2.2.3 Unauthorized use

The product must not be put into operation if:

- it is visibly damaged (e.g. after being transported)
- it was stored under adverse conditions for a lengthy period of time (storing conditions, see chapter 7 TECHNICAL DATA).

#### 2.3 User qualification

Carrying out photometric determinations with the aid of test sets frequently requires the handling of hazardous substances.

We assume that the operating personnel know how to handle hazardous substances due to their professional training and experience. The operating personnel must particularly be able to understand and correctly implement the safety labels and safety instructions on the packages and inserts of the test sets.

#### 2.4 Handling of hazardous substances

When developing test sets, WTW carefully sees that the tests can be carried out as safely as possible. Some hazards by dangerous substances, however, cannot always be avoided.

If self-produced tests or solutions are used, the responsibility concerning any risks caused by those tests or soulutions lies with the user (personal responsibility).



#### WARNING

Improper handling of certain reagents can cause damage to your health.

In any case follow the safety labels on the packing and the safety instructions of the package insert. Protective measures specified there have to be followed exactly.

Safety datasheets

The safety datasheets of the chemicals comprise all instructions on safe handling, occurring hazards, preventive actions and actions to take in hazardous situations. Follow these instructions in order to work safely.

## 3 Commissioning

#### 3.1 Scope of delivery

- Spectrophotometer photoLab® 7600 UV-VIS
- Power pack with connection cable
- Buffer batteries 4 x AA alkaline manganese (Mignon)
- Zero cell (16 mm, round)
- Short instructions
- CD-ROM with
  - detailed operating manual
  - Analysis instructions
  - SpectralTransfer software

Packing This photometer is sent out in a protective transport packing.



#### CAUTION

Keep the original packing including the inner packing to protect the instrument against hard shocks if it has to be transported.

The original packing is also required for the proper return of the instrument if it has to be repaired.

Note that damage caused by improper transport voids all warranty claims.

#### 3.2 General notes on handling

The photoLab<sup>®</sup> 7600 UV-VIS photometer is an optical precision meter. Therefore, it should always be handled with care, especially in mobile use. Always protect the meter from conditions that could damage the mechanical, optical and electronic components. Heed the following points especially:

- The temperature and humidity during operation and storage must be within the limits specified in chapter 7 TECHNICAL DATA.
- The following influences always have to be avoided with the meter:
  - Extreme dust, moisture and wetness
  - Exposure to intensive light and heat
  - Fumes that are corrosive or contain high concentrations of solvents.
- For measuring, the meter must be placed on a flat surface.
- Spilled liquid or other material should be removed immediately (see section 5.2 CLEANING).
- If a cell has broken in the cell shaft, the cell shaft should be cleaned immediately (see section 6.1 ACTIONS IN THE CASE OF A BROKEN CELL).

- The cell shaft should always be closed when the photometer is not used.
- For the instrument to be transported the cell shaft has to be empty.
- For mobile used we recommend using the suitable transport case FC spectral 6000 (see section 8.1 ACCESSORIES).

#### 3.3 Initial commissioning

Perform the following activities:

- Insert the buffer batteries (see section 3.3.1)
- Connect the power supply (see section 3.3.2)
- Switch on the photometer (see section 3.3.3)
- Set the language (see section 3.3.4)
- Set the date and time (see section 3.3.5)
- Carry out a zero adjustment (see section 4.4)




When you set the language, date and time according to the mentioned sections of this operating manual you will quickly become familiar with the simple operation of the 3.3.43.3.5photoLab<sup>®</sup> 7600 UV-VIS. More detailed instructions on operation are given in section 4.2 GENERAL OPERATING PRIN-CIPLES.

#### 3.3.1 Inserting the buffer batteries

The buffer batteries supply the integrated clock while the photometer is switched off. Four alkaline manganese batteries (type AA or Mignon) separately included in the scope of delivery are used as the buffer batteries.

Insert the batteries as follows:



- **1** Turn the photometer upside down and place it on a soft surface.
- 2 Open the lid of the battery compartment (1).
- 3 Insert the four batteries in the battery compartment. Make sure that the poles of the batteries are in the correct position.

The  $\pm$  signs on the batteries must correspond to the  $\pm$  signs in the battery compartment.

4 Close the lid of the battery compartment.

**Battery service life** The power consumption of the clock is very low. The lifetime of high quality batteries is at least 5 years.

#### 3.3.2 Connecting the power supply


The power is supplied via the enclosed plug-in power pack. The power pack supplies the photometer with low voltage (12 VDC).



#### CAUTION

The line voltage of the usage location must fulfill the specifications stated on the power pack (the specifications are also given in chapter 7 TECHNICAL DATA). Always use the supplied 12 V original power pack only.

Before plugging in the power cable check whether it is undamaged. If the power cable is damaged, the instrument must not be operated.



Operation with a mobile 12 V power source You can also operate the photoLab<sup>®</sup> 7600 UV-VIS on the move and independent of the local power supply.

To do so, a 12 V power supply unit such as a commercial 12 V power station or a 12 V car battery and a 12 V vehicle connection cable is required (see section 3.4.6).

#### 3.3.3 Switching on the photometer for the first time

During the initial commissioning, the photometer automatically guides you through the setting of the meter language, date and time after switching on (see following sections).

Gen	eral setup	16.01.15 9:52
	German	
1	English	
	Français	
	Español	
	Italiano	
	Bulgarian/Български	
	Česko	
	Simplified Chinese/ 中	
	Traditional Chinese/ 繁	
	Greek/Ελληνικά	
		1

1 Press <ON/OFF>.

The photometer is switched on.

The display switches to the setting of the language (see section 3.3.4).

After the setting of the language the photometer carries out the self-test.

When the initial commissioning is completed, the photometer displays the *Home* menu each time after it is switched on and after the self-test (see section 4.1).

#### 3.3.4 Setting the language

During the initial commissioning the photometer automatically guides you to the setting of the meter language after switching on.



- **1** Select a language with  $< \Delta > < \nabla >$ .
- 2 Confirm the selected language with <**START·ENTER**>.

The language has been set. The currently selected language is marked by a check.

The display switches to the setting of the *Date* and *Time* (see section 3.3.5).

After the initial commissioning, you can change the language in the *General* setup / Language menu at any time (see section 4.2.4).

#### 3.3.5 Setting the date and time

During the initial commissioning, the instrument automatically guides you to the setting of the time and date after the setting of the language.



Date/Time	16.01.15 9:52
Date Time	16.01.2015 9:52:09
Time	9:52:09
Date	
16 .01.2015	
	ОК

 Date/Time
 16.01.15
 9:52

 Date
 16.01.2015

 Time
 9:52:09

 09
 : 52
 : 09

The Date/Time menu is open.

Using <▲><▼>, select a menu item and confirm or open it with <**START-ENTER**>.

- Select and confirm *Date*.
   The input field for the current date pops up.
- 2 Enter the current date with <0...9> and confirm.

The input field closes. The date is accepted.

- 3 Select and confirm *Time*. The input field for the current time pops up.
- 4 Enter the current time with <0...9> and confirm.

The input field closes. The time is accepted.

After the initial commissioning, you can change the date and time in the *General setup / Date/Time* menu at any time (see section 4.2.4).

#### 3.4 Connecting optional accessories

3.4.1 Communication interfaces

#### Connections




Fig. 3-1 Communication interfaces on the rear panel

You can connect the following accessories to the photometer:

- PC (see section 3.4.2)
- Printer (see section 3.4.2)
- USB storage media (see section 3.4.3)
- USB-PC keyboard (see section 3.4.4)
- Barcode reader (see section 3.4.5)
- 12 V vehicle charging cable (see section 3.4.6)



If you want to connect several USB devices such as a USB-PC keyboard and a USB memory device to the meter, you can increase the number of USB-A sockets by a commercially available USB-2 hub with separate power supply.

#### 3.4.2 PC/printer

PC and printer can be connected to the photometer as follows:

Interface	PC	Printer	Functions
USB-A		1	The data is printed out with <b><print></print></b> .
USB-B	1	-	Enables the direct connection of photom- eter and PC. With this you can transmit measurement data to the PC (see sec- tion 4.13 and section 4.15) or update the photometer software (see section 4.21.1).
			The direct connection with the PC is established with the aid of the "Spectral- Transfer" program. The program is pro- vided on the supplied CD-ROM.
			More instructions on how to establish the connection are given in the operating manual of the "SpectralTransfer" pro- gram (see CD-ROM).



Suitable are all printers that can interpret the PCL-3 printer control language.

#### 3.4.3 USB memory device

Using a USB memory device (such as a USB flash drive), you can

- update the meter software and method data (section 4.21)
- transmit data to the USB memory device (section 4.12 and section 4.13).

USB memory devices are connected to the USB-A interface.



Please follow the instructions on using USB memory devices (see section 4.12.2).

#### 3.4.4 PC keyboard

With the PC keyboard it is possible to enter letters, e.g. to assign names for identification (ID).

In addition, the following keys of the PC keyboard are assigned with the following functions of the photometer:

Photometer	PC keyboard
<start·enter></start·enter>	Enter
<esc></esc>	Esc
<f1> to <f4> (function keys)</f4></f1>	F1 to F4
< <b>▲</b> >< <b>▼</b> >< <b>♦</b> > (arrow keys)	Arrow keys
<home></home>	F5
<pre><print></print></pre>	F6
<store></store>	F7
<zero·blank></zero·blank>	F8
<timer></timer>	F9
<on off=""></on>	F12
Symbols and characters according to the operating manual	Corresponding key on the key- board
09	09
-	-

The USB-PC keyboard is connected to the USB-A interface.

#### 3.4.5 Barcode reader

The barcode reader enables the simplified entering of alphanumerical character strings and can be used in all operating situations that require the entry of text or numerals. The barcode reader is connected to the USB-A interface.



Before using the bar code reader with the photometer, configure the bar code reader to use code 39 (see operating manual of your bar code reader).

#### 3.4.6 Operation with a 12 V vehicle connection cable

With a commercial 12 V vehicle connection cable, you can operate the photoLab<sup>®</sup> 7600 UV-VIS spectrophotometer on the move and independent of the local power supply.

To do so, a 12 V power supply such as a commercial 12 V portable power source or a 12 V car battery is required.

	12 V vehicle connection cable	
	12 V power supply unit (e.g. portable power source or car battery)	photoLab <sup>®</sup> 7600 UV-VIS
Safety	battery. Make sure the power source is su	ttery, follow the safety instructions of the uitable for operation of the spectrophotome- er source and technical data of the spectro-
Operating time with battery	<ul> <li>The maximum operating time dep</li> <li>Battery (e.g. nominal capacity)</li> <li>Operating mode of the spectro ments)</li> </ul>	
	<ul><li>Photometer (instrument type)</li></ul>	
Example	Operating time with a type 12 V / $^{-1}$ 16 h	19 Ah battery in optimum condition: approx.
	standby mode. With battery operation	er consumes electricity even while it is in , we recommend you disconnect the vehi- vhile you are not using the photometer.
12 V vehicle connection cable	Vehicle connection cables with th	ne following characteristics are suitable:

Voltage	12 V
Amperage	8 A
Barrel connector	2.5 x 5.5mm
Positive pole at the inner contact	

## 4 Operation

4.1 Switching on or off the photometer

#### Switching on

Starting the

Self test

Self test

Then press <START/ENTER>

Self test							16.01.15	9:!
Please ma	ke sur	e no cel	lisin	serted	and ti	ne c	over is clo	sed
Then pres	s <st< td=""><td>ART/ENT</td><td>ER&gt;</td><td></td><td></td><td></td><td></td><td></td></st<>	ART/ENT	ER>					
		,	,					
Setup							Info	)
						П		
Login							16.01.15	9:5
		er name						
Ent	er use							
	minist							

Please make sure no cell is inserted and the cover is closed.

1 Switch the photometer on with **<ON/OFF>**.

The display shows

- the *Self test* dialog (if the user management is not active).

or

 the Login dialog (if the user management is active).

With activated user management:

2 Login

Enter user name and password or register as a guest (see section 4.17.4).

Then the photometer displays the *Self test* dialog.

- **3** Remove all cells and close the cell shaft cover.
- 4 Start the self-test with <**START·ENTER**>.

The photometer carries out the self-test.

**Self test** During the self-test, all cells must be removed and the cell shaft cover closed. The self-test runs in the background and may take some minutes.

16.01.15 9:52

Self test			16.01.15	9:52
Кеер соче	er closed			
System te Filter test Lamp test Waveleng				

The self-test includes:

- the test of the memory, processor, internal interfaces, filter and lamp
- a calibration for each wavelength

After the self-test is completed, the main menu is displayed.




The result of the self-test can be viewed and printed with the *[Info]* function key (see section 4.19).

Automatic wavelength	With the automatic wavelength calibration, the photometer checks and calibrates the accuracy of the wavelengths created (by the monochromator).
calibration	The wavelength calibration of the photometer takes place regularly after the photometer was switched on (within the framework of the self-test) and is automatically repeated during operation after 15, 30, 60, 120 and 240 minutes.
	While the automatic wavelength calibration is running on the photometer, a note is displayed. The automatic wavelength calibration only starts if the cell shaft is empty.
	If a cell is in the cell shaft the wavelength calibration is carried out only after the cell was removed.
Switching off	To switch the photometer off, keep the <b><on off=""></on></b> key depressed until the

**Switching off** To switch the photometer off, keep the **<ON/OFF>** key depressed until the photometer is switched off.

#### 4.2 General operating principles

#### 4.2.1 Navigating with function keys and menus



Use of the function keys The function keys F1 to F4 are below the display. Their functions change depending on the operating situation and mode. The current functions are displayed in the function key menu at the bottom edge of the display.

Apart from navigation, the function keys are also used for other operations:

- Opening a selection list or input field
- Executing a command (directly or with intermediate query)
- Switching over between two display options, such as absorbance ↔ transmission

Navigation with arrow keys (<▲><▼><◀><►>) and <START·ENTER> These operating elements are used to select an item from a menu or list. The current selection is displayed in reverse video. Pressing of **<START·ENTER>** confirms the selection.

Apart from navigation, the **<START·ENTER>** key is also used for other operations:

- Opening a selection list or input field
- Confirming a selection
- Confirming entries of text and numerals
- Executing a command (directly or with intermediate query)
- Activating an item in a selection list ( $\checkmark$  = active)


#### 4.2.2 Display of navigation paths in short form

In this operating manual, the introductory navigation steps leading to individual menus or dialogs are clearly shown in a gray box. The box indicates a section of the menu tree.

Starting point of the description is always the main menu, which can be reached with the **<HOME>** key from any operating situation. From there navigation takes place downward.

The following example shows the elements of the menu tree with the relevant operating steps:

Operating example: Navigation to the setting menu for the language



Further navigation options:

- The **<ESC>** key moves one level up in the menu tree.
- The **<HOME>** key directly calls up the main menu.



If you are "lost" in a menu, press **<HOME>** and restart navigating from the main menu.

#### 4.2.3 Entry of numerals, letters and characters

Numerals, letters, punctuation marks and special characters are entered with the alphanumeric keypad of the meter or using an external keyboard. Entries are required in operating situations such as the following:

- Entering the date and time
- Entering an ID e.g. when storing measurement data

- Selecting a method with the [Search] function
- Programming user-defined methods
- Entering user name and password
- Administrating users

**Character set** The following characters are available:

- Numerals 0 ... 9
- Letters A ... Z and a ... z
- Punctuation marks. -
- Special characters ° / + <sup>2</sup> <sup>3</sup> # %

#### **Operating principle** Entering characters is always possible if there is an input field on the display.



The numerals and characters (expect for the small letters) assigned to the keys of the alphanumeric keypad are printed on the keys. Example: With the **<7/PQRS>** key you can enter the following characters: 7, P, Q, R, S, p, q, r, s.

Select the required character by pressing the key several times (similar to a mobile phone). When pressing a key that is assigned to several characters once, the respective numeral appears first. To enter a numeral, one keypressing is always sufficient.

When pressing the key for the first time a line pops up that displays all characters possible with this key. The currently selected character is highlighted.

A character is taken over in the input field if

- the character is highlighted for more than one second,
- the character is confirmed with <START.ENTER>,
- another alphanumeric key is pressed.



During mere number entries (such as entering a wavelength), the keys of the alphanumeric keypad are assigned to the respective numeral only. Each keypressing directly enters the numeral (like a pocket calculator).

Special characters Sp

Special characters are entered with the <1/\*> key.


Operating example: Entering the ID The *Enter ID* input field appears if you press the **<STORE>** key while the storing symbol is visible. In the following example a measurement dataset with the ID "Test" is stored.

Enter ID 8	1	Press < <b>8/TUV</b> > several times until "T" appears in the input line.
8 TUVtuv Enter 1D T 8 T TUVtuv		Below the input field, a selection line pops up with all characters that are available for this key, e.g. B T U V t u v. The currently selected character is highlighted.
		After approx. one second the character is taken over and the selection line closed.
Enter ID Test_	2	Complete the ID with <b><a9></a9></b> and confirm.

#### Correcting incorrect entries

Using <<>>, erase all characters until you have reached the incorrect digit and repeat the entry from there.

#### 4.2.4 Detailed operating example: Changing the language



General setup	16.01.15	9:52						
Language								
Date/Time								
Display settings								
User management								
Measured value memory								
Software/methods update								
Reset								
Data transfer/Printer								
Exchange methods/profiles								
Save data to USB memory device	Save data to USB memory device							
Unlock application packages								

- 1 Call up the main menu with the **<HOME>** key.
- 2 Open the *General setup* menu with the F1 function key [Setup].
- 3 Using <**▲**><**▼**>, select the *Language* menu item and open with <**START·ENTER**>.

The *Language* menu shows a list with the available languages. The currently active language is marked by a check.

Lan	guage	16.01.15	9:52						
	German								
<	English								
	Français								
	Español								
	Italiano								
	Bulgarian/Български								
	Česko								
	Chinese/ 中文								
	Traditional Chinese/ 繁體中文								
	Greek								
	Indonesian								

4 Select the required language from the list with <▲><▼> and confirm with <**START**·ENTER>.

The selected language is taken over immediately. The photometer moves up one menu level.

#### 4.3 Photometer settings and system administration

The general photometer settings are done in the **<HOME>** -> *General setup* menu. These comprise:

- Language (see section 4.3.1)
- Date/time (see section 4.3.2 and section 4.2.4)
- Display characteristics (see section 4.3.3)
- User management (see section 4.17)
- Administration of the measurement data memory (see section 4.12)
- Software and method update (see section 4.21)
- Reset of the settings to default values (see section 4.18)
- Settings for data transmission (see section 4.15.2)

#### 4.3.1 Language


The complete list of the available instrument languages is given in the *Language* chapter and in chapter 7 TECHNICAL DATA menu of the photometer.



How to set the language is described in detail in the operating example in section 4.2.4.

#### 4.3.2 Date/Time

The date format is set automatically with the language setting. According to the locally usual version, the date format is displayed in the order, Day.Month.Year (*DD.MM.YY*) or Month/Day/Year (*MM/DD/YY* or MM.DD.YY).



#### 4.3.3 Display settings

Here you can adjust the display contrast to the lighting conditions.

<home></home>	
[General setup]	
<ul> <li>Display settings</li> </ul>	

Display	/ setting	s		16	5.01.15	9:52
Brightr	ie ss			50	) %	
			 		_	
	1				1	

- Select and confirm *Brightness*.
   A slide control for the display brightness appears.
- 2 Set and confirm the display brightness with <**4**><**>**>.

#### 4.4 Zero adjustment

A valid zero adjustment is required for the calculation of measured values in the modes, *Concentration, Absorbance / % Transmission, Special / Multi wavelengths* and *Kinetics*. With a zero adjustment, the absorbance of a cell filled with distilled water ("zero cell") is measured and stored as the zero value.

Factory zero adjustment for concentration measurements For all measurements with WTW test sets (*Concentration* mode), a factory zero adjustment is available in the delivery condition. We recommend replacing it with a zero adjustment of your own. If a zero adjustment exists already for a method, the date and time of the last zero adjustment are displayed in the top right area of the display.

Please select method for me	[ZERO	15.01.15 11:11]
Please select method for me		
or insert a barcoded cell or insert AutoSelector.	easuring	

Zero adjustment for absorbance measurements In the *Absorbance* mode, the zero adjustment has to be carried out separately for each cell type and each used wavelength. If a zero adjustment exists already for the inserted cell type at the selected wavelength, the date

and time of the last zero adjustment are displayed in the top right area of the display.

Absorbance	16.01.15 9:52
[ZERO	15.01.15 11:11]
To start measurement, insert cell or press <start enter=""></start>	
525 nm	10 mm
Setup Wavelength Transmission	Reference

If no zero adjustment is available, the photometer will prompt you to carry out a zero adjustment.



The cells must be absolutely clean and free of scratches. Always use a cell of the same type for zero adjustment and measurement of the sample.

#### Notes on zero adjustment

Zero adjustment with round cells:

- Only use clean, scratch-free round cells with distilled water. The minimum filling level is 20 mm. A ready zero cell is included in the scope of delivery of the photometer and PhotoCheck (see chapter 8 ACCESSORIES AND OPTIONS).
- A ready zero cell can, in principle, be used any number of times. We recommend, however, to regularly check the zero cell for visible contamination and scratches and refill or exchange it if necessary (at least every 24 months).

Zero adjustment with rectangular cells:

- For rectangular cells, the zero adjustment must be carried out with the same cell type (manufacturer and glass type [e.g. optical glass, quartz glass, plastic]) that is used for measurement. This is important because cells of different manufacturers have a different absorption behavior. When changing the cell type repeat the zero adjustment with the new type.
- Prior to zero adjustment, clean the rectangular cell and fill it with distilled water. The minimum filling level is 20 mm.
- Rectangular cells always have to be inserted in the cell shaft with the same orientation for measurement and zero adjustment (e.g. cell printing on the left side ).



Ordering information is given in chapter 8 ACCESSORIES AND OPTIONS. The cells listed in the chapter 8 ACCESSORIES AND OPTIONS are especially adapted to the WTW test set program. General requirements of the cells are given in chapter 7 TECHNICAL DATA. Note that the spectral transparency of the cell must be suitable for the intended application (example, quartz cell for UV range).

#### Carrying out a zero adjustment


The zero adjustment takes place similarly in the *Concentration, Absorbance* / % *Transmission, Special / Multi wavelengths* and *Kinetics* modes.

Concen	ntration			16.	01.15 9:52
	Adjust				
	Zero adj	ustment			
	Reagent	: blank value			
	Calibrate	e the method			
3: A6/2	25				NH4-N
16 m m	l			0.20	- 8.00 mg/l
Set	tun	Method list	Citation form		Unit

Zero adjustme	nt		16.01.15	9:52
	sert zero celi ( <start ente<="" th=""><td>distilled water R&gt;</td><td>)</td><td></td></start>	distilled water R>	)	

- 1 In the respective mode, press the <**ZERO·BLANK**> key.
- 2 In *Concentration* mode only: Select and confirm *Zero adjustment*.

The zero adjustment window pops up.



Zero adjustment		16.01.15	9:52
	Zero adjustment successful		
		1	0 mm
		OK	

- **3** Close the inner turn-up lid.
- 4 Depending on the cell type, insert the zero cell as follows:

#### Round cell:

Insert the round cell in the round cell shaft so that it touches the bottom.

If the inner turn-up lid is opened too wide, a message prompts you to close the inner turn-up lid.

#### Rectangular cell:

Open the inner turn-up lid.

Insert the rectangular cell vertically so that it touches the bottom and left edge of the cell shaft. The opaque sides of the rectangular cell must point to the front and back.

The photometer has an external light recognition. If there is too much external light, a message prompts you to close the cell shaft cover.

The photometer automatically starts the zero adjustment and subsequently stores the value.

5 After a successful zero adjustment switch to measurement with [OK].

#### Validity of the zero adjustment

The data of the zero adjustment is stored in the photometer separately for each cell type. As long as the data is valid, it is automatically used again after a temporary change to a different cell type. The validity depends on the respective mode:

Mode	Validity of the zero adjustment
Absorbance / % Transmission	<ul> <li>Till the next zero adjustment with the same wavelength *</li> </ul>

Mode	Validity of the zero adjustment
<i>Concentration</i> (user-defined methods) and	<ul> <li>Till the next zero adjustment for the same method *</li> </ul>
Special / Multi wavelengths	
Kinetics	Till another kinetic profile is loaded
	• Till the <i>Kinetics</i> mode is exited or the photometer is switched off

\* The photometer displays that a zero adjustment is available and the time it was carried out. You can then decide whether to use this zero adjustment or carry out a new zero adjustment.

When to repeat the zero adjustment?

We recommend to repeat the zero adjustment in the following cases:


- If the photometer was subject to mechanical stress such as strong shock or transport
- If the ambient temperature changed by more than 5 °C since the last zero adjustment
- At least once per week
- If a new cell type (different manufacturer, different glass type is used)
- Basically each time you want to measure with the highest possible accuracy.


# 4.5 Measuring in *Concentration* mode

4.5.1 Measuring cell tests with barcode

# <HOME>







Concentration			16.01.15 9:52
		40	
		49	mg/l
18: C3/25 (44)	5 nm)		COD
16 m m			10 - 150 mg/l
Setup	Method list	Citation form	Unit

Inserting a cell with barcode starts a measurement.

- **1** Open the cell shaft cover.
- 2 Close the inner turn-up lid.

If the inner turn-up lid is open, a message prompts you to close the inner turn-up lid.

3 Insert the barcoded round cell in the round cell shaft so that it touches the bottom. When doing so, align the line mark with the notch at the front of the round cell shaft.

The photometer selects the method based on the bar code and automatically starts measurement.

- 4 Further options:
  - Select a different citation form with [*Citation form*], (e.g. NH₄ <-> NH₄-N).
  - Select a different measuring unit with [Unit], (e.g. mg/l <-> mmol/l).
  - Make further settings such as dilution or blank value measurements with [Setup] (see section 4.5.6).

Display if the measured value is not within the measuring range (see section 4.5.4).


## 4.5.2 Measuring reagent tests with AutoSelector

On one side, the AutoSelector has a bar code encoding a method number. When the AutoSelector is inserted in the round cell shaft, the bar code is automatically read and the relevant method selected by the photometer.


## <HOME> Concentration

Concentration			16.01.15	9:52
or insert or insert	lect method fo a barcoded ce AutoSelector.			
Setup	Method list	Last method	New Met	:ho d

**1** Open the cell shaft cover.



- 2 Insert the AutoSelector in the round cell shaft so that it touches the bottom. When doing so, align the line mark with the notch at the front of the round cell shaft.
  - The photometer selects the correct method with the aid of the barcode.



- **3** Open the inner turn-up lid.
- 4 Insert the rectangular cell vertically so that it touches the bottom and left edge of the cell shaft. The opaque sides of the rectangular cell must point to the front and back.

The correct measuring range is automatically selected when the rectangular cell (1, 2, 5 cm) is inserted.

The photometer has an external light recognition. If there is too much external light, a message prompts you to close the cell shaft cover.

The photometer starts measuring automatically.

- 5 Further options:
  - Select a different citation form with [*Citation form*], (e.g. NH₄ <-> NH₄-N).
  - Select a different measuring unit with [Unit], (e.g. mg/l <-> mmol/l).
  - Make further settings such as dilution or blank value measurements with [Setup] (see section 4.5.6).

Display if the measured value is not within the measuring range (see section 4.5.4).

### 4.5.3 Measuring reagent-free tests and user-defined methods

16.01.15 9:52

NH<sub>4</sub>-N

0.05 - 3.00 ma/

Unit

1.92 ma/l

Method list Citation form

User-defined methods and reagent-free methods normally do not have a barcode and therefore, no automatic method recognition. In such a case, select the method manually:


<HOME> Concentration

Concentration

18: 14752

Setup

10 m m



inner turn-up lid 1 Select the method manually (see section 4.5.5).

The photometer is ready to measure.


2 Depending on the type, insert the cell as follows:

Round cell:

Close the inner turn-up lid.

Insert the round cell in the round cell shaft so that it touches the bottom.

If the turn-up lid is opened too wide, a message prompts you to close the inner turn-up lid.



Concentration			16.01.15	9:52
	$\cap$	<b>` ^ ^</b>		
	U.C	529	ma/	
			51	
15: Color436 (4	136 nm)			DFZ
50 m m			0.1 - 50.0	) m <sup>-1</sup> l
Setup	Method list	Citation form	Unit	

### Rectangular cell:

Open the inner turn-up lid.

Insert the rectangular cell vertically so that it touches the bottom and left edge of the cell shaft. The opaque sides of the rectangular cell must point to the front and back.

The photometer has an external light recognition. If there is too much external light, a message prompts you to close the cell shaft cover.

- 3 Further options:
  - Select a different citation form with [Citation form], (e.g. NH₄ <-> NH₄-N).
  - Select a different measuring unit with [Unit], (e.g. mg/l <-> mmol/l).
  - Make further settings such as dilution or blank value measurements with [Setup] (see section 4.5.6).

Display if the measured value is not within the measuring range (see section 4.5.4).

### 4.5.4 Exceeding the upper or lower limits of the measuring range

Measured values outside the limits of the measuring range are displayed in red.

Measured value display if the measured value is outside the measuring range:

Ra	nge	Display	Example: MR: 10 - 150 mg/l
	LL < MV < UL	Measured value	128 mg/l
1	UL < <b>MV</b> < UL + 10%	Upper limit of measuring range exceeded by up to 10% and measured value	> 150 157 mg/l
	LL - 50% < <b>MV</b> < LL	Lower limit of measuring range undercut by up to 50% and measured value	< 10 7 mg/l
2	<b>MV</b> > UL + 10%	Upper limit of measuring range exceeded by more than 10%	> 150 mg/l
	<b>MV</b> < LL - 50%	Lower limit of measuring range undercut by more than 50%	< 10
3	Invalid measured value	Bars	
	e.g. <b>MV</b> < 0		mg/l

MR = Measuring range

UL = Upper limit value of the measuring range

LL = Lower limit value of the measuring range

MV = Measured value

### 4.5.5 Selecting a method manually

Selecting the method last used

<HOME> Concentration [ [Last method] The method last used is immediately selected.

# Selecting a method from the *Method list*

Sele	ct method	d (all)	16.01.15 9:52
Sear	ch:		
1	C3/25	COD	10 - 150 mg/l
2	C4/25	COD	25 - 1500 mg/l
3	A6/25	NH4-N	0.20 - 8.00 mg/l
4	N2/25	NO <sub>3</sub> -N	0.5 - 25.0 mg/l
5	N5/25	NO <sub>2</sub> -N	0.010 - 0.700 mg/l
6	P6/25	PO <sub>4</sub> -P	0.05 - 5.00 mg/l
7	P7/25	PO <sub>4</sub> -P	0.5 - 25.0 mg/l
160	09717	Pb	0.10 - 5.00 mg/l
67	14834	Cd	0.025 - 1.000 mg/l
39	14552	Cr	0.05 - 2.00 mg/l
La	st used		

The list of methods is displayed. The methods are ordered by method numbers. The arrows ▼ or ▲ on the right edge indicate that the list comprises more methods further up or down.

The method last selected is highlighted.

#### Select the method:

- Select the required method with
   <▲><▼>. The active selection is displayed in reverse video.
- 2 Accept the selection with **START·ENTER**>.

## Narrowing down the method list

You can narrow down the method list and thus make the search easier:

- Using [Last used], you can restrict the method list to the ten methods last used.
- With the search function you can search certain character strings in the list. The search takes place as a full-text search of the entire list contents. Thus you can search for a method number or certain citation form.

### **Search function**

Sele	ect method	l (last used	) 16.01.15 9:52
CO			
14	14540	COD	10 - 150 mg/l
23	14541	COD	25 - 1500 mg/l
All	methods		

Search for a character string:

Enter the character string to be searched for in the search window with **<A...9>**.

The list appearing below shows all hits containing the character string. The hit list is updated with each character that is entered.



Note the case sensitivity when searching. It is not required or possible to enter inferior characters. When searching for chemical formulas, inferior characters are treated as normal characters. Example: The search for "NH4" shows all hits that contain "NH4" as well as "NH<sub>4</sub>".

# 4.5.6 Settings for *Concentration* mode

Prior to measuring, check the settings for the selected method.

<home> Concentrat Select a me</home>	ethod	
Concentration		16.01.15 9:52
		10.01.15 5.52
Sample blank va	due	
User-defined bla		
Turbidity correct		
, Display absorba		
AQA		
Edit method		
New method		
Measurement da	ata memory	

The menu shows an overview of all settings.

Active settings are marked by a check.

Overview of the	Menu item	Explanation
settings	Dilution	Here you can set the dilution prior to measuring if you want to use a diluted sample.
		In the measured value display, the dilution is indicated in the form $[1 + x]$ (parts sample + parts distilled water).
		For more detailed information on dilution, see section 4.5.7.
	Sample blank value	Here you can measure while taking a sample blank value into account.
		In the measured value display, measurements with sample blank value are marked by [SB] (Sample blank).
		For more detailed information on sample blank value, see section 4.5.8.

Menu item	Explanation
User-defined blank value	If available, a user-defined reagent blank value is used.
	In the measured value display, measurements with a user-defined reagent blank value are marked by [BV/Lot number].
	For more detailed information on reagent blank value, see section 4.5.9.
Turbidity correction	Activates/deactivates the automatic turbidity correction.
	In the measured value display, measurements with automatic turbidity correction are marked by [TURB].
	For more detailed information on the automatic turbidity correction, see section 4.5.11.
Display absorbance	Activates/deactivates the display of the absorbance value in addition to the main measured value.
AQA	Here you can view and change the AQA settings without discarding the current measurement.
Edit method	Here you can edit user-defined methods.
New method	Here you can create user-defined methods.
Measurement data memory	Here you can view the measurement data memory.

### 4.5.7 Measuring diluted samples

If the concentration of a sample exceeds the measuring range of a method, you can specifically dilute the sample so that the concentration of the diluted sample is in the measuring range of the method. Thus a valid measurement is possible.

After entering the factor for the dilution the meter converts the concentration to that of the undiluted sample.



Optimum measurement results are achieved if the concentration of the diluted sample is in the middle of the measuring range of the method after diluting.

## Setting the dilution

Concentration

<HOME>

Concentration	16.01.15	9:52
Please select method for measuring or insert a barcoded cell or insert AutoSelector.		
Setup Method list Last method	New Me	thod
Concentration	16.01.15	9:52
To start measurement, insert cell or press <start enter=""></start>		
3: A6/25	Ν	IH4-N
16 mm	0.20 - 8.00	mg/l
Setup Method list Citation form	Unit	
	Ι	
Concentration	16.01.15	9:52

Inserting a cell with barcode starts a measurement.

If a cell without barcode is used: Select the method manually (see section 4.5.5).

The photometer is ready to measure.

- 1 Open the setting menu with [Setup].
- 2 Select and confirm *Dilution*. The input field for the dilution pops up.
- 3 Enter and confirm the dilution (<0...9>).

The entered dilution is taken into account with the next measurement.

The entered value for the dilution factor is valid for the selected method only. The dilution factor is erased if

NH<sub>4</sub>-N 0.20 - 8.00 mg/l

Unit

• the photometer is switched off

Method list Citation form

3: A6/25

Setup

16 m m

- a different method is selected
- the factor 0 is entered in the *Dilution* menu.

If a dilution factor is active, it is indicated on the display during measurement in the form [1 + x].

### 4.5.8 Sample blank value

By measuring and using a sample blank value, measurement errors due to coloring and turbidity of the sample matrix can be eliminated to a large extent.

The sample blank value is a characteristic of the sample (coloration) to be currently determined. It is determined by measuring the blank sample.

The sample blank value is determined with the same procedure as the corresponding analysis but without the coloring reagent. The sample blank values required are explained in detail in the relevant analysis specification.

**Validity** The sample blank value applies only to the next measurement and has thus to be remeasured for each new measured value determination.

Single and multiple<br/>determinationThe sample blank value can be determined by single or multiple<br/>determination. With multiple determination, the sample blank value is<br/>calculated as the median from the individual measured values.

Measuring the sample blank value

<HOME> Concentration

Concentration	16.01.15	9:52
Please select method for measuring or insert a barcoded cell		
or insert AutoSelector.		

Inserting a cell with barcode starts a measurement.

If a cell without barcode is used: Select the method manually (see section 4.5.5).

Concentration	16.01.15	9:52
To start measurement, insert cell or press <start enter=""></start>		
3: A6/25	Ν	IH <sub>4</sub> -N
16 mm	0.20 - 8.00	mg/l
Setup Method list Citation form	Unit	

Sample blank value	16.01.15	9:52
To start measurement, insert cell or press <start enter=""></start>		
3: A6/25 16 mm	N 0.20 - 8.00	IH <sub>4</sub> -N   mg/

Sample blank value			16.01.15 9:52	
	Last measure	d absorbance		
	0.115	0.115		
	Median			
	0.115 (1 N	1easuremen	t(s))	
3: A6/25 16 mm			NH <sub>4</sub> -N 0.20 - 8.00 mg/l	
Next meas.	Discard		Apply	

Concentration			16.01.15 9:52
			[SB]
	neasurement,   or press <st< td=""><td>ART/ENTER&gt;</td><td></td></st<>	ART/ENTER>	
3: A6/25			NH <sub>4</sub> -N
16 m m			0.20 - 8.00 mg/l
Setup	Method list	Citation form	Unit

The photometer is ready to measure.

- 1 Open the setting menu with [Setup].
- 2 Select and confirm *Sample blank value*.
- **3** Insert the cell with a suitable blank sample.

The first single measurement for the sample blank value takes place.

The following data is displayed as the result:

- The measured absorbance from the (last) single measurement.
- The median from all single measurements carried out up to now.
- 4 If necessary, carry out further single measurements for the formation of the median with [Next meas.] or discard the last single

measurement with [Discard].

5 To accept the median value, press [Apply].

The photometer is ready to measure.

The use of the sample blank value is indicated by [SB] in the top right corner of the display.

### 4.5.9 Reagent blank value

The evaluation of the photometric measurement always refers to the comparison value of a test sample without the substance to be determined (reagent blank value). Thus the influence of the basic absorbance of the reagents on photometric measurement, e.g. the inherent coloring, is taken into account.

In practice, the reagent blank value is measured with the same amount of deionized water instead of sample.

Factory and user-<br/>defined reagent<br/>blank valuesWith photometric concentration determination, the reagent blank value is a<br/>constant. The method data for all measurements with WTW test sets<br/>(Concentration mode) include an exactly determined reagent blank value.<br/>This value is overwritten if you measure the reagent blank value yourself<br/>(setting, User-defined blank value, see section 4.5.6).



You can increase accuracy if you determine the reagent blank value with a test of a new lot and use the reagent blank value for all further measurements with this lot. This is especially recommended for measurements in the vicinity of the lower limit of the measuring range. To be able to attribute the reagent blank value in the measured value documentation later, you can enter the lot number of the reagent package (*Lot number*) during the blank value determination.

Validity	The factory blank values always remain stored in the meter and can be activated at any time. The reagent blank values you measured yourself also remain stored in the meter until they are overwritten by a new blank value measurement.
Single and multiple determination	The reagent blank value can be determined with single or multiple determination. With multiple determination, the reagent blank value is calculated as the median from the individual measured values.

User-defined	For user-defined methods, you can activate the reagent blank value function
methods	as follows only:

Entry type	Function type	Reagent blank value possible?
Entry of a function	Linear	Yes
(with and without entering the ordinate intercept)	Nonlinear	No
Entry of value pairs or measurement and storage of standard solutions	Linear	Yes
	Parabola (second-order function)	Yes
(with entering/measuring and storing E0)	Polygon line	No
Entry of value pairs or	Linear	Yes
measurement and storage of standard solutions (without entering/measuring and storing E0)	Parabola (second-order function) Polygon line Polygon line through zero	No



If no value for E0 is stored during the entry of value pairs or the measurement and storing of standard solutions for a nonlinear function (parabola or polygon line), the message, *No blank value correction is intended for this method.* appears when the *User-defined blank value* function is activated. The blank value (E0) can be entered later by editing the method.

# Measuring the reagent blank value

Concentration

<HOME>

Concentration			16.01.15	9:52
Please se	lect method fo	or measuring		
	a barcoded ce			
or insert	Auto Selector.			
Setup	Method list	Last method	New Met	hod

Inserting a cell with barcode starts a measurement.

If a cell without barcode is used: Select the method manually (see section 4.5.5).

Concentration		16.01.15 9:52
		10.01.15 9:52
To start measurement insert cell or press <s< td=""><th></th><td></td></s<>		
3: A6/25		NH <sub>4</sub> -N
16 m m		0.20 - 8.00 mg/l
Setup Method list	Citation form	Unit

Concentration			16.01.15 9:52
Adjust			
Zero ad	justment		
Reagent	blank value		
Calibrat	e the method		
3: A6/25			NH <sub>4</sub> -N
16 m m			0.20 - 8.00 mg/
Setup	Method list	Citation form	Unit

Reagent blank	value		16.01.15	9:52
	neasurement, I or press <st< th=""><td>ART/ENTER&gt;</td><td></td><td></td></st<>	ART/ENTER>		
3: A6/25 16 mm				IH <sub>4</sub> -N ma/l
16 mm			0.20 - 8.00	

The photometer is ready to measure.

1 Using <**ZERO·BLANK**>, open the *Adjust* selection list.

or

Open the setting menu with [Setup].

2 Select and confirm *Reagent blank*.

The window for the measurement of the reagent blank value pops up.

The data of the last measurement appears in the measured value display.

**3** Insert the cell with the blank sample.

The first single measurement for the reagent blank value takes place.

The following data is displayed as the result:

- The measured absorbance from the (last) single measurement.
- The median from all single measurements carried out up to now.

Reagent blank value	16.01.15 9:52	4	If necessary, carry out further
Last measured absorbance 0.600 Median 0.600 (1 Measuremer	nt(s))		single measurements for the formation of the median with [Next meas.] or discard the last single measurement with [Discard].
3: A6/25 16 mm Next meas. Discard	NH <sub>4</sub> -N 0.20 - 8.00 mg/l	5	To accept the median value, press <i>[Apply]</i> .
Discalu	Арріу		The <i>Lot number</i> entry field pops up.
		6	Enter and confirm the <i>Lot number</i> ( <b><a9>)</a9></b> . The blank value measurement is completed.
Reagent blank value	16.01.15 9:52		The photometer is ready to measure.
To start measurement, insert cell or press <start enter=""></start>			The use of the reagent blank value is indicated by [BV/Lot number] in the top right corner of the display.
3: A6/25 16 mm Setup Method list Citation form	NH <sub>4</sub> -N 0.20 - 8.00 mg/l Unit		

### 4.5.10 User calibration (standard adjustment)

With some methods for concentration measurement, there is the possibility to optimize with a user calibration the original calibration stored with the method.

This makes sense, for example, if the original calibration of the method has changed due to the lot.

When creating a user-defined method you can also allow a user calibration (see section 4.5.12).

A user calibration is only valid if it deviates from the original calibration by at least 30%.

The absorbance for a user calibration can be measured as a single measurement or multiple measurement. With multiple measurement, the absorbance is calculated as the median from the individual measured values.

If a method is called up for which a user calibration is possible, a query

appears asking whether the user calibration should be utilized. If a method is called up for which a user calibration is required, measuring is only possible with a valid user calibration.

The implementation of the user calibration is documented together with the measured value and indicated in the measured value display by [Cal].

**Validity** A user calibration is always stored for the method just called up. A user calibration is only erased if

- a new user calibration is carried out
- the original calibration is selected for measurement
- the user calibration is erased manually
- the photometer is reset to the default condition.

# Performing a user calibration

# <HOME>

Concentration

Concentration			16.01.15	9:52
Please select method for measuring or insert a barcoded cell or insert AutoSelector.				
Setup	Method list	Last method	New Met	:hod

Select the method manually (see section 4.5.5).

If data for zero adjustment, reagent blank value or a user calibration are already existing, the photometer informs you of this. You can accept or discard the existing values.

If no zero adjustment is available, the photometer will prompt you to carry out a zero adjustment.

Concentration			16.01.15	9:52
Adjust				
Zero adj	ustment			
Reagent	blank value			
Calibrate	the method			
7302: Cu-1				Cu
20 m m			0.04 - 5.00	mg/l
Setup	Method list	Citation form	Unit	

Calibrate the method		16.01.15 9:52
	Target value (Cu)	Absorbance
E0	0.00 mg/L	
1	4.00 mg/l	
Ba	ck	Next

Calibrate the method		16.01.15	9:52	
To start measurement, insert cell or press <start enter=""></start>				
7302: Cu-1 20 mm			0.04 - 5.00	Cu   mg/l

1 Using <**ZERO·BLANK**>, open the *Adjust* selection list.

or

Open the setting menu with [Setup].

2 Select and confirm *Calibrate the method*.

If data for a user calibration are already existing, the list for all standard solutions includes the calibration data of the last user calibration.

If no data of a user calibration are available, the list to measure the *Absorbance* for all calibration standards required pops up.

3 In the *Target value* column, enter the nominal values of the individual standard solutions.

For E0 (reagent blank value) the nominal value is fixed and cannot be changed. The corresponding absorbance has to be measured.

4 Select an absorbance value and confirm with <**START·ENTER**>.

The measurement window pops up.

5 Insert the cell with the corresponding standard or reagent blank value (for E0).

The first single measurement for the calibration is carried out.

The following data is displayed as the result:

- The measured absorbance from the (last) single measurement.
- The median from all single measurements carried out up to now.

Calibrate the method			16.01.15 9:52
	Last measure	d absorbance	
	0.600		
	Median		
	0.600 (1 M	leasuremen	t(s))
7302: Cu-1			Cu
20 m m			0.04 - 5.00 mg/l
Next meas.	Discard		Apply

Calibra	te the method	16.01.15 9:52
	Target value (Cu)	Absorbance
E0	0.00 mg/L E	0.000
1	4.00 mg/l 1	0.600
Ba	ck	Next

Calibrate the r	Calibrate the method		16.01.15	9:52
The cali	pration has bee	n successfully co	ompleted.	
			•	
Protocol ID:	2			
Date:	16.01.2012			
User:	Admin			
Curve type:	Straight line			
Correction:	105%			
7302: Cu-1				Cu
Cancel	Calibration		Appl	,
Cance	Calibration		Apply	y

6 If necessary, carry out further single measurements for the formation of the median with [Next meas.] or discard the last single

measurement with [Discard].

7 To accept the median value, press [Apply].

The list of the standards required for this method pops up. The measured absorbance is entered for the standard or reagent blank value (E0).

8 Select all fields consecutively in the column *Absorbance* and start the corresponding measuring procedure with

# <START.ENTER>.

When <u>all values</u> have been measured (including the reagent blank value E0):

9 Accept the values with Next.

The result of the calibration pops up.

If necessary, display the list with the value pairs nominal value and absorbance with *Calibration data*.

If necessary, display the calibration line with *Graphic* in the window with the value pairs.

**10** Accept the calibration with *Apply*.

Calibrate the r	nethod	4E	16.01.15 9:52
User calibration	n:		
Protocol ID: Date:	2 16.01.2012		
User:	Admin Straight line		
Curve type: Correction:	Straight line 105%		
7302: Cu-1			Cu
End	Calibration	Delete	New

Calibrate the method 16.01.15 9:52 User calibration: Protocol I Lot number Date: User: Admin Curve type: Straight line 105% Correction 7302: Cu-1 Cu End Calibration Delete New Calibrate the method 16.01.15 9:52 [Cal][BV/Lot number][10.01.15 8:32] To start measurement, insert cell or press <START/ENTER>

To start measurement, insert cell or press <START/ENTER> 7302: Cu-1 Cu 20 mm 0.04 - 5.00 mg/l Setup Method list Unit

Viewing the data of the user calibration

<HOME> Concentration If necessary, display the list with the value pairs nominal value and absorbance with *Calibration data*.

If necessary, display the calibration line with *Graphic* in the window with the value pairs.

If necessary, erase the user calibration with *Delete*.

If necessary, carry out a new user calibration with *New measurement*.

**11** Use *End* to finish calibration.

The input field *Lot number* pops up. Here you can input the *Lot number* of the reagent blank value (E0).

12 Enter the Lot number of the reagent blank value (<A...9>) and confirm.The user calibration is completed.

The photometer is ready to measure.

If the user calibration is utilized, [Cal] appears on the display.

Note: calibration is not successful if a new value deviates by more than 30 % from the value of the old calibration.

Concentration	16.01.15	9:52
Please select method for measuring or insert a barcoded cell or insert AutoSelector.		
Setup Method list Last method	New Met	:hod

Concentration			16.01.15	9:52
Adjust				
Zero ad	justment			
Reagen	t blank value			
Calibrate	e the method			
7302: Cu-1				Cu
20 m m			0.04 - 5.00	mg/l
Setun	Method list	Citation form	Unit	

Calibrate the r	nethod	4 E	16.01.15 9:52
User calibration	<b>1</b> :		
	2		
Protocol ID:			
Date:	16.01.2012		
User:	Admin		
Curve type:	Straight line		
Correction	105%		
7302: Cu-1			Cu
End	Calibration	Delete	New

Select the method manually (see section 4.5.5).

If data for zero adjustment, reagent blank value or a user calibration are already existing, the photometer informs you of this. You can accept or discard the existing values.

1 Using <**ZERO·BLANK**>, open the *Adjust* selection list.

or

Open the setting menu with [Setup].

2 Select and confirm *Calibrate the method*.

The *Calibrate the method* window opens.

The data of the last measurement appears in the window.

If necessary, display the list with the value pairs nominal value and absorbance with *Calibration data*.

If necessary, display the calibration line with *Graphic* in the window with the value pairs.

If necessary, erase the user calibration with *Delete*.

If necessary, carry out a new user calibration with *New measurement*.

If necessary, use *End* to finish the calibration.

# Measuring with user calibration

<HOME> Concentration

Concentration	ncentration 16.01.15 9:5					
	[Cal][BV/2c][ZERO 10.01.2012 11:08]					
User calib	ration					
	A calibration dated xxx is available for this method. Should it be used?					
Yes						
No						
7302: Cu-1				Cu		
12 mm			0.04 - 5	5.00 mg/l		
Setup	Method list	Citation form		Unit		

Select the method manually (see section 4.5.5).

If data for zero adjustment, reagent blank value or a user calibration are already existing, the photometer informs you of this. You can accept or discard the existing values.

If the existing user calibration should not be utilized, a query with further options pops up:

- Use default calibration The existing user calibration is erased. The following measurements are done with the original calibration stored with the method.
- *Carry out user calibration* The existing user calibration is erased. The process for a new user calibration is started.
- *Cancel* The existing user calibration is retained. The previous query is displayed.

The photometer is ready to measure when all the required data have been confirmed or remeasured.

Concentration			16.01.15 9:52		
	[Cal][BV/2c][ZERO 10.01.2012 11:08]				
	neasurement, I or press <st< th=""><th>ART/ENTER&gt;</th><th></th></st<>	ART/ENTER>			
7302: Cu-1 12 mm			Cu 0.04 - 5.00 mg/l		
Setup	Method list	Citation form	Unit		

# 4.5.11 Automatic Turbidity correction

The *Turbidity correction* function activates the automatic recognition and compensation of the light absorption caused by turbid substances.

After being activated the function remains permanently switched on. Measured values that were measured with *Turbidity correction* are labeled with [TURB] (turbidity correction) on the display and in the documentation (printout and memory). The *Turbidity correction* function is not active in the delivery condition.



The setting for automatic turbidity correction is used with all methods where the automatic turbidity correction makes sense. The photometer automatically decides whether or not to use the function.

Switching on the turbidity correction is activated and deactivated in the setting menu of the concentration measurement (see section 4.5.6 SETTINGS FOR CONCENTRATION MODE).

### 4.5.12 Programming / modifying user-defined methods

**Overview** For *Concentration* mode, you can yourself develop and store user-defined methods under the method numbers 1001 to 1100. The photometer software supports you when creating the methods.

**Calibration data and calibration function in** photometry, the calibration function describes the dependency between the measured parameter (e.g. concentration) and the photometric measurement result (e.g. absorbance) of a sample. The knowledge of this dependency is a prerequisite for the development of a photometric method. The calibration function is usually determined by means of a series of measurements with standard solutions of known concentrations (nominal value), e.g. a 10-point calibration.



In measuring operation, the reverse calibration function is used to output the measured absorbance as a concentration value.

**Line types** The dependency between the nominal value and absorbance is often linear in a wide range as shown in the following example:




Fig. 4-1 Example of a linear calibration function after a 10-point calibration

In the case of a linear dependency, the calibration function is determined by

means of linear regression. The slope and axis intercept (E0) are the characteristics of the calibration line.

In the case of a nonlinear dependency, the points of the measuring ranges can be connected to each other as a polygon line or approximated as a parabola:






Fig. 4-2 Example of a polygon line calibration function after a 10-point calibration





Determining the calibration function

You have the following options to create a method:

#### Measure and store:

Carry out a series of measurements with the following sample solutions while at the same time the photometer takes over the values:

- Blank sample to determine the reagent blank value (with deionized water instead of sample, see section 4.5.9)
- at least one, up to ten, standard solutions in different concentrations.

The photometer stores nominal value/absorbance value pairs of the individual measurements and determines the resultant characteristics of the calibration. When doing so, you can select the following line types: *Polygon line, Straight line* or *Parabola*.

### Enter as value pairs:

Entry of the value pairs, Nominal value (concentration) / Measured absorbance of an <u>already available</u> test series with the following sample solutions:

- Blank sample to determine the reagent blank value (with deionized water instead of sample, see section 4.5.9)
- at least one, up to ten, standard solutions in different concentrations.

Based on the entered value pairs, the photometer determines the characteristics for the calibration. For this you can select the following line types: *Polygon line, Straight line* or *Parabola*.

### • Enter a function:

Entry of a function to calculate the concentration from the absorbance (reverse calibration function). You can enter on the photometer the coefficients of a polynomial equation of the following type:

 $c = a0 + a1 \cdot A + a2 \cdot A^2 + a3 \cdot A^3 + a4 \cdot A^4 + a5 \cdot A^5$ 

with:

С	Measurement result, e.g. concentration
a0 to a5	Coefficients (input range 0.000 to 1000.000)
А	Absorbance



Entering the formula is especially simple if you measure with a commercial test set for which the manufacturer has given the value for the coefficients a1. It is often called the "Factor" and corresponds to the reciprocal value of the slope of the straight line of the calibration function.

If a <u>linear</u> function (straight line) should be entered, it is necessary to enter the coefficients a0 and a1 to receive correct measured values. If the exact value for a0 is not known at the time the formula is entered, it is sufficient to enter the coefficient a1. In this case, the *User-defined blank value* function (in the *Concentration / Setup* menu) has to be activated to measure with this method.

Prior to measuring with this method, a blank value measurement has to be carried out. This procedure determines the value for a0, which then replaces the value from the programming of the method.

If the *User-defined blank value* function is not activated, the photometer uses the value zero for the coefficient a0.

More information on the entry of the formula (determination of coefficients)

Linear function	If the value for a1 (slope of the reverse calibration function) is unknown, you can very simply program the method in the photometer by measuring/storing or entering the value pairs (see above).
	For entry as a formula, you can determine the coefficients of the reverse calibration function by linear regression. The concentration has to be on the Y axis and the absorbance on the X axis.
	In the case of a linear function, the coefficients of the reverse calibration function can also be determined from the determined reagent blank value and the slope (m) of the calibration function (Y axis = absorbance, X axis = concentration). Proceed as described below.
	Explanation of the coefficients of the formula:
	<ul> <li>a0 = - E0*a1         [E0 = reagent blank value (absorbance at concentration 0)]     </li> </ul>
	<ul> <li>a1 = 1/m Reverse value of the slope of the calibration function (often referred to as "Factor") m = slope of the calibration function</li> </ul>
	<ul> <li>a2, a3, a4, a5 = further coefficients (when entering a linear function: zero)</li> </ul>
Nonlinear function	The coefficients of the reverse calibration function are determined by multiple regression. When doing so, the concentration has to be on the Y axis and the absorbance on the X axis.

## Further method data

Input field	Possible entries
Number *	1001 1100
Designation	Any name (max. 18 characters)
Version	Any version designation (max. 18 characters)
Wavelength *	Freely selectable (in nm)
Cell *	16 (round), 10, 20 or 50 mm
Citation form	e.g. PO4-P (max. 18 characters)
Unit **	e.g. mg/l (max. 18 characters)
Resolution *	0.001, 0.01, 0.1 or 1
Lower and upper limit of the measuring range *	Any value between zero and the highest concentration of the used standard solutions
Timer 0 to 3	Up to four analysis timers freely adjustable
AQA2 target value	Any value within the measuring range
AQA2 tolerance	Any
Required	1 or greater
measurements	Number of measurements after which a measured value is documented. With more than one measurement, the documented value is the median from all measurements.
Reagent blank value required	Yes/No
User calibration possible	Yes/No
User calibration required	Yes/No

\* necessary inputs

\*\* default: mg/l



If a nonlinear calibration line is programmed for a method, it may be the case that the default setting of the following menu items cannot be changed:

- Reagent blank value required
- User calibration possible
- User calibration required

How to program user-defined methods

<home> Concentration</home>	
Edit method	16.01.15 9:52
	10.01.15 9.52

Number			1001
Designation			Nitrite
Version		-	01
Wavelength			525
Cell			10 mm
Citation form			NO2-N
Unit			mg/l
Resolution			0.001
Calibration curv	e	Measure sta	andard solutions
Method list	Delete		Next

Enter the general method data 1 here. The next available method number is already entered as the number.


You have the following options when filling out the input fields:

- Fill out all empty input fields one after the other
- Using [Method list], select an already existing method as a model, give it a new method number and adjust the entries
- Using [Method list], select an existing method in order to change it (without changing the number).
- You can delete the method completely with [Delete].
- 2 Select the menu item, *Calibration* curve. Select the method for the determination of the calibration line. The following variants can be selected:
  - Measure standard solutions
  - Enter value pairs
  - Enter formula
- **3** Using *[Next]*, accept all entries on the page and switch to the next page.



During the following proceeding, you can return to the previous page at any time with [Back], e. g. if you want to correct entries, add further value pairs or eliminate outliers.

### Variant 1: Measure standard solutions



Edit method			16.01.15 9:52
	Target value		Absorbance
E0	0.000		
1			
Back	Add	Delete	Next

Edit meth	od		16.01.15 9:52
	Target value		Absorbance
E0	0.000		
1	0.300		
2	0.600		
3	1.000		
Back	Add	Delete	Next

Edit meth	nod	16.01.15 9:52
	Target value	Absorbance
E0	0.000	
1	0.300	
2	0.600	
3	1.000	
Back	Add Delet	e Next

- 1 Select and confirm *Measure standard solutions*.
- 2 Enter and confirm details of the standard solutions (optional).
- **3** Using *[Next]*, accept all entries on the page and switch to the next page.

The table for the measurement of standard solutions pops up.

In the first two lines of the table, the two value pairs (measuring points) that are at least required for a calibration are already prepared (reagent blank value E0 and any further nominal value).

4 Create further values pairs with [Add] as necessary.

You can delete a highlighted value pair with [Delete].

5 In the *Target value* column, enter the nominal values of the individual standard solutions.

Measuring the standard solutions:

6 Using the arrow keys <▲><▼> and <◀><▶>, navigate to the relevant input field in the *Absorbance* column and press <START-ENTER>.

Absorbance E0	16.01.15 9:52		The measurement display
			appears.
To start measurement,			Insert the cell with the respective standard.
insert cell or press <start enter<="" td=""><td>,</td><td></td><td>The absorbance is measured. The result of the first single measurement is displayed.</td></start>	,		The absorbance is measured. The result of the first single measurement is displayed.
525 nm	16 mm		
Absorbance E0	16.01.15 9:52	8	If necessary, carry out further
Last measured absorbance 0.009	2		single measurements for the formation of the median with [Next meas.]
Median			or
0.009 (1 Measurem	ent(s))		discard the last single measurement with [Discard].
525 nm	16 mm	9	To accept the median value, press
Next meas. Discard	Apply		[Apply].




If the zero standard concentration (reagent blank value E0) is not measured and stored, the photometer calculates the calibration line without this value. If the function *User-defined blank value* (in the *Concentration / Setup* menu) is activated for measuring with this method, the value for a0 is determined and replaces the calculated axis intercept from the programming of the method (a0 see page 62).

Edit methe	bd		16.01.15 9:52
	Target value		Absorbance
E0	0.000		0.009
1	0.300		0.664
2	0.600		1.292
3	1.000		2.178
Back	Add	Delete	Next

- **10** Repeat the steps 6 to 9 until all input fields in the *Absorbance* column are filled out.
- **11** Using *[Next]*, accept all entries on the page and switch to the next page.

The value pairs are displayed in a diagram (standard: Polygon line).



Edit method

Timer 0

Timer 1

Timer 2

Timer 3 AQA2 target value

AQA2 tolerance

Back

Required measurements

User calibration required

Reagent blank value required User calibration possible The related formula f(x) and correlation coefficient  $R^2$  are displayed above the diagram.

- **12** If required, select a different line type for the line adjustment with *[Curve type]*.
  - Polygon line
  - Straight line
  - Parabola
- **13** If required, enter different measured value limits with [Meas. range].
  - Lower limit
  - Upper limit
- **14** Using *[Next]*, complete the editing of the calibration line and proceed to the next page.

The timers and AQA2 data linked to the method are displayed.

- **15** If necessary, enter intervals for up to 4 timers.
- **16** If necessary, enter the *AQA2 target value* and *AQA2 tolerance*.
- 17 If necessary, set the number from how many single measurements the documented measured value will be calculated.
- **18** If necessary, specify whether a reagent blank value is required.
- **19** If necessary, set whether a user calibration is possible and/or required.
- 20 Complete the programming of the method with [Complete].

The method is programmed and selected for measuring.

### Variant 2: Enter value pairs

Unlike variant 1, the fields of the *Absorbance* column are filled out manually here. Accordingly, the steps 6 to 10 are not applicable here. Apart from that, the proceeding is identical to variant 1.

16.01.15 9:52

00:00:00

00:00:00

00:00:00 00:00:00

1.00 mg/l

0.10 mg/l

Complete

No

No

No

Variant 3: Enter formula

Edit method			16.01.15 9:52	
$c = a0 + a1 \cdot A + a2 \cdot A^2 + a3 \cdot A^3 + a4 \cdot A^4 + a5 \cdot A^5$				
a0			0.605	
a1			2	
a2				
a3				
a4				
a5				
Lower limit of m	neasuring rang	e	1,000 mg/l	
Upper limit of measuring range			3.000 mg/l	
Method list	Delete		Next	

- Select and confirm *Enter formula*.
   Input fields for the coefficients (a0 ... a5) of the formula are displayed.
- 2 Enter and confirm the factors.

If no value is entered for a coefficient the photometer automatically uses the value 0.



Entering the formula is especially simple if you measure with a commercial test set for which the manufacturer has given the value for the coefficients a1. It is often called the "Factor" and corresponds to the reciprocal value of the slope of the straight line of the calibration function.

If a linear function (straight line) should be entered, it is necessary to enter the coefficients a0 and a1 to receive correct measured values. If the exact value for a0 is not known at the time the formula is entered, it is sufficient to enter the coefficient a1. In this case, the *User-defined blank value* function (in the *Concentration / Setup* menu) has to be activated to measure with this method. Prior to measuring with this method, a blank value measurement has to be carried out. During this procedure the value for a0 is determined and replaces the previous value.

- **3** Enter and confirm the measuring range limits.
- 4 Complete the entering of the formula with [Next].

The timers and AQA2 data linked to the method are displayed.

Edit method	16.01.15 9:52		
Ti	00.00.00		
Timer 0	00:00:00		
Timer 1	00:00:00		
mer 2 00:00			
Timer 3	00:00:00		
AQA2 target value	2.000 mg/l		
AQA2 tolerance	0.200 mg/l		
Required measurements	1		
Reagent blank value required	No		
User calibration possible	No		
User calibration required	No		
Back	Complete		

- 5 If necessary, enter intervals for up to 4 timers.
- 6 If necessary, enter the AQA2 target value and AQA2 tolerance.
- 7 If necessary, set the number from how many single measurements the documented measured value will be created.
- 8 If necessary, specify whether a reagent blank value is required.
- **9** If necessary, set whether a user calibration is possible and/or required.
- **10** Complete the programming of the method with [Complete].

The method is programmed and selected for measuring.

### 4.5.13 The *IQ-LabLink* procedure

The IQ-LabLink procedure enables the secure data exchange between the sensors of the IQ SENSOR NET online measuring system and a photometer (such as the photoLab<sup>®</sup> 6000 series) with the aid of a commercial USB memory device.

During the matrix adjustment of the sensor, the measurement data of the IQ SENSOR NET sensor are adjusted to the photometrically determined reference data.

Previously, the reference data had to be entered manually during the matrix adjustment. With the aid of the IQ-LabLink procedure, the reference data can now be directly transferred to the IQ SENSOR NET sensors by means of a USB memory device, without manual input and with no risk of confusion.

- IQ SENSOR NET:
  - Terminal/controller with USB-A interface and software for the IQ-LabLink procedure (e.g. MIQ/TC 2020 XT)
  - Online sensor with software for the IQ-LabLink procedure (e.g. VARiON<sup>®Plus</sup> 700 IQ)
- Photometer:
  - Photometer with software for the IQ-LabLink procedure (e.g. photoLab<sup>®</sup> 7xxx)

System requirements for the IQ-LabLink procedure Course of the IQ-LabLink procedure

### Step 1 on the IQ SENSOR NET terminal:

Automatic creation of a job file on the USB memory device with current sensor values, parameters, designation of the measuring location and automatic allocation of a job reference number for clear identification.

### Step 2 on the photometer:

Automatic recognition of the job files, menu-guided measurement of all required parameters, storage of the determined data in the job file.

### Step 3 on the IQ SENSOR NET terminal:

Automatic recognition of the job files, complete reading of all data required for the matrix adjustment on keypressing.



To carry out the matrix adjustment of an online sensor with the IQ-LabLink procedure requires operating steps on both instruments: IQ SENSOR NET system and photometer of the photoLab<sup>®</sup> 6000 series.

The detailed description of the cross instrument operating steps for the matrix adjustment with the IQ-LabLink procedure on the IQ SENSOR NET and the photometer is given in an additional operating manual. This operating manual can be downloaded from the Internet at www.wtw.com.

### 4.6 Measuring the Absorbance / % Transmission

### 4.6.1 General information

The absorbance or transmission respectively is measured without the use of any methods or profiles. All settings are configured during measurement.

Measuring against the Reference absorbance The absorbance or transmission can alternatively be measured against the absorbance of the zero adjustment or against a *Reference absorbance* determined by yourself (see section 4.6.3 MEASURING AGAINST THE REFERENCE ABSORBANCE).

### 4.6.2 Measuring the absorbance or transmission


### <HOME> Absorbance / % Transmission

 Absorbance
 16.01.15
 9:52

 To start measurement, insert cell or press <START/ENTER >

 300 nm

 Setup
 Wavelength
 Transmission
 Reference




The settings of the last measurement are active.

- 1 Using [Wavelength], change the wavelength as necessary.
- 2 Using [Absorbance] <-> [Transmission], you can switch over between absorbance and transmission measurement.
- 3 If necessary, use or measure a reference measurement with [*Reference*] (see section 4.6.3).
- 4 Depending on the type, insert the cell as follows:

### Round cell:

Insert the round cell in the round cell shaft so that it touches the bottom.

If the inner turn-up lid is opened too wide, a message prompts you to close the inner turn-up lid.



Absorbance		5E	16.01.15	9:52	
0.860					
489 nm			1	0 mm	
Setun	Wavelength	Transmission	Refere	nce	



Rectangular cell:

Open the inner turn-up lid.

Insert the rectangular cell vertically so that it touches the bottom and left edge of the cell shaft. The opaque sides of the rectangular cell must point to the front and back.

The photometer has an external light recognition. If there is too much external light, a message prompts you to close the cell shaft cover.

The photometer starts measuring automatically.

5 Using [Absorbance] <-> [Transmission], switch over the display from Absorbance to Transmission or vice versa.

# 4.6.3 Measuring against the Reference absorbance

Each time the photometer is switched on, the absorbance or transmission is measured against the absorbance of the zero adjustment as a basis. You can, however, also determine a *Reference absorbance* and use it as the basis.

The *Reference absorbance* refers to the adjusted wavelength. The measured value remains stored until

• the photometer is switched off

- the cell type is changed
- the wavelength is changed
- a new reference value is measured
- it is deleted manually ([Reference] / Delete).
- the Absorbance / % Transmission measuring mode is exited

# Single and multiple determination

The Reference absorbance can be determined with single or multiple determination. With multiple determination, the mean value is calculated as the median from the individual measured values.

#### Reference absorbance Measuring the

Absorbance / % Transmission

<HOME>

Absorbance	16.01.15 9:52
To start measurement, insert cell or press <start enter=""></start>	
489 nm	10 mm
Setup Wavelength Transmission	Reference

Reference absorbance			16.01.15	9:52
	neasurement, I or press <st< td=""><td>ART/ENTER&gt;</td><td></td><td></td></st<>	ART/ENTER>		
489 nm			1	0 mm
			1	

The settings of the last measurement are active.

1 Start the reference measurement with [Reference].

If a value for the reference absorbance is already stored, it can be deleted or overwritten by a new reference measurement.

After the reference absorbance value has been deleted, the photometer measures against the absorbance of the zero adjustment.

2 Insert the cell with the reference sample.

The first single measurement for the Reference absorbance is carried out.

The following data is displayed as the result:

- The measured absorbance from the (last) single measurement.
- The median from all single measurements carried out up to now.

Reference absorbance 16.01.15 9:52		16.01.15 9:52
	Last measured absorbance	
	0.232	
	Median	
	0.232 (1 Measuremer	nt(s))
489 nm		10 mm
Next meas.	Discard	Apply
Absorbance		16.01.15 9:52
		Reference
	t measurement, cell or press <start enter=""></start>	
mserre	Lei of press < 3TART/LIVTER>	
489 nm		10 mm
Setup	Wavelength Transmission	Reference
- a ca p	, and the second second	

- 3 If necessary, carry out further single measurements for the formation of the median with [Next meas.] or discard the last single measurement with [Discard].
- 4 To accept the median value, press [Apply].

The photometer is ready to measure.

The reference absorbance is displayed in the top right corner during absorbance or transmission measurement.

#### 4.7 Special / Multi wavelengths methods

### 4.7.1 Basic information on Special / Multi wavelengths measurements

In the Special / Multi wavelengths mode of the photoLab<sup>®</sup> 7600 UV-VIS, you can carry out measurements with special methods and functions.

You can use the following functions for these methods:

- Measurements at different wavelengths
- Multiple measurements at one wavelength (e.g. before and after adding a reagent)
- Use of procedure variables.
   Procedure variables provide a value that has to be entered prior to each measurement on the photometer (e.g. volume, pH value or temperature)
- Check whether a value meets a condition. With a condition you can check a value for validity (e.g. absorbance value, procedure variable or the result of a formula).
- Formula editor for the convenient programming of any user-defined methods

**Special methods** The method list in the Special / Multi wavelengths mode comprises:

- preprogrammed multi wavelengths methods
- preprogrammed special methods
- special methods programmed by the user



If you program any special methods yourself, you can use all extended functions of the Special / Multi wavelengths mode.

### 4.7.2 Programming / modifying the Special / Multi wavelengths methods



For multi wavelength methods, you can use the method numbers 2001 to 2499. All special methods can also be selected in the method list of the concentration mode.

The creation of a user-defined method is done in the following steps:

- Enter the general method data Method number, method name, unit etc.
- Enter the wavelengths for absorbance measurements (A<sub>x nm</sub>)
   Minimum 1, maximum 10
- Define the procedure variables (K<sub>x</sub>) (optional)
   Procedure variables are used to take into account any influence quantities that cannot be measured by the photometer.

The values for these procedure variables have to be entered for all measurements with the method, e.g. the temperature or pH value.

• Enter the formula to calculate the measurement result Enter the formula with which you want to calculate the measurement result in the formula editor.

#### • Enter an additional condition (optional)

Conditions are used to check the measurement result for validity.

The condition is entered with the formula editor.

Example: Determination of chlorophyll a according to Nusch The chlorophyll determination is based on two measurements (before and after adding an acid) of the optical density (= absorbance) of the extract of an aqueous sample at 665 nm.

Chlorophyll a ( $\mu$ g/l) = 29.6 \* ( $A_{(before) 665 nm} - A_{(after) 665 nm}$ )\*( $V_{Extract}/V_{Sample}$ )

with.	

${oldsymbol{\mathcal{A}}}_{( ext{before})}$ 665 nm	1st absorbance measurement at 665 nm (before adding the acid)
$oldsymbol{A}_{(after) \ 665 \ nm}$	2nd absorbance measurement at 665 nm (after adding the acid)
$V_{\scriptscriptstyle Extract}$	Volume of the extract (in ml)
$V_{\scriptscriptstyle {Sample}}$	Volume of the aqueous sample (in ml)

**Converted equation** For entry on the photometer, assign names that you can enter in the formula editor on the photometer to the variables of the equation.

 $R = 29.6 * (A_{665nm} - A_{665nm_2})*(K_1/K_2)$ 

R = result (concentration chlorophyll A in $\mu$ g/l)
Variables for absorbance. These values are measured by the photometer. Here: Two measurements at the same wavelength, at different points in time.
With several measurements (e.g. before and after adding the acid), the names of the variables differ in the index _y added after an underscore (e.g. $A_x$ mm, $A_x$ mm_2, $A_x$ mm_3, etc.).
Procedure variables
K1 = Volume of the extract (in ml)
K2 = <i>Volume of the aqueous sample (in I)</i> Freely selectable numerical values

Edit method	16.01.15 9:52
Number	2001
Name	Chlorophyll a
Version	1.0
Citation form	Chl a
Unit	µg/l
Resolution	0.1
Cell	10 mm
Lower limit of measuring	0 µg/l
Upper limit of measuring	1000 µg/l
Method list	Delete Next

Wavelength 16.01.15 9:52 665 nm Wavelength 1 Delete Back Add Next 16 01 15 9:52 Procedure variables Procedure variables are variables whose current numerical values have to be entered during the course of the measurement (e.g. weighted sample or dilution). If a procedure variable is required to calculate the result: Create a procedure variable (K) with <Add> Back Next

1 Enter the general method data here. The next available method number is already entered as the number.

You have the following options when filling out the input fields:

- Fill out all empty input fields one after the other
- Using [Method list], select an already existing method as a model, give it a new method number and adjust the entries
- Using *[Method list]*, select an existing method in order to change it (without changing the number).
- You can delete the method completely with [Delete].
- 2 Using [Next], accept all entries on the page and switch to the next page.

Enter the wavelengths for the absorbance measurements  $(A_{x \text{ nm}})$ .

**3** Add another wavelength with *[Add]*.

Delete a highlighted wavelength with [Delete].

4 Using *[Next]*, accept all entries on the page and switch to the next page.

Create all required procedure variables.

**5** Create a procedure variable required for the formula with [*Add*] and enter a designation, e.g. the measured parameter.

or

Using *[Next]*, accept all entries on the page and switch to the next page.

Procedure variat	les		16.01.15 9:52
K 1 K 2			V (extract) V (sample)
Back	Add	Delete	Next
<u>.                                    </u>			
Formula entry			16.01.15 9:52
Use the <operators> softkey to select an operation, function or constant (e.g.: +, -, *, tan, log, e, Pi).</operators>			
Use the <variables> softkey to select an absorbance at a certain wavelength or a procedure variable.</variables>			
Enter numerals via the keyboard.			
You can erase the last entry with <◀>.			

Operators Variables

Next

Back

 Formula entry
 16.01.15
 9:52

 Variables
 4(665 nm)
 4(665 nm)

 K1 (V extract (ml))
 K2 (V sample

 Back
 Operators
 Variables

 Next
 Back

6 Add another procedure variable with *[Add]*.

or

Delete a highlighted procedure variable with [Delete].

7 Using *[Next]*, accept all entries on the page and switch to the next page.

Enter the formula.

8 Enter any numbers with <0...9>.

Use [Operators], <**▲**><**▼**> <**◀**><**▶**> and <**START**·**ENTER**> to enter an operator, a function or a constant.


Use [*Variables*], <**▲**><**▼**> <**◀**><**▶**> and <**START**·**ENTER**> to select a variable.

The formula is displayed after each step.

Using <<>> you can delete the last element of the formula.

Use [Back] to quit the formula editor.

9 Use [Variables], <▲><♥>
 <◀><▶> and <START·ENTER>
 to select and confirm a variable.
 The current version of the formula
 is displayed.



Formula entry	16.01.15 9:52
R = 29.6 * (A665nm - A665nm_2)	
Back Operators Variables	Next
Condition	16.01.15 9:52
Here you can enter a formula for a condition	The measured
value is only valid if this condition is met.	. The measureu
value is only value if this condition is met.	
Back Operators Variables	Next

**14** Using *[Next]*, accept all entries on the page and switch to the next page.

If an error is in the formula, an error message appears. The formula editor is only exited once the error is eliminated.

If necessary, enter the formula for a condition.

15 Enter any numbers with <0...9>.

Use [Operators], <▲><▼> <◀><►> and <START·ENTER> to enter an operator, a function or a constant.

Use [Variables], <▲><♥> <◀><▶> and <**START**·ENTER> to select a variable.

The condition is displayed after each step.

Using << > you can delete the last element of the condition.

Use [Back] to quit the formula editor.

- **16** Complete the condition.
- 17 Complete the programming of the method with [Next].

Condition	16.01.15 9:52
A <sub>665 nm</sub> < 2	
b5	
Back	Next

Edit method Sequence Measure ment Measure ment	16.01.15 9:52 Designation	If the formula includes several measurements at the same wavelength (measurement sequence), you can assign names to the individual measurements of the sequence.
		<b>18</b> Enter the names for the individual measurements of a sequence.
Back	Next	
Edit method	16.01.15 9:52	19 Complete the programming of the
Sequence	Designation	method with [Next].
Measurement	before	
Back	Next	
Special / Multi wavelengths	16.01.15 9:52	The method is programmed and selected.
V extract (ml)		The photometer is ready to measure.
Press <start enter=""> to ente</start>	er the value	
2001:Chl a 10 mm	Chlorophyll a	
Setup Method list Citation	n form Unit	

### 4.7.3 Selecting a Special / Multi wavelengths method

To select a method for Special / Multi wavelengths measurements, proceed as follows:

<HOME> Special / Multi wavelengths [Method list]

Selec	16.01.15 9:52			
2001	Chl a	Chlorophyll a	µg/I	
2002 2503	Chl a NO 3 UV	Chlorophyll a NO <sub>3</sub> -N	mg/ mg/	
Las	st used			

The list of methods is displayed. The methods are ordered by method numbers.

#### Select the method:

- Select the required method with
   <▲><▼>. The active selection is displayed in reverse video.
- 2 Accept the selection with **START·ENTER**>.

The photometer is ready to measure.

### Narrowing down the method list

If the list is very long, you can narrow down the method list and thus make the search easier as follows:

- Using [Last used], you can restrict the method list to the ten methods last used.
- With the search function you can search certain character strings in the list. The search takes place as a full-text search of the entire list contents. Thus you can search for a method number or certain citation form.

Search function	Select method	(last used)	16.01.15 9:52
	Chl_		10.01.13 5.32
	2001 Chla	Chlorophyll a	µlg/l
	2002 Chla	Chlorophyll a	mg/m3
	All methods		

#### Search for a character string:

Enter the character string to be searched for in the search window with **<A...9>**.

The list appearing below shows all hits containing the character string. The hit list is updated with each character that is entered.



Note the case sensitivity when searching.

#### 4.7.4 Carrying out Special / Multi wavelengths measurements

### <HOME>

Special / Multi wavelengths

Special / Multi wavelengths	16.01.15 9:52	1
Please select method for measu	iring!	
Setup Method list Citation	form Unit	
Special / Multi wavelengths	16.01.15 9:52	
V extract (ml)		
Press < START/ENTER> to enter	r the value	•
,		2
2001 CH	Chlanashulla	
2001:Chl a 10 mm	Chlorophyll a 0.00 - 1000.00 µg/l	
Setup Method list Citation	form	
Setup Method list Citation	form	
Setup Method list Citation Special / Multi wavelengths	form	3
Special / Multi wavelengths		3
		3
Special / Multi wavelengths		3 4
Special / Multi wavelengths		3 4
Special / Multi wavelengths V extract (ml) 10 ml V sample (ml) 100 ml		3 4
Special / Multi wavelengths V extract (ml) 10 ml V sample (ml) 100 ml		3 4
Special / Multi wavelengths V extract (ml) 10 ml V sample (ml) 100 ml Proceed with <start enter=""></start>	16.01.15 9:52	3 4
Special / Multi wavelengths V extract (ml) 10 ml V sample (ml) 100 ml		3 4 5
Special / Multi wavelengths V extract (ml) 10 ml V sample (ml) 100 ml Proceed with <start enter=""> 2001:Chl a</start>	16.01.15 9:52	3 4 5
Special / Multi wavelengths V extract (ml) 10 ml V sample (ml) 100 ml Proceed with <start enter=""> 2001:Chl a 10 mm</start>	16.01.15 9:52 Chlorophyll a	3 4 5
Special / Multi wavelengths V extract (ml) 10 ml V sample (ml) 100 ml Proceed with <start enter=""> 2001:Chl a 10 mm</start>	16.01.15 9:52 Chlorophyll a	3 4 5
Special / Multi wavelengths V extract (ml) 10 ml V sample (ml) 100 ml Proceed with <start enter=""> 2001:Chl a 10 mm</start>	16.01.15 9:52 Chlorophyll a	3 4 5
Special / Multi wavelengths V extract (ml) 10 ml V sample (ml) 100 ml Proceed with <start enter=""> 2001:Chl a 10 mm</start>	16.01.15 9:52 Chlorophyll a	3 4 5
Special / Multi wavelengths         V extract (ml)       10 ml         V sample (ml)       100 ml         Proceed with <start enter="">         2001:Chl a         10 mm         Setup         Repeat</start>	Chlorophyll a	3 4 5
Special / Multi wavelengths V extract (ml) 10 ml V sample (ml) 100 ml Proceed with <start enter=""> 2001:Chl a 10 mm Setup Repeat</start>	Chlorophyll a	3 4 5
Special / Multi wavelengths         V extract (ml)       10 ml         V sample (ml)       100 ml         Proceed with <start enter="">         2001:Chl a         10 mm         Setup         Repeat</start>	Chlorophyll a	3 4 5
Special / Multi wavelengths         V extract (ml)       10 ml         V sample (ml)       100 ml         Proceed with <start enter="">         2001:Chl a         10 mm         Setup         Repeat</start>	Chlorophyll a	3 4 5

Chlorophyll a 0.00 - 1000.00 µg/l

Unit

2001:Chl a 10 mm

Setup Method list Citation form

I Select the required method with [Method list] (see section 4.7.3).

For the description of the measuring procedure, the selfprogrammed method "Chl a" is selected as an example.

For methods with procedure variables: Enter the values of all procedure variables one after the other.

- 2 Go on to the next step with **START·ENTER>**.
- **3** Follow the instructions on the display.
- 4 Enter the volumes of the sample and extract.


If necessary, repeat the last step with [Repeat].


5 Go on to the next step with **<START·ENTER>**.

The photometer is ready to measure.

If necessary, carry out a zero measurement.

Special / Multi	16.01.15 9:52			
Measurement 1				
To start measurement, insert cell or press <start enter=""></start>				
2001:Chl a Chlorophyll a 10 mm				
Setup	Method list	Citation form	Unit	





The photometer is ready to measure.

6 Go on to the next step with <**START·ENTER**>.

7 Depending on the type, insert the cell as follows:

#### Round cell:

Insert the round cell in the round cell shaft so that it touches the bottom.

If the inner turn-up lid is opened too wide, a message prompts you to close the inner turn-up lid.

#### Rectangular cell:

Open the inner turn-up lid.

Insert the rectangular cell vertically so that it touches the bottom and left edge of the cell shaft. The opaque sides of the rectangular cell must point to the front and back.

The photometer has an external light recognition. If there is too much external light, a message prompts you to close the cell shaft cover.

8 Close the cell shaft cover.

Special / Multi wavelengths	16.01.15 9:52	9	Start the measurement.
Measurement 1			
To start measurement,			
insert cell or press <start enter=""></start>			
2001:Chl a 10 mm	Chlorophyll a		
Setup	Cancel		
Special / Multi wavelengths	16.01.15 9:52		An intermediary result is displayed
V extract (ml) 10 ml V sample (ml) 100 ml Measurement 1 A(665 n) = 0.600			if there are several measurements.
Proceed with <start enter=""></start>			The photometer is ready for the next measurement.
2001:Chl a	Chlorophyll a	10	Go on to the next step with <b><start·enter></start·enter></b> .
10 mm	enieropiiji a		The photometer is ready to
Setup Repeat	Cancel		measure.
Special / Multi wavelengths	16.01.15 9:52	11	Start the measurement.
Maaaning and D	1		
Measurement 2			
To start measurement,			
insert cell or press <start enter=""></start>			
2001:Chl a 10 mm	Chlorophyll a		
Setup Repeat	Cancel		
Special / Multi wavelengths	16.01.15 9:52		The result is displayed.
V extract (ml) 10 ml V sample (ml) 100 ml			If an entered condition is not met,
Measurement 1 $A(665 n) = 0.600$ Measurement 2 $A(665 n) = 0.000$			no measured value is displayed.
	ma (ml	12	If necessary, start a new measurement with the method.
1.78	mg/ml		measurement with the motiou.
Start new analysis with < START/ENTER>			
Setup	Cancel		

### 4.8 OptRF measurement

#### 4.8.1 General information

The OptRF measurement menu contains WTW methods for the optical reagent free measurement of some parameters in the UV range. The results - like those of WTW online measurements - are calculated by measuring and evaluating over a UV range between 200 and 390 nm. Therefore, a quartz cell is required for this measurement.

At the moment, the parameters are COD total, COD dissolved, NO<sub>3</sub> and NO<sub>2</sub> in the outlet of municipal waste water treatment plants, in standard operation without malfunctions. To measure, no photometric test set and no digestion is required. To measure COD dissolved, no photometric test set and no digestion <u>or filtering</u> of the sample is required.



<u>COD dissolved</u>: Always carry out the OptRF measurement with the unfiltered sample, even if the sample is filtered for the reference measurement.

To achieve the highest accuracy possible and to exclude any interferences caused by the sample matrix, you should always carry out a user calibration to optimize the system for the relevant municipal waste water treatment plant.

To do so, the measured value is adjusted with reference values by means of suitable reference methods. This is done before the first usage and then at regular intervals or as necessary (user calibration, see section 4.8.5).



Basically, the methods can only be used for the specified application because the sample matrix and effects of the turbidity affect he measurement. For more detailed information, please refer to the application instructions for measurement in the analysis specifications.

4.8.2 Selecting OptRF measurement

Selecting the method last used

The method last used is immediately selected.

# Selecting a method from the *Method list*

3001	CODt_H_Outlet_10	COD	0.0 - 75.0 mg/l
3002	CODs_H_Outlet_10	COD	0.0 - 75.0 mg/l
3003	NO3_H_Outlet_10	NO 3-N	0.0 - 3.0 mg/l
3004	NO2_H_Outlet_10	NO <sub>2</sub> -N	0.00 - 4.00 mg/

The list of methods is displayed. The methods are ordered by method numbers.

The method last selected is highlighted.

Select the method:

- Select the required method with
   <▲><▼>. The active selection is displayed in reverse video.
- 2 Accept the selection with **START-ENTER**>.

# Narrowing down the method list

You can narrow down the method list and thus make the search easier:

- Using *[Last used]*, you can restrict the method list to the ten methods last used.
- With the search function you can search certain character strings in the list. The search takes place as a full-text search of the entire list contents. Thus you can search for a method number or certain citation form.

#### Search function

Select method (last used)	16.01.15 9:5	
NO3		
3003 NO3_H_Outlet_10	NO3-N	0.0 - 3.0 mg/
3002 CODs_H_Outlet_10	COD	0.0 - 75.0 mg
All methods		

#### Search for a character string:

Enter the character string to be searched for in the search window with **<A...9>**.

The list appearing below shows all hits containing the character string. The hit list is updated with each character that is entered.



Note the case sensitivity when searching. It is not required or possible to enter inferior characters. When searching for chemical formulas, inferior characters are treated as normal characters. Example: The search for "NH4" shows all hits that contain "NH4" as well as "NH<sub>4</sub>".

4.8.3	Carrying	out measurements
-------	----------	------------------

<HOME> OptRF measurement

OptRF measurement	16.01.15	9:52
Please select method for measuring!		
Setup Method list Citation form		

 OptRF measurement
 16.01.15
 9:52

 Zero measurement required!
 Press <ZERO/BLANK>

 3003:NO3\_H\_Outlet
 NO3-N

 10 mm
 0.0 - 3.0 mg/l


OptRF measurement			16.01.15 9:52		
Adjust	Adjust				
Zero adjustment					
User cal	User calibration				
3003:NO3_H_Outlet NO: 10 mm 0.0 - 3.0 m					
Setup	Method list	Citation form			

1 Select the required method with *[Method list]* (see section 4.7.3).

- 2 Using <**ZERO·BLANK**>, open the *Adjust* selection list.
- 3 Carry out a zero adjustment.
- 4 If necessary, perform a user calibration.

OptRF measurement [ZERC	16.01.15 9:52 16.01.15 9:51]	The photometer is ready to measure.
		5 Further options:
Start measurement with $<$ START/ENT	TER>	<ul> <li>Select a different citation form with [Citation form] (e.g. NO<sub>3</sub> &lt;-&gt; NO<sub>3</sub>-N).</li> </ul>
3003:NO3_H_Outlet 10 mm Setup Method list Citation form	NO <sub>3</sub> -N 0.0 - 3.0 mg/l	<ul> <li>Make further settings such as dilution or user calibration with [Setup] (see section 4.5.6).</li> </ul>

Menu item	Explanation
Dilution	Here you can set the dilution prior to measuring if you want to use a diluted sample.
	In the measured value display, the dilution is indicated in the form $[1 + x]$ (parts sample + parts distilled water).
	For more information, see section 4.8.4.
User calibration	Here you can enter the value pairs (raw value [#] and corresponding reference value) for the method.
	For more information, see section 4.8.5.
Measurement data	Access to the measurement data memory
memory	For more information, see section 4.12.



OptRF measure	ment		16.01.15 9:52
		[ZERC	) 16.01.15 9:51] [MQ 1]
1	Raw value	<b>2.3</b> e: #2.0 mg/	mg/I 
3003:NO3_H_O	utlet		NO <sub>3</sub> -N
10 mm			0.0 - 3.0 mg/l
Setup	Method list	Citation form	

<u>COD dissolved</u>: Always carry out the OptRF measurement with the unfiltered sample, even if the sample is filtered for the reference measurement.

6 Insert the cell:

Rectangular cell:

Open the inner turn-up lid.

Insert the rectangular cell vertically so that it touches the bottom and left edge of the cell shaft. The opaque sides of the rectangular cell must point to the front and back.

The photometer has an external light recognition. If there is too much external light, a message prompts you to close the cell shaft cover.

The result is displayed.

If no user calibration was carried out, the measured value (raw value) and the result of the measurement are identical. After a user calibration, the measured raw value and the result adjusted according to the calibration are displayed.

Rating number for the measured value quality (MQ):

- 1 3: Excellent sufficient
- 4: Measured value invalid
- 7 Further options:
  - Make further settings such as dilution or user calibration with [Setup].
- 8 If necessary, start a new measurement with the method.

Possible errors	Remedy
Absorption too high	Dilution

Possible errors	Remedy
MQ value high (MQ = 2 or MQ = 3)	User calibration
MQ value too high (MQ = 4)	Use a matrix that is more suitable



A user calibration does not affect or change the evaluation number for the measured value quality. The user calibration just optimizes the calculated concentration values.

#### 4.8.4 Setting the dilution

If the concentration of a sample exceeds the measuring range of a method, you can specifically dilute the sample so that the concentration of the diluted sample is in the measuring range of the method. Thus a valid measurement is possible.

After entering the factor for the dilution the meter converts the concentration to that of the undiluted sample.



<HOME>

The displayed measuring range of the method changes with the entered dilution factor. The measuring range is no longer displayed after the dilutation factor was entered.

Optimum measurement results are achieved if the concentration of the diluted sample is in the middle of the original measuring range of the method (without entering a dilution factor).

#### Setting the dilution

OptRF measurement

OptRF measure	ement		16.01.15 9:52
		[ZERC	0 16.01.15 9:51]
Start mea	Start measurement with <start enter=""></start>		
3003:NO3_H_0	Dutlet		NO <sub>3</sub> -N
10 mm			0.0 - 3.0 mg/l
Setup	Method list	Citation form	

1 Select a reagent free method and carry out the zero measurement (see section 4.8.3)

The photometer is ready to measure.

OptRF measurement	16.01.15 9:52
Sample + distilled water	
1+_	
3003:NO3_H_Outlet	NO <sub>3</sub> -N
10 mm	0.0 - 3.0 mg/l
Setup Method list Citation form	
	L.

- 1 Open the setting menu with [Setup].
- 2 Select and confirm *Dilution*. The input field for the dilution pops up.
- 3 Enter and confirm the dilution (<0...9>).

The entered dilution is taken into account with the next measurement.

The entered value for the dilution factor is valid for the selected method only. The dilution factor is erased if

- the photometer is switched off
- a different method is selected
- the factor 0 is entered in the *Dilution* menu.

If a dilution factor is active, it is indicated on the display during measurement in the form [1 + x].

#### 4.8.5 User calibration

For optical reagent free methods, a user calibration should be carried out for the required parameters so the system is optimally adapted to the matrix of the relevant waste water treatment plant.

**Calibration points** and value pairs The adjustment is carried out at one or two points within the measuring range. A value pair is determined at each point. Each value pair consists of the raw value of the sample and the corresponding reference value.

**Selection of calibration points** If the required accuracy is low, the default setting of the first value pair can be adopted unchanged. The first value pair is equivalent to the zero point in this case. The second value pair should be close to the working point of the application. This method results in a single-point calibration.

If the application has to cover a wide value range with high measuring accuracy, both value pairs should be used for calibration. When doing so, the calibration value pairs should be determined and input at the limits of the application range.



The measurement results will be ever more accurate the closer the composition of the test sample corresponds to the status at the time of the calibration. If there is a fundamental change of the characteristics of the sample, a new calibration is required. The implementation of the user calibration is documented together with the measured value and indicated in the measured value display by [Cal].

- Validity A user calibration applies to the method currently called up. A user calibration is only erased if
  - a new user calibration is carried out
  - the original calibration is selected for measurement
  - the user calibration is erased manually
  - the photometer is reset to the default condition.

#### **Performing a user calibration** Select a sample for the user calibration that represents the real application conditions as well a possible. Determine the concentration of the required parameters of this sample promptly with OptRF measurements and with a suitable reference measurement.

The user calibration comprises the following steps:

- Carry out the reference measurement: Determine the concentration of the required parameter for this sample with a reference method.
- Carry out the *OptRF measurement*. Determine the concentration of the required parameter for this sample with the *OptRF measurement* method.
- Enter the value pairs for *OptRF measurement* (raw value) and corresponding reference value on the photometer.


Reference measurement

 Measure your sample with a reference method of your choice. Measure your sample in the lower and upper concentration range to be expected and note down the measured values (reference values).

Then measure the sample with a method in the menu *OptRF measurement*.

OptRF measurement

<**HOME>** OptRF measurement 2 Select the method and carry out the zero measurement (see section 4.8.3)



The photometer is ready to measure.

**3** Use the same samples as for the reference measurement and note down the raw values measured.

<u>COD dissolved</u>: Always carry out the OptRF measurement with the unfiltered sample, even if the sample is filtered for the reference measurement.

- 4 Open the setting menu with [Setup].
- 5 Select and confirm *User calibration*.

The data currently used (factory calibration or the calibration values last entered) are displayed as value pairs (Raw value and reference value) for two calibration points.

- 6 Select a value with <▲><▼> and open it for editing with <START·ENTER>
- 7 Enter a value (<0...9>) and confirm with <START.ENTER>.
- 8 Complete all four values.
- **9** Confirm the values of the user calibration with *[Apply]*.

# Entering the value pairs

OptRF measur			
	Raw value	Referer	nce value
Value pair 1	0.00 #mg/l	0.0	00 #mg/l
Value pair 2	5.00 #mg/l	6.0	00 #mg/l
Reset all			

OptRF measurement	16.01.15 9:52
[ZER0	0 16.01.15 9:51]
Start measurement with <start en<="" td=""><td>TER &gt;</td></start>	TER >
3003:NO3 H Outlet	NO 3-N
10 mm	0.0 - 3.0 mg/l
Setup Method list Citation form	Unit

The photometer is ready to measure.

#### 4.9 Spectrum

#### 4.9.1 General information

With the *Spectrum* function, the *Absorbance* or *Transmission* in dependency of the wavelength is measured and recorded. The wavelength range can be freely selected within the measuring range of the photometer. The increment is 1 nm.

A spectrum is recorded without using any methods or profiles. All settings are configured during measurement.

- **Baseline** Before a spectrum is recorded, a baseline has to be recorded using a suitable zero cell, e.g. with deionized water. The baseline has to cover at least the wavelength range of the spectrum to be recorded. Once the baseline is measured, it remains stored in the photometer until
  - a new baseline is recorded
  - the wave length range is extended in the [Setup] menu
  - the Spectrum mode is exited or the photometer is switched off
- Settings You can record a spectrum with standard settings without opening the setting window.

The following settings are possible for a	a spectrum:	
---	-------------	--

Input field	Possible entries
Wavelength start	190* 1100 nm
Wavelength stop	190 1100* nm
Mode	Absorbance* or Transmission
Smoothing	Yes* or No
Color of graph	Color selection for the curve
Scaling	Auto* or Manual
Scaling: Auto*	During measurement, the instrument adjusts the axis scaling (minimum and maximum value of the axis) to the measured values. The entire curve is always visible.
Scaling: Manual Y-axis min Y-axis max	The axis scaling (minimum and maximum value of the axis) is set manually.

\* default setting



You can store the current settings as a profile with [Save]. You can load a stored profile with [Open]. Profiles for spectra have the file extension, ".profil".

Spectrum	
Spectrum	16.01.15 9:52
4.0 Spectrum 2 You have to record a baseline firs ( <zero>). 1 Adjustment of wavelength range <general setup="">. 0.0 -1.0 400 Wavelength [nm]</general></zero>	it
Setup	Open
Spectrum	16.01.15 9:52
Wavelength start Wavelength stop Mode Smoothing Color of graph Scaling	1100 nm Absorbance Yes Blue Auto
	Apply

#### 


### 4.9.2 Recording a Spectrum


<HOME>


A message containing operating instructions is displayed.

- 1 Open the setting menu with [Setup].
- 2 If necessary, change the standard settings for the spectrum.
  - Wavelengths for start point and end point of the spectrum to be recorded
  - Representation mode (Absorbance / Transmission)
  - Smoothing of the curve (Yes/ No)
  - Color of the curve
  - Scaling of the Y-axis
     Auto: (entire value range)
     Manual: (selected value range)
- **3** Accept all entries with [Apply].

A message containing operating instructions is displayed.







Recording the baseline:

- 4 Press the <ZERO·BLANK> key. The photometer records the baseline.
- 5 Wait until the baseline is completely recorded.

The photometer is ready to measure after the baseline has been recorded.

Recording the spectrum:


6 Depending on the type, insert the cell as follows:


Round cell:

Insert the round cell in the round cell shaft so that it touches the bottom.


- 7 Close the inner turn-up lid.
- 8 Close the cell shaft cover.
- 9 Start the measurement with <**START·ENTER**>.

After the spectrum has been recorded, the following message appears: *Recording of spectrum is completed.* 





Cursor information



Rectangular cell:

Open the inner turn-up lid.

Insert the rectangular cell vertically so that it touches the bottom and left edge of the cell shaft. The opaque sides of the rectangular cell must point to the front and back.

- 10 Close the cell shaft cover.
- 11 Start the measurement with <START·ENTER>.
- **12** Wait until the spectrum is completely recorded.

At the end of the recording the following message appears: *Recording of spectrum is completed.* 

**13** Confirm the message with **<START·ENTER>**.

The cursor appears at the absolute maximum of the spectrum.

- 14 You have the following options:
  - Immediately edit the spectrum (see section 4.9.3)
  - With **<PRINT>**, you can output the spectrum to a connected printer as a graphic.
  - You can save the spectrum as a \*.csv file with <STORE>. As the storage location, you can select the photometer (*Internal DataB folder*) or a USB memory device connected to the USB-A connection (*USB memory*). Stored spectra can be recalled and edited at any time (see section 4.9.3).

#### 4.9.3 Loading/editing a spectrum

A spectrum can be edited immediately after measurement. Stored spectra

can be loaded and edited as well.

The following tools are available for editing:

- Cursor function for incremental moving along the curve with indication of the x and y values
- Zoom function to scale up a section
- Mathematical functions for various evaluating and calculating operations. The functions are described from page 102.

Loading a stored spectrum

<HOME> Spectrum – [Open]


Open (Inte	rnal DataB folder)	16.01.15 9:52
I		
26.02.07	Holmium.csv	
23.02.07	K2Cr2O7_340nm.csv	
Location	Delete	
		1


The list with the spectra stored in the exchange memory is displayed.

- 1 If necessary, you can select a different memory location for the spectrum with *[Location]* (USB memory device at the USB-A connection).
- 2 Select the required spectrum. The original view of the curve is displayed.

The cursor consists of a horizontal and vertical line that cross each other on a point of the curve. A box names the x and y values of the point of the curve.

Move the cursor along the x axis (wavelength) with  $< \P > < \triangleright >$ . You can scan and evaluate the curve point after point.





#### ba77094e01 06/2015

- Extreme values (zoomed area) Highlights the extreme values (minimum and maximum values) of the displayed spectrum.
- Mark points
  - Opens an edit mode where you can highlight individual points of the spectrum.

With the *[Mark]* function key you can highlight individual points. The wavelength and measured valued are displayed at the highlighted point.

With the [Delete] function key you can remove individual points.

- Delete all marks Erases all highlighted points in the spectrum.
- Original Displays the original, unedited spectrum.
- Integral

Calculates the area between the zero line and curve within a freely selectable wavelength interval [X1,X2].

• Derivative

Calculates the derivative of the total spectrum. To calculate the second and third derivative, the function can be carried out several times.

• Compare spectrum

Loads a second spectrum into the same diagram for direct comparison. The second spectrum is displayed in the color magenta.

- Add spectrum Adds a stored spectrum to the current spectrum.
- Subtract spectrum Subtracts a stored spectrum from the current spectrum.
- Divide spectrum (ratio)
   Divides the absorbance or % transmission values of the current spectrum
   by the values of a stored spectrum
- Add fixed value Adds a constant absorbance or % transmission value to the current spectrum.
- *Multiply fixed value* Multiplies the absorbance or % transmission values of the current spectrum by a constant value.




The addition, subtraction and division of two spectra always applies to the common wavelength range of both spectra only.

#### 4.9.4 Saving / exporting a spectrum

The saving of a spectrum saves both the edited and the original spectrum. Consequently, the original spectrum can be restored from each stored spectrum.

#### Storing in memory



1 Record a spectrum (see section 4.9.2) or Load a stored spectrum (see

section 4.9.3).

- 2 If necessary, connect a USB memory device to the USB-A interface.
- 3 Open the save dialog with **<STORE>**.
- If necessary, change the storage location with [Location]: Internal DataB folder.
   Exchange folder in the instrument or USB memory: USB memory device connected at the USB-A connection.
- If necessary, change the file name.
   The photometer automatically suggests an unambiguous file name consisting of the
- wavelength range, date and time.6 Save the file with
  - START·ENTER>.

**Export to a PC** To export a stored spectrum to a PC: see section 4.13.3

#### 4.10 Kinetics

The Kinetics function enables the temporal tracing of the absorbance or transmission of a sample at a certain wavelength.

The photometer automatically calculates the slope between two adjacent measuring points from the available measurement data. The catalytic activity can also be determined and displayed if required.

To record the kinetics, the photometer carries out single measurements at regular intervals (measuring interval) and stores the measured values as a time function.

All settings for a recording are administrated as a profile. Profiles can be created, stored, edited and deleted. Each measurement requires a respective profile.

#### 4.10.1 Creating/editing profiles for Kinetics recordings



Profiles for Kinetics records are stored under the numbers 4001 to 4020.

In the delivery condition, a profile is stored for demonstration purposes.

A profile for a Kinetics recording comprises the following data:

Input field	Possible entries
Number *	4001 4020
Name	Any name (max. 18 characters)
Mode *	Absorbance or Transmission
Wavelength *	Freely selectable (in nm)
Duration *	Total duration in the format hh:mm:ss (hours:minutes:seconds)
Interval *	Measuring interval = time interval between two successive single measurements in the format hh:mm:ss (hours:minutes:seconds)
	Exception: With the setting, <i>Measurements/interval</i> : <i>Max/</i> <i>interval</i> the interval is defined differently (see below).
Delay	Time between the start of the recording and the start of the first single measurement

Input field	Possible entries
Scaling	Auto or Manual
Scaling: Auto **	During measurement, the instrument adjusts the axis scaling (minimum and maximum value of the axis) to the measured values. The entire curve is always visible.
Scaling: Manual Y-axis min Y-axis max	The axis scaling (minimum and maximum value of the axis) is set manually.
Measurements/interval	1/interval or Max/interval
	Here you define how many measurements are carried out per interval.
	This setting has an impact on the calculation of the slope of the individual intervals (see section 4.10.6).
Catalytic activity	Yes or No
(only with <i>Mode</i> : <i>Absorbance</i> )	Here you determine whether the catalytic activity should be calculated.
	The catalytic activity is a measure for the amount of substance that is converted per time unit. To accelerate the substance conversion, a catalyst or enzyme (biological catalyst) is used in most cases.
	Carry out this measurement at room temperature.
Catalytic activity: Yes	
Factor Unit Resolution	The catalytic activity or enzymatic activity is calculated from the slope of the curve.
	Cat. A. = mean value <i>Slope</i> [Δ/min] * Factor
	Here you can enter the value for <i>Factor</i> .
	The calculated value for the catalytic activity is displayed in the menu, [Edit] / Slope & catalytic activity, together with the unit and resolution selected here.
* necessary inputs ** default: <i>Auto</i>	

### Creating/editing a profile


<home></home>
Kinetics
– [Setup]
<ul> <li>Edit profile</li> </ul>

Edit profile (1 c	of 2)	16.01.15 9:52
Number		4001
Name		NADH
Mode		Absorbance
Wavelength		340 nm
Duration		02:00:00
Interval		00:00:30
Delay		00:01:00
Scaling		Auto
Profile list	Delete	Next

1 Enter the data for the profile here. The next available profile number is already entered as the number.

You have the following options when filling out the input fields:

- Fill out all empty input fields one after the other
- Using [*Profile list*], select an already existing profile as a model, give it a new profile number and adjust the entries
- Using [*Profile list*], select an existing profile in order to change it (without changing the number).
- You can delete the profile completely with [Delete].
- 2 With [Next] you can switch to further settings.
- **3** Enter further data for the profile here.
- 4 Accept all entries with [Complete].
  - The profile is created and selected. The photometer is ready to measure.





The *Catalytic activity* function is only available if the Absorbance mode was selected.

<home> Kinetics  - [Profi</home>	le list]			
Select profile (all)		16.01.15 9:52		The li The p to the
4001 NADH 4002 A740	Absorbance Absorbance			<u>Selec</u>
			1	Selec < <b>≜</b> >< displa
Last used			2	Accep < <b>STA</b>
	1	1		

#### 4.10.2 Loading a profile for kinetic recording

To load a profile for kinetic recording, proceed as follows:

The list of profiles is displayed. The profiles are ordered according to the profile number.

Selecting a profile:

- Select the required profile with
   <▲><▼>. The active selection is displayed in reverse video.
- 2 Accept the selection with **START·ENTER**>.

The photometer is ready to measure.

## Narrowing down the list of profiles

If the list is very long, you can narrow down the profile list and thus make the search easier as follows:

- Using [Last used], you can restrict the profile list to the ten profiles last used.
- With the search function you can search certain character strings in the list. The search takes place as a full-text search of the entire list contents. Thus you can search for a profile number or name.

#### Search function

bsorbance	
bsorbance	

#### Search for a character string:

Enter the character string to be searched for in the search window with **<A...9>**.

The list appearing below shows all hits containing the character string. The hit list is updated with each character that is entered.




Note the case sensitivity when searching.


<home> Kinetics</home>			
KINEUCS			
Kinetics	16.01.15 9:52	profi	cessary, select a different le with <i>[Profile list]</i> (see on 4.10.2).
Zero measurement requirec Press <zero blank=""></zero>	11		t the zero measurement wit <b>RO·BLANK&gt;</b> .
4001: Demo Setup Profile list	Absorbance Open		
Kinetics	16.01.15 9:52	3 Carr	y out the zero measuremer
Please insert zero cell (disti or press <start enter=""></start>	lled water)		
10 mm Setup Profile list	0.0 - 22.1 mg/l		
Kinetics	16.01.15 9:52 [ZERO 16.01.15 9:51]		photometer is ready to sure.
Close the cell shaft cover. To start measurement, pres ENTER>.	ss <start <="" td=""><td></td><td></td></start>		
10 mm	0.0 - 22.1 mg/l		

Open

# 4.10.3 Recording the Kinetics

Setup Profile list





4 Depending on the type, insert the cell as follows:

### Round cell:

Insert the round cell in the round cell shaft so that it touches the bottom.

If the inner turn-up lid is opened too wide, a message prompts you to close the inner turn-up lid.

### Rectangular cell:

Open the inner turn-up lid.

Insert the rectangular cell vertically so that it touches the bottom and left edge of the cell shaft. The opaque sides of the rectangular cell must point to the front and back.


- 5 Close the cell shaft cover.
- 6 Start the measurement with <START·ENTER>.

The photometer starts recording automatically.

7 Wait until the recording is finished.

Stopping the recording:

- Use [Stop] to terminate the recording prematurely. The curve recorded up to this point can be stored and edited (see section 4.10.6).
- Use **<ESC>** to completely cancel measurement. The curve recorded up to this point is discarded.




8 After the specified *Duration* has expired, the cursor appears.

You have the following options:

- You can move the cursor along the curve and have the measurement data for each point displayed (see section 4.10.6)
- With **<PRINT>**, you can output the kinetic curve to a connected printer as a graphic.
- You can store the kinetic curve with **<STORE>** (see section 4.10.4).
- Execute further functions to edit the kinetic record (see section 4.10.6)
- Close the kinetic record with **<ESC>**.

Storing in memory



# 4.10.4 Saving / exporting a Kinetics record

1 Carry out the kinetic recording (see section 4.10.3) or

Load a stored kinetic record (see section 4.10.4).

- 2 If necessary, connect a USB memory device to the USB-A interface.
- 3 Open the save dialog with <STORE>.
- If necessary, change the storage location with [Location]: Internal DataB folder.
   Exchange folder in the instrument or USB memory: USB memory device connected at the USB-A connection.
- 5 If necessary, change the file name.
- 6 Save the file with <START.ENTER>.

Export to a PC Export a stored kinetic record to a PC: see section 4.13.3 Example of a 6|4001|1|1|525|1280913092|59|5|1|0.000|0.301|0|1.000|µkat|2 kinetic recording (\*.csv file) Device: Serial number:Software: User: photoLab 7600 UV-VIS 09130512 1.30-WTW-1.60 Administrator Start time Wavelength [nm] 04.08.2010 11:11525 Time [s] Absorbance 0 0,092 5 0,077 10 0,073 15 0,069 . . . . . . .

#### Line 1 - explanations:

Column	Value	Explanation
1	6	Version of the file format for the CSV file
2	4001	Profile number
3	1	Measurement of absorbance (0) or transmission (1)
4	1	Measurement once per interval (0) or as often as possible (1)
5	525	Wavelength (in nm)
6	1280913092	Start time (internal data format)
7	59	Duration (in sec)
8	5	Interval time (in sec)
9	1	Scaling automatic (0) or manual (1)
10	0.000	Minimum for manual scaling
11	0.301	Maximum for manual scaling
12	0	Enzymatic activity Off (0) or On (1)
13	1.000	Factor for enzymatic activity
14	µkat	Unit of enzymatic activity
15	2	Decimal points for enzymatic activity

### 4.10.5 Loading a kinetic record

You can load and view stored kinetic records.


Loading a stored Kinetics record <HOME> Kinetics – [Open]

		16.01.15	9:52
26.02.07	Enzyme kinetics.csv		
24.02.07	A740.csv		
24.02.07	kinetics_4002_070224_1410.csv		
Location	Delete		
Location	Delete		

The list with the stored Kinetics records is displayed (*Internal DataB folder*).

- 1 With [Location] select the memory location of the kinetic record (Internal DataB folder or USB memory for a USB memory device at the USB-A connection).
- 2 Select the required Kinetics record.

The curve is loaded.



You have the following options:


- You can move the cursor along the curve and have the measurement data for each point displayed (see section 4.10.6)
- With **<PRINT>**, you can output the kinetic curve to a connected printer as a graphic.
- You can store the kinetic curve with **<STORE>** (see section 4.10.4).
- Execute further functions to edit the kinetic record (see section 4.10.6)
- Close the kinetic record with **<ESC>**.

# 4.10.6 Editing a kinetic record

The following functions are available for kinetic records:

- Moving along the curve with the cursor
- Displaying a list with the slopes of the curve for each interval
- Scaling the Y-axis of the diagram
- Combined display of two kinetic records in one graphic
- Display of the difference of two kinetic records

### Cursor Cursor information



The cursor consists of a horizontal and vertical line that cross each other on a point of the curve. A box names the x and y values of the point of the curve.

Move the cursor along the x axis (time axis) with  $< \P > < \triangleright >$ . You can scan and evaluate the curve point after point.

Slope of the curve & catalytic activity

The function, *Slope & catalytic activity* indicates the slope of the kinetic curve in the individual intercepts (intervals) of the curve. An intercept corresponds to the *Interval* entered in the profile.

A 16 01 15 9:52

cat

Time

5 s

10 s

15 s

20 s

25 s

30 s

1 Indicate the slope of the kinetic curve in the individual intercepts (intervals) of the curve with [Edit]/ Slope & catalytic activity.

If the calculation of the catalytic activity was selected when the profile was created it is displayed here together with the slope.



Back

Kinetics

Interval 1

2

3 4 5

6

0.63

Slope [ $\Delta$ /min] ( $\Delta$ /

0.000

0.000

0.000

0.000

0.000

0.000

The *Slope & catalytic activity* function is only available if the kinetic recording was done in the Absorbance mode.

The displayed slope for an interval is determined as follows, depending on the slope:

Measurements/interval	Slope
1/interval	Slope, converted to the interval "1 minute"
Max/interval	Slope of the straight line determined by linear regression in an interval, converted to the interval, "1 minute"

Scaling of the Y-axis You can manually determine the scaling of the Y-axis with [Setup]/Scaling/ Manual.

*Compare kinetics* For direct comparison, you can load a second kinetic record into the same diagram with *[Edit] / Compare kinetics.* 



The *Compare kinetics* function can only be carried out if both kinetic records were made in the Absorbance mode.

*Subtract kinetics* You can subtract a stored kinetic record from the current kinetic record with [*Edit*] / *Subtract kinetics*.



The *Subtract kinetics* can only be carried out if both kinetic records were made with the following settings:

- Mode: Absorbance
- Measurements/interval: 1/interval
- Equal interval

# 4.11 Timer

You can use the timers to remind you with an audio signal of a time interval that has expired.

The photometer has two types of timers:

- The User defined timer is a timer that can be freely assigned. The interval and name can be freely set. Only one freely assignable timer is available. It cannot be erased (see section 4.11.1).
- Analysis timer are timers permanently stored in the photometer. The names and intervals of the analysis timers are stored in the method data of a measuring method (*Concentration* mode). The number of available analysis timers corresponds to the number of reaction times prescribed in the analysis instructions of the programmed methods (see section 4.11.2).

The photometer administrates all timers in the timer overview.

The timer overview (the *Timer* menu) is opened with the **<TIMER>** key. The *Timer* menu can be opened from any operating situation. Operation of the timer does not disturb any other functions. The timer overview can be exited with the **<ESC**> key.

When the *Timer* menu is opened for the first time, only the user-defined timer is in the timer overview. You can include analysis timers into the list or remove them according to your requirements (see section 4.11.2).

The timer overview displays the status of each timer and, of a started timer, the remaining time of the specified time interval.

All timers are started manually.

As soon as one single timer has been started the timer symbol appears on the display in all operating modes.

As soon as a timer is started, it is given the timer status Active.

When the specified time interval has expired, the timer status changes from *Active* to *Expired* and an audio signal sounds.

In the timer status *Expired* the audio signal sounds until the timer is stopped manually.

After the stop, the timer status changes to *Inactive* and the audio signal is switched off.

#### 4.11.1 User defined timer

If you want to manually enter time intervals, use the *User defined timer* function.

# <TIMER>

Timer		16.01.15 9:52
Designation	<b>T</b> i	Chathing
Designation	Time	Status
User defined timer	00:15:00	Inactiv
A6/25 - 1	00:15:00	Inactiv
Start Sto	op Edit	Add

The *Timer* menu is open.

- 1 Highlight the User defined timer.
- 2 If necessary, change the name and time of the timer with [Edit].
- **3** Start the highlighted timer with *[Start]*.

The status of the timer is *Active*. When the specified time interval has expired, and audio signal sounds and the timer status changes to *Expired*.

4 Stop the highlighted timer with [Stop].

The status of the time changes to *Inactive*. The audio signal is switched off.

# 4.11.2 Analysis timer

Between the individual steps of a measurement, reaction times often have to be observed. The length of the reaction time is defined in the relevant analysis instructions.

For all required reaction times, the analysis timers with the corresponding time intervals are stored in the instrument. The names of the analysis timers include the method name and a current number so several timers within a method can be distinguished from each other.

To be able to use an analysis timer for a method, you have to load it first in the timer overview.

To do so, first select the required method and then add the available analysis timers to the timer overview so they can be started as necessary.

The timer overview always comprises the free timer and the selected analysis timers.

1 Select the required method in the *Concentration* mode.

Manual selection of the method (see section 4.5.5).

# <TIMER>

Timer			
Designation		Time	Status
User defined tim	ier	00:15:00	Inactive
A6/25 - 1		00:15:00	Inactive

2 Open the Timer menu.

The *Timer* menu is open.

3 If necessary, add a new timer to the list with [Add].

#### Note:

The *[Add]* function key is only displayed if a method is selected for which analysis timers were programmed but are not yet displayed in the list of timers.

- 4 Highlight an analysis timer.
- 5 If necessary, remove the analysis timer from the list with [*Remove*].
- 6 Start the highlighted timer with [Start].


The status of the timer is *Active*. When the specified time interval has expired, an audio signal sounds and the timer status changes to *Expired*.

7 Stop the highlighted timer with [Stop].

The status of the time changes to *Inactive*. The audio signal is switched off.



4.12.1 Overview



Measured data	Save, back up, export
Concentration, Absorbance / % Transmission	Measurement datasets of these measuring modes are first stored in the measured value memory of the photometer (5000 memory locations) with <b><store></store></b> or <i>AutoStore</i> .
<i>Special / Multi wavelengths</i> <i>OptRF measurement</i>	The measured value memory is available from the <i>Measurement data memory</i> menu. Here you can view, filter and export into a PC- readable file (*.csv) the stored measurement datasets ( <b>STORE</b> >).
	Csv files of these measuring modes cannot be reimported to the photometer.
	Measurement datasets of these measuring modes can also be stored to a pdf file (see section 4.12.11).
Spectrum Kinetics	You can store and export measurement data of these measuring modes directly as a PC-readable file (*.csv) with <b><store></store></b> .
	Csv files of these measuring modes can be reimported and displayed on the photometer.
	Measurement data of these measuring modes can also be stored to a pdf file (see section 4.12.11).
AQA records	You can store and export measurement data of these measuring modes directly as a PC-readable file (*.csv) with <b><store></store></b> .
	Csv files of records cannot be reimported to the photometer.
	Measurement data of these measuring modes can also be stored to a pdf file (see section 4.12.11).
User-defined methods / profiles	Method data and profile data are stored and exported with the <i>Exchange methods/profiles</i> function in the <b><home></home></b> / <i>General setup</i> menu.

For each export procedure you can select the location where the PCreadable files (\*.csv, \*.pdf) should be stored: either to the photometer (*Internal DataB folder*) or an external memory (*USB memory*). On an external memory, the data are stored in the directory "DataB photoLab 7....".

The files stored in the photometer (*Internal DataB folder*) can later be transferred to a connected PC or to an external memory (*USB memory*).

# 4.12.2 Instructions on using USB memory devices

The safety of data stored on USB memory devices depends on the quality of the memory device and the data transmission. Data is stored partly or not at all if for example:

- The power supply of the external memory device is interrupted during the write process, or
- The external memory device is prematurely disconnected from the photometer during the data backup.

To prevent a data loss we recommend the following:

- Save all data internally in the photometer first.
- After performing a backup leave the USB memory device connected to the photometer for some time.
- Check whether the stored data is complete, e.g. on a PC.
- Use the USB memory device for data transport but not for permanent data storage.

#### 4.12.3 Measurement datasets

Elements of a	A complete measurement dataset consists of:
measurement dataset	<ul> <li>Consecutive number (is automatically assigned by the photometer)</li> </ul>
Galasei	Date/time
	<ul> <li>Identification (e.g. ID or "AutoStore")</li> </ul>
	• User name
	<ul> <li>Measured parameter, e.g. method number, dilution, wavelength (depending on the measuring mode)</li> </ul>
	<ul> <li>Measured value with unit and, if necessary, citation form</li> </ul>
<b>Operations with</b>	Measurement datasets can be
measurement datasets	<ul> <li>stored (see section 4.12.4)</li> </ul>
Galaseis	<ul> <li>displayed and printed (see section 4.12.6)</li> </ul>
	<ul> <li>filtered, i.e. selected or hidden based on certain criteria (see section 4.12.7 and section 4.12.8)</li> </ul>
	<ul> <li>deleted (see section 4.12.9).</li> </ul>
If the memory is full	You can erase measurement datasets (see section 4.12.9), or overwrite the oldest dataset with the next storing procedure. A security prompt appears before a dataset is overwritten. To backup the measurement data, you can transmit the measurement datasets from the measurement data memory to the internal DataB folder or a USB memory device connected to the USB-A

connection and archive them further from there (see section 4.13.3).

### 4.12.4 Saving measurement datasets manually

After each measurement, you can store the measurement data manually with the **STORE**> key. The data are stored in the measurement data memory. The memory symbol in the header indicates that the measurement data displayed on the screen are ready to be stored. With the measuring modes, *Concentration, Absorbance / % Transmission,* and *Special / Multi wavelengths* you have the additional option to automatically store all new measured values at the time of the measurement (*AutoStore,* see section 4.12.5).

**Storing with identification (ID)** When storing manually, an input field for the identification (ID) appears after pressing the **<STORE>** key. Here you can enter an individual combination of alphanumeric characters for later easier identification of the measurement datasets. 30 digits are available for this.

The following measurement data are stored in the measured value memory automatically (see section 4.12.5) or manually (with the **<STORE>** key, see section 4.12.4):

- Concentration
- Multi wavelength
- OptRF measurement
- Absorbance / % Transmission

The data stored in the measured value memory can be filtered with filter criteria and then exported to the PC-readable \*.csv format.

The photometer automatically offers a file name during the storage procedure.

Example: Saving data from the measured value memory

#### <HOME>

Concentration, Absorbance / % Transmission, Special / Multi wavelengths or OptRF measurement [Setup] Measurement data memory

		2	<pre>STORE&gt;.</pre>
			The photom proposes th <i>DataB folde</i>
Save (Internal DataB folder)	16.01.15 9:52	3	If necessary with [Location of the second se
		4	If necessary proposed file
		5	Save the me
Location			The data are If the photor <i>folder</i> ) is sel

- 1 If necessary, set the filter criteria with [Setup].
- 2 Open the save dialog with

neter automatically ne location Internal erand a file name.

- y, change the location ion] (USB memory).
- y, change the le name.
- easurement data with NTER>.

re stored. meter (Internal DataB elected as the location. the data can then be copied to a USB memory device (see section 4.13.1).

#### 4.12.5 Saving measurement datasets automatically

For the measuring modes, Concentration, Absorbance / % Transmission, and Special / Multi wavelengths you can record every measured value automatically (AutoStore). The AutoStore function is active in the default condition.

All automatically stored measurement datasets are given the ID "AutoStore". The "AutoStore" ID is overwritten if the same measured value is manually stored afterwards (**<STORE>**).

This ensures that every measurement dataset is stored in the data memory only once.

Activate or deactivate the AutoStore function as follows:

AutoStore Activating or deactivating the function

### <HOME>

Concentration. Absorbance / % Transmission, Special / Multi wavelengths or **OptRF** measurement

- [Setup]

Measurement data memory

The available functions are displayed.

1 Select and confirm AutoStore.

The *AutoStore* function is active  $(\checkmark)$  or inactive (no checkmark).



The AutoStore setting applies to the measuring modes Concentration, Absorbance / % Transmission, Special / Multi wavelengths and OptRF measurement in general.

#### 4.12.6 Displaying measurement data memory

Depending on the operating situation, you can recall the measured value memory as follows:

 From the main menu
 <HOME>

 [Setup],

 Measurement data memory

 From a measuring mode
 Concentration,

 Absorbance / % Transmission,
 Special / Multi wavelengths or

OptRF measurement - [Setup] - Measurement data memory

Each of these options displays the contents of the measurement data memory as a list as follows.

Measurement data memory	<u>a</u> e	16.01.15 9:52
15.01.15 14:00 3.50 mg/l Ni 15.01.15 14:05 3.64 mg/l Ni 15.01.15 14:10 3.69 mg/l Ni 15.01.15 14:15 3.72 mg/l Ni 15.01.15 14:20 3.72 mg/l Ni 15.01.15 14:25 3.75 mg/l Ni 15.01.15 14:35 3.73 mg/l Ni 15.01.15 14:35 3.80 mg/l Ni	Administrator Administrator Administrator Administrator Administrator Administrator Administrator Administrator Administrator	AutoStore AutoStore AutoStore AutoStore AutoStore AutoStore AutoStore AutoStore AutoStore AutoStore
Filter ✓ Memory space usage: 9/5000 Setup Single value	Delete	

If there are more datasets available than can be displayed, the arrows  $\blacktriangle$  and  $\blacktriangledown$  are displayed additionally.

*Filter* ✓ indicates that the filter settings are active. In this case, only those datasets are displayed that correspond to the selected filter criteria (see section 4.12.7).

#### Options

Measurement datasets can be

- displayed in short form as a list or in details as individual values ([List] <->
  [Single value])
- filtered (see section 4.12.7 and section 4.12.8)
- deleted (see section 4.12.9).
- with <STORE>, you can store the entire displayed list as a \*.csv file in the internal DataB folder or on a USB memory device connected to the USB-A connection. The filter settings apply to the storing process. You can freely select the file name. Thus you can, e. g. store in a separate file and systematically archive measurement data of a certain period.
- with **<PRINT>**, the entire displayed list can be printed. The filter settings apply to the print process.

### 4.12.7 Filtering measurement datasets

The functions to display, delete and download stored measurement datasets refer to all stored measurement datasets that correspond to the specified filter criteria.

### Filter criteria The following filter criteria can be set:

- *Mode* (measured parameter)
- User
- *ID* (identification)
- Date (date from ... to ...)
- *Method* (for the measured parameters, *Concentration* and *Multi wavelength*)

<home></home>
Concentration,
Absorbance / % Transmission,
Special / Multi wavelengths or
OptRF measurement
– [Setup]
<ul> <li>Measurement data</li> </ul>
memory
– Setup

Filter
--------

- Mode (Concentration, etc.)
- User (**<A...9**>)
- ID (<**A...9**>)
- Date (from ... to ...)
  - Method

Measurement	data memory	8 E	16.01.15 9:52
$\begin{array}{c} \textbf{15.01.15} \ \textbf{14:00} \\ \textbf{15.01.15} \ \textbf{14:00} \\ \textbf{15.01.15} \ \textbf{14:10} \\ \textbf{15.01.15} \ \textbf{14:11} \\ \textbf{15.01.15} \ \textbf{14:12} \\ \textbf{15.01.15} \ \textbf{14:22} \\ \textbf{15.01.15} \ \textbf{14:22} \\ \textbf{15.01.15} \ \textbf{14:23} \\ \textbf{15.01.15} \ \textbf{14:33} \\ \textbf{15.01.15} \ \textbf{14:34} \\ \textbf{15.01.15} \ \textbf{14:44} \end{array}$	5 3.64 mg/l Ni 5 3.69 mg/l Ni 5 3.72 mg/l Ni 5 3.72 mg/l Ni 5 3.75 mg/l Ni 5 3.73 mg/l Ni 5 3.80 mg/l Ni	Administrator Administrator Administrator Administrator Administrator Administrator Administrator Administrator	AutoStore AutoStore AutoStore AutoStore AutoStore AutoStore AutoStore AutoStore AutoStore
Filter 🗸 Memory space Setup	usage: 9/5000 Single value	Delete	

The filter setting menu is displayed.

- **1** Set the filter criteria.
- 2 If necessary, deactivate any selected filter criteria with [Reset entry].
- **3** Confirm the filter selection with *[Apply]*.

The *Measurement data memory* list is displayed.

The following information is displayed additionally:

- Current memory occupancy
- Active filter criteria (*Filter* ✓)



Alternatively, you can <u>hide</u> measurement datasets that meet the specified filter criteria with the *Selected values: invert selection* function (see section 4.12.8).

### 4.12.8 Inverting filters

With the *Selected values: invert selection* function you can <u>hide</u> all measurement datasets that correspond to the specified criteria of the filter (see section 4.12.7).



You can use this function to select and delete measurement datasets no longer used.

<HOME> Concentration, Absorbance / % Transmission, Special / Multi wavelengths or OptRF measurement [- [Setup] - Measurement data memory - Setup

Measurement	data memory	4e	16.01.15 9:52
15.01.15 14:00	57	Administrator	AutoStore
15.01.15 14:05	51	Administrator	AutoStore
15.01.15 14:10	3.69 mg/l Ni	Administrator	AutoStore
15.01.15 14:15	3.72 mg/l Ni	Administrator	AutoStore
15.01.15 14:20	3.72 mg/l Ni	Administrator	AutoStore
15.01.15 14:25	3.75 mg/l Ni	Administrator	AutoStore
15.01.15 14:30	3.73 mg/l Ni	Administrator	AutoStore
15.01.15 14:35	3.80 mg/l Ni	Administrator	AutoStore
15.01.15 14:40	3.78 mg/l Ni	Administrator	AutoStore
Filter 🗸 Memory space	usage: 9/5000		
Setup	Single value	Delete	

The *Measurement data memory* list is displayed. All measurement datasets corresponding to the filter criteria are hidden.

# 4.12.9 Erasing stored measurement datasets

If you no longer need any stored measurement datasets, you can erase them individually or altogether.

# <HOME>

Concentration, Absorbance / % Transmission, Special / Multi wavelengths or OptRF measurement

– [Setup]

Measurement data memory

15.01.15 14:00 3.50 mg/l Ni	Administrator	AutoStore
15.01.15 14:05 3.64 mg/l Ni	Administrator	AutoStore
15.01.15 14:10 3.69 mg/l Ni	Administrator	AutoStore
15.01.15 14:15 3.72 mg/l Ni	Administrator	AutoStore
15.01.15 14:20 3.72 mg/l Ni	Administrator	AutoStore
15.01.15 14:25 3.75 mg/l Ni	Administrator	AutoStore
15.01.15 14:30 3.73 mg/l Ni	Administrator	AutoStore
15.01.15 14:35 3.80 mg/l Ni	Administrator	AutoStore
15.01.15 14:40 3.78 mg/l Ni	Administrator	AutoStore
Filter 🗸		
Memory space usage: 9/5000		

The *Measurement data memory* list is displayed.

The filter settings used last are active.

#### Erasure functions The follo

The following erasure functions are available.

- Erasing an individual measurement dataset
   Highlight a measurement dataset
   Remove the highlighted measurement dataset with [Delete].
- Erasing all measurement datasets of the displayed list
   1 Open the setting menu with [Setup].

Erasing all measurement datasets

 2 Select and confirm Delete memory (selected values only).

All measurement datasets corresponding to the current filter criteria are erased.

or

Select and confirm *Delete memory (all values)*.

All measurement datasets are erased.

# 4.12.10 Saving kinetic recordings, spectra and AQA files

After the following measurements, the *Save* dialog opens and prompts you to save the data in a \*.csv file:

- Kinetics
- Spectrum
- AQA3/MatrixCheck

If the data are not saved in \*.csv format, they are lost when the measuring mode is terminated.



During a kinetic recording, the current measurement is always saved in the file, "KineticsBackup.csv" for safety reasons.

### 4.12.11 Saving data as a pdf file

All data that can be printed (printer symbol on the display) can also be saved as a pdf file. The pdf file contains the data that are also output to a USB printer. Kinetic recordings and spectra are stored in the pdf file as a graphic.

Saving as a pdf file and printing is done with the **<PRINT>** key. Prerequisite is that pdf printing is set as the printer in the menu **<HOME>**/*General setup*/ *Data transfer/Printer/Function of PRINT key*.

Subsequently, enter a file name and select the storage location (internal folder DataB or USB memory device).

# 4.13 Saving / exporting files

If you want to back up or process measurement data files outside the photometer, you can copy them to external media.



Please follow the instructions on using USB memory devices (see section 4.12.2).

#### 4.13.1 Copying all measurement data files to a USB memory device

Even if no PC is directly connected to the photometer, you can very simply transfer all measurement data files from the photometer (*Internal DataB folder*) to a connected USB memory device.

<ho< th=""><th>ME&gt;</th></ho<>	ME>
[Setı	ıp]
	Save data to USB memory device

When the data saving procedure is finished, a message appears.

1 Confirm the message with **<STORE>**.

All measurement data files from the photometer (*Internal DataB folder*) have been transferred to the USB memory device.

The complete folder structure from the photometer is created on the USB memory device. The individual measurement data files are stored in subfolders sorted by measurement data types.

# 4.13.2 Copying user-defined methods / profiles to a USB memory device

### <HOME>

[Setup] – Exchange methods/profiles / Store to USB memory

device

A list is displayed that includes all user-defined methods and profiles available on the photometer. All methods and profiles are checked off with a checkmark. All methods and profiles checked off are saved.

1 If necessary, select individual methods/profiles with <▲><▼> and remove the checkmark with <START·ENTER>.

These methods/profiles will not be saved.

2 Start the save process with [Store].

A message appears when the data have been saved.

3 Confirm the message with **<START**.**ENTER**>.

The save process is completed. The data are stored in the *Exchange\_Method\_Profile* folder on the USB memory device. The individual files with the methods/ profiles are in subfolders.

Already existing files with identical names are overwritten without confirmation prompt.

# 4.13.3 Copying files to a PC

You can copy from the photometer to a PC the following data:

- Measured data
- Spectra
- Kinetic recordings
- AQA records
- User-defined methods
- Profiles

After saving measurement data in \*.csv or \*.pdf format, you can copy them to a PC. Measurement data in csv format can be directly imported to and processed in spreadsheets such as  $Microsoft^{®}$  Excel<sup>®</sup>.



Depending on the country variant, some spreadsheet programs require a certain decimal separator for the correct import of numerical values (comma or point). The decimal separator can be selected in the following menu:

<HOME> -> General setup -> Data transfer/Printer -> Decimal separator for csv-Files.

Files containing measurement data can be copied to a PC in the following ways:

- By using a USB memory device as a temporary storage (see section section 4.13.2 and section 4.13.1). Subsequently, you can connect the USB memory device to a PC and read out the data.
- By means of the "SpectralTransfer" program (see operating manual of the "SpectralTransfer" program). The "SpectralTransfer" program and the corresponding operating manual is provided on the enclosed CD-ROM.
- Using the "pHotoLab Data spectral" program (see operating manual of the "pHotoLab Data spectral" program). The "pHotoLab Data spectral" program is available as an accessory.

# 4.14 Importing files

You can import to a pHotoLab<sup>®</sup> 7xxx spectrophotometer the data that were created with the same or another pHotoLab<sup>®</sup> 7xxx spectrophotometer, and were then saved to a USB memory device or a PC.

You can import the following data:

- Spectra
- Kinetic recordings
- User-defined methods
- Profiles

# 4.14.1 Importing spectra or kinetic recordings from a USB memory device

You can import to the photometer any spectrum or kinetic recording by opening an externally stored spectrum or kinetic recording with the Open function of the photometer.

# 4.14.2 Importing methods / profiles from a USB memory device



When importing methods make sure that your photometer supports the wavelengths of the imported methods.

### <HOME>

[Setup]

Exchange methods/profiles
 / Import from USB memory device

A list is displayed including all user-defined methods and profiles stored in the corresponding subfolders of the Exchange directory on the USB memory device.

All methods and profiles are checked off with a checkmark. All methods and profiles checked off are imported.

 If necessary, select individual methods/profiles with <▲><▼> and remove the checkmark with<START·ENTER>.

These methods / profiles are excluded from importing.

2 Start the import with [Import].

A confirmation prompt appears before any data on the photometer are overwritten.

A message appears when the data have been imported.

3 Confirm the message with **<START**.**ENTER**>.

The import is completed. The imported methods / profiles are available on the photometer.

# 4.14.3 Importing files from a PC

You can import files from the PC to the photometer in the following ways:

- By means of the "SpectralTransfer" program (user-defined methods only) (see operating manual of the "SpectralTransfer" program). The "SpectralTransfer" program and the corresponding operating manual is provided on the enclosed CD-ROM.
- Using the "pHotoLab Data spectral" program (see operating manual of the "pHotoLab Data spectral" program). The "pHotoLab Data spectral" program is available as an accessory.

# 4.15 Printing the data (USB)

#### 4.15.1 Printer and terminal programs

**Usable printers** Data can be printed with standard printers (ink-jet or laser) connected to the USB-A interface. Suitable are all printers that can interpret the PCL-3 printer control language.

The printer symbol indicates that the display contents can be printed. To print, press **<PRINT>**.

**pdf file** As an alternative, you can also output the print data to a pdf file.

In den following paragraphs, "Print" means:



- output to a USB printer
- output to a pdf file.

#### 4.15.2 Settings for data transmission

Settings are possible for the data transmission to a printer or PC.

**Decimal separators** for CSV files For the output of CSV files you can select either a comma or a point as the decimal separator. The setting is made in the following menu:

**<HOME>** -> General setup -> Data transfer/Printer -> Decimal separator for csv-Files -> Comma (12,34) or Point (12.34).

**Short and long version** When printing measurement datasets, you can select a short or long version with different information contents. The setting is made in the following menu:

<HOME> -> General setup -> Data transfer/Printer -> Data format (print) -> Short or Extended.

**Printer** Here you can set which function is assigned to the **<PRINT>** key:

- Output to a USB printer
- Output as pdf file

The setting is made in the following menu:

<HOME> -> General setup -> Data transfer/Printer -> Function of PRINT key -> USB printer or PDF file.

### 4.15.3 Printing measurement datasets

This section describes how to print measurement datasets of the measuring modes, *Concentration*, *Absorbance / % Transmission*, and *Special / Multi wavelengths*.

By means of sample printouts, the printed information is described below:

Concentration and Special / Multi wavelengths mode

```
21 05.06.07 14:05:41 C4/25 844 mg/l COD Inlet
Administrator 0.005 02.06.07 11:02:13 2 AQA1: 9 AQA2: 14
```

Structure of the lines from left to right:

1st line:

[Consecutive no.] [Date] [Time] [Method name] [Measured value] [Unit] [Citation form] [Dilution] [ID or "AutoStore"]

2nd line (long version only):

[User] [Reagent blank value] [Date of blank value measurement] [Time of blank value measurement] [Lot ID of blank value measurement] [AQA1: label] [AQA1: record no.] [AQA2: label] [AQA2: record no]



Optional elements (e.g. dilution or ID) are output only if they were really used for measurement or storage.

Absorbance / % Transmission mode

14 05.06.07 11:25:01 445 nm 0.609 Absorbance AutoStore Administrator 0.133 02.06.07 09:59:01 AQA1: 9

Structure of the lines from left to right:

1st line:

[Consecutive no.] [Date] [Time] [Wavelength] [Measured value] ["Absorbance" or "Transmission" mode ] [ID or "AutoStore"]

2nd line (long version only):

[User] [Value of reference absorbance] [Date of reference absorbance] [Time of reference absorbance] [AQA1: label] [AQA1: record no.]



Optional elements (e.g. ID or reference absorbance) are output only if they were really used for measurement or storage.

# 4.15.4 Printing spectra



If you output a spectrum to a USB printer or as a pdf file, the current graphic representation is shown on the display.

# 4.15.5 Printing Kinetics records



If you output a kinetic record to a USB printer or as a pdf file, the current graphic representation is shown on the display.

# 4.16 Analytical quality assurance (AQA)

### 4.16.1 General information

The target of the analytical quality assurance (AQA) is to secure correct and precise measurement results.



Settings for AQA checks are only available for users of the user group, administrator. Every registered user can carry out the AQA check (see also section 4.17.1).

Analytical quality assurance (AQA) can be carried out in two steps independent of each other:

- AQA1: Monitoring of the photometer
- AQA2: Monitoring of the total system. It comprises the photometer, the test used, the accessories and the user's way of working.

The monitoring includes a check procedure that has to be successfully repeated by the user within a certain period (AQA interval).



The AQA monitoring is not active in the delivery condition.

AQA in measured value documentation All values that are measured after a passed check and within the AQA interval are given the addition *Protocol ID* in the measured value documentation. This addition is used to identify the relevant AQA inspection record.

# 4.16.2 Photometer monitoring (AQA1)

At least one set of test standards such as Spectroquant<sup>®</sup> PhotoCheck, Certipur<sup>®</sup> or other commercial test equipment (e.g. filter) is required for the photometer monitoring.

The administrator specifies which test standard has to be used as the minimum requirement for the AQA1 monitoring.

The extent of the monitoring can be enlarged with further test standards.



Settings for AQA checks are only available for users of the user group, administrator.

Every registered user can carry out the AQA check (see also section 4.17.1).

### Spectroquant<sup>®</sup> PhotoCheck

The PhotoCheck consists of 12 test standards in duplicate, 2 zero cells and 2 cells to check the barcode reader. Each PhotoCheck package contains a

lot dependent test certificate with all nominal values (absorbances) and tolerances of the test standards. These values are entered in the photometer during the configuration of the AQA1 check.

**Certipur® test** standards Each Certipur<sup>®</sup> standard is provided with a lot dependent test certificate with all nominal values (absorbances) and tolerances of the test standards. The values were preset in the factory.



Observe the shelf life of the test standards. The values in the photometer always have to be checked when a new package of test standard is used. If necessary, adjust the values at the photometer.

Overview of the	
photometer	
monitoring	

Photometer monitoring (AQA1) consists of the following parts:

- Configuring settings in the *AQA1 setup* menu.
  - Activate AQA1
  - Specify AQA1 Interval
  - Activate/deactivate the meter lock for missing or expired AQA1 check
  - Define the extent of the AQA1 monitoring by activating or deactivating the individual test standards.
  - Enter the nominal values, tolerances and lot numbers for the individual test standards
- Carrying out the AQA1 check. The photometer compares the results with the nominal values while taking into account the tolerances.

The steps are described in detail below.

Activating AQA1 The AQA1 monitoring is activated in the *Mode* menu:

<home></home>
[AQA]
├ AQA1 setup ├ Mode
– Mode

r				
AQA1 setup			16.01.15	9:52
Mode Interval Lock if Config AQA1 ina	ctive	Weeks		
Weeks				

Select and confirm *Weeks*. AQA1 is active. The *Interval* setting indicates *Weeks* as the interval unit.

# Defining the AQA1 Interval

The AQA1 Interval defines the interval between two AQA1 checks. When an interval has expired, the following consequences become effective:

- Warning and loss of the AQA1 labeling
- Locking of the photometer against all measurements (if activated).

<home></home>
[AQA]
⊢ AQA1 setup ⊢ Interval
– Interval

AQA1 setup		16.01.15 9:52
Mode Interval Lock instrument	Weeks	
Config Interval 2 Weeks		

Enter a numeric value (2 to 52 weeks) (<0...9>) and confirm
 The *Interval* defined for the AQA1 check is active.

# Configuring the lock of the photometer

Here you configure whether or not the photometer will be locked against all measurements if there is no valid AQA1 check or the interval for the AQA1 check has expired.

<home></home>
[AQA]
⊢ AQA1 setup ⊢ Lock instrument
– Lock instrument

AQA1	setup	16.01.15	9:52
Mode	Weeks		
Interv Lock i			
Config S	Should the instrument be locked for measurements if AQA1 check is inva has expired?		
	No		
	Yes		

1 Select and confirm *Yes*. The photometer is locked against all measurements if the AQA1 check is invalid or the AQA1 interval has expired.

Configuration of tests ...

⊢ AQA1 se ⊢ Config	etup Iuration o	f tests
		1
AQA1 setup		16.01.15 9:52
PhotoCheck	Active	
Certipur UV-VIS 1	Inactive	
Certipur UV-VIS 1A	Inactive	
Certipe Certipur UV-VIS 1		
Certipi General setup		
Certipi Activate		
L	1	

Apply

<HOME>

PhotoCheck		16.01.15 9:52
Lot number: Use by		HC616115 16.04.2008
	Target value	Tolerance
445/1	0.196	± 0.020
445/2	0.500	± 0.030
445/3	0.998	± 0.040
445/4	1.508	± 0.050
525/1	0.197	± 0.020
525/2	0.495	± 0.030
525/3	0.992	± 0.040
525/4	1.496	± 0.050
		Apply

All possible test standards or test standard sets are listed.

- 1 Select and confirm a test standard or test standard set.
- 2 Adjust and confirm the extent of the monitoring with *Activate* or *Deactivate*.
- **3** Confirm the test standard (set) once again.
- 4 Switch to the adjustment of the nominal values and tolerances with *Setup*.

# Example PhotoCheck:

- 5 Using <▲><♥> and <◀><▶>, select the Lot number, Target value or Tolerance entries and open them for editing with <START.ENTER>.
- 6 Enter and confirm the required value (<0...9>)
- 7 Accept all values with [Apply].

#### Carrying out the AQA1 check (example: PhotoCheck)

The AQA1 check comprises the check with all test standards activated in the menu *AQA menu / AQA1 setup / Configuration of tests ...* for AQA1 (see page 142).

First, a zero adjustment for all wavelengths takes place. Subsequently, the first individual checks with the selected test standards take place (e.g. PhotoCheck).

< <b>HOME</b> [AQA]  - AC	> QA1 che	ck				
PhotoCheck			16.01.15	9:52		
Reference measurement Please insert zero cell (distilled water).						
		1	1			

PhotoCheck 445/1			16.01.15	9:52
Please ins	ært PhotoCheci	k 445/1		

The photometer is ready for the zero adjustment.

1 Insert the zero cell. The cell is automatically recognized and the zero adjustment is started for all wavelengths.

> After the successful zero adjustment, the photometer is ready to measure the PhotoCheck test standard 445/1.

2 Insert the cell. The cuvette is automatically recognized and the measurement started.

After measuring, the measurement result, Target value, Tolerance and an evaluation (OK or failed) are displayed.

The photometer offers to repeat the measurement if the check failed.

If the check was successful, the measurement of the next PhotoCheck test standard, e.g. 445/2, appears on the display.

 Measure all test standards in the same way. After all test standards are successfully measured, the check is passed.

# **Test record** A test record is displayed after the check. It can be printed and stored as a file (in the internal DataB folder or USB memory device at the USB-A

connection, see section 4.12.1).

Sample printout:

```
photoLab 7600 UV-VIS 09130512 1.30-WTW-1.60 Administrator
AQA1
         ΟK
Protocol ID
                   9
Executed by:Administrator
                   16.01.2015
Executed
Valid until:
                  16.02.2015
PhotoCheckOC4790940K
445-1 0.200 +- 1000.192
445-2
        0.500 +- 2000.511
        1.000 +- 2001.006
445-3
        1.500 +- 2001.526
445-4
          0.200 +- 2000.247
525-1
          . . . . . . . . . . . . . . . . .
. . . . .
          . . . . . . . . . . . . . . . . .
(etc.)
```



Afterwards you can view the last AQA1 test record under AQA1 info.

# 4.16.3 Total system monitoring (AQA2)

For the total system monitoring, standard solutions with a defined analyte content are required (preferably certified Spectroquant<sup>®</sup> CombiCheck standards).



Settings for AQA checks are only available for users of the user group, administrator.

The AQA check can be carried out by any registered user.

Spectroquant<sup>®</sup> CombiCheck

Spectroquant<sup>®</sup> CombiCheck standards are multiparameter standards ready to use, i. e. they can be used for several test sets (methods).

In addition to the CombiCheck standards, one parameter standard solutions can also be used. They are prepared by dilution to the respective end concentration. The end concentration should be in the middle of the measuring range.




The suitable CombiCheck standards and one parameter standards are listed in the WTW catalog or on the Internet.

# Overview of the total system monitoring

Total system monitoring (AQA2) consists of the following parts:

- Configuring the general settings in the AQA2 setup menu.
  - Activate AQA2
  - Select the AQA2 interval unit (Weeks or Measurements)
  - Activate/deactivate the measurement lock for missing or expired AQA2 check. The measurement lock is effective for all methods that were activated for AQA2 monitoring
- Selecting the method to be activated for AQA2
- Configuring the method-specific settings in the AQA2 setup menu.
  - Activate AQA2
  - Specify AQA2 Interval
  - Enter the nominal value, tolerance and designation (standard ID) for the test standard
- Carrying out the AQA2 check. During the check the test is carried out with the standard solution as the sample while the other conditions are the same. The photometer compares the result with the nominal value while taking the tolerance into account.

The steps are described in detail below.



[AQA]	
– AQA2 setu	ıp

AQA2 setup		16.01.15 9:52
Mode	Weeks	
Lock methods Method	Yes	
Method list		

- 1 Select and confirm *Mode*. The *Mode* selection field pops up.
- 2 Select and confirm *Weeks* or *Measurements*.

AQA2 is active. For all methods, the AQA2 intervals are entered either in weeks or number of measurements.

**3** Accept the general settings with *[Apply]*.



When the mode (*Weeks* or *Measurements*) is changed, all AQA2 intervals are reset to the preset values.


**Locking the method** Here you configure whether or not a method will be locked against measurement if there is no valid AQA2 check or the interval for the AQA2 check has expired.

<HOME> [AQA] - AQA2 setup AQA2 setup 16.01.15 9:52 Mode Weeks Lock methods Yes Method ...

- 1 Select and confirm *Lock methods*.
- 2 Select and confirm Yes.

The method lock is enabled.

Each method will be locked if the AQA2 check is invalid or the AQA2 interval has expired.



Method list

Method list

- 1 Select a method (see section 4.5.3).
- 2 Select and confirm AQA2.
- **3** Select and confirm *AQA2 active*. AQA2 is active for this method.

Defining the AQA2 Interval, nominal value and tolerance The AQA2 Interval defines the interval between two AQA2 checks. When an interval has expired, the following consequences become effective:

• Warning and loss of the AQA2 labeling

• Locking of the method against measurement (if activated).

Setting range:

- 1 to 12 weeks (default: 12 weeks) or
- 1 to 10000 measurements (default: 200 measurements)




The unit of the AQA2 interval (Weeks or Measurements) is defined in the line, Mode (see page 145).

AQA2 setup	16.01.15 9:52
Method	3: A6/25
AQA2	AQA2 active
Interval	12 Weeks
Target value	4.00 mg/⊨NH <sub>4</sub> -N
Tolerance	0.50 mg/I NH <sub>4</sub> -N
Standard ID	
Method list	

- 4 Select the Interval and enter the AQA2 interval.
- 5 If necessary, adjust the values for Target value and Tolerance.
- 6 Option: Select Standard ID and enter a designation. The designation is recorded in the AQA2 documentation.

Repeat the steps 1 to 8 if you want to configure further tests for AQA2.



If the check was successful, the AQA2 check function is finished.

**Test record** 

A test record is displayed after the check. It can be printed and stored as a file (in the internal DataB folder or USB memory device at the USB-A

connection, see section 4.12.1).

Sample printout:

```
photoLab 7600 UV-VIS 09130512 1.30-WTW-1.60 Administrator
AQA2 OK
Protocol ID32
Executed by:Administrator
Executed 16.01.2015
Valid until:13.03.2015
Method 6: P6/25 PO4-P
Standard IDCC10 OC557775
Target value0.80 +- 0.08 mg/l
Measured value0.84 mg/l
```



Later you can view the last AQA2 test records for all methods monitored with AQA2 under *AQA2 info*.

### 4.16.4 AQA3/MatrixCheck

The *MatrixCheck* is used to check if the photometric determination is disturbed by other substances present in the sample (sample matrix). The MatrixCheck can be carried out by spiking or diluting:

The photometer enables a simplified MatrixCheck with the aid of the Spectroquant<sup>®</sup> CombiCheck R-2 addition solution. The MatrixCheck can be carried out immediately. The volumes required for the sample and standards are displayed on the screen. The MatrixCheck is then carried out with a single spiking.

For the MatrixCheck with a standard of your own, however, you can enter the number of spikings or dilutions yourself (max. 3).



If the user administration is active, only users of the user group *Administrator* are allowed to change the settings for AQA checks. The AQA check can be carried out by any user.

MatrixCheck by spiking

For the MatrixCheck by spiking, the photometric determination is repeated after a defined amount of analyte has been added to the test sample in the form of standard solutions.

The nominal value for the determination is calculated from the added amount

of analyte, provided that there is no disturbance due to the sample matrix.
After the photometric determination the measured value is compared to the
nominal value and the recovery rate is calculated. A matrix disturbance is
likely if the recovery rate is less than 90 % or more than 110 %.

# MatrixCheck by<br/>dilutingFor the MatrixCheck by dilution, the photometric determination is repeated<br/>after the test sample has been diluted with distilled water.

The nominal value for the determination is calculated from the dilution, provided that there is no disturbance due to the sample matrix. After the photometric determination the measured value is compared to the nominal value and the recovery rate is calculated. A matrix disturbance is likely if the recovery rate is less than 90 % or more than 110 %.

# Practical instructions After evaluating the measured value of the sample the photometer suggests for the MatrixCheck to spike or dilute the sample and standard with suitable volumes. You can change the suggested values of the volumes for the sample and standard. The photometer checks your entries and informs you of errors (e.g. if a nominal value is outside the measuring range of the test). For each spiking or dilution, the relevant nominal concentration value is displayed.

- To be able to reliably recognize matrix effects by <u>spiking</u>, the <u>volume</u> <u>increase</u> after spiking should be <u>small</u>.
- To be able to reliably recognize matrix effects by <u>diluting</u>, the <u>dilution factor</u> should be <u>high</u>.
- You can carry out the MatrixCheck as a series of measurements, consisting of up to three determinations with different spiking volumes or dilutions respectively.
- Prepare all test sample solutions simultaneously at the beginning of the series of measurements.

The MatrixCheck consists of the following parts:


- Configuring settings in the AQA3/MatrixCheck setup menu.
  - Specify the maximum deviation from the nominal value after spiking or diluting (default setting: 10%)
- Carrying out the AQA3 / MatrixCheck

Specifying the maximum deviation from the nominal value

Overview of the

AQA3/MatrixCheck

The assessment of the recovery rate is determined with the maximum deviation from the nominal value. The assessment of the recovery rate is displayed next to the recovery rate after the check has been carried out.



AQA3/MatrixCheck setup	16.01.15 9:52
Maximum difference	10%
Maximum difference	
10.0 %	
	1
	L

- 1 Enter and confirm a numerical value.
  - The setting is active.
- 2 Exit the menu with **<ESC>**.

### Carrying out the MatrixCheck

Concentration		8B	16.01.15 9:52
		45	mg/l
18: C3/25 (44 16 mm	5 nm)		COD 10 - 150 mg/l
Setup	Method list	Citation form	Unit

- 1 Measure the original sample without spiking or diluting it (see section 4.5.1 to 4.5.3).
- 2 The measured value is displayed.
- **3** Open the setting menu with *[Setup]*.
- 4 Select and confirm AQA.
- 5 If necessary, check the settings in the menu, *AQA3/MatrixCheck setup*.
- 6 Select and confirm AQA3/ MatrixCheck.

The display for the MatrixCheck opens up.

MatrixChec	k (Spike)		16.01.15 9:52
Method Sample con	centration	1: C3/25 45 mg/ICOD	
Standard ID Standard co		0 0 mg/ICOD	
Sample [ml]	Standard [ml]	Target value [mg/l]	
10	0	45	
10	0	45	
10	0	45	
Dilute		Delete	Next

If the spiking with the standard values of the CombiCheck R-2 suggested by the photometer would cause the measuring range to be exceeded, the MatrixCheck by diluting is automatically suggested.



The following description shows the proceeding for the MatrixCheck by spiking. To switch to the MatrixCheck by dilution, use the *[Dilute]* function key. The proceeding is similar there, but the entry of the Standard ID and Standard concentration is not applicable.

MatrixChec	k (Spike)		16.01.15 9:52
Method Sample con	centration	1: C3/25 45 mg/ICOD	
Standard ID Standard co	) oncentration	COD 1500 400 mg/ICOD	
Sample [ml]	Standard [ml]	Target value [mg/l]	
10	0.5	62	
10	1	77	
10	1.5	91	
Dilute		Delete	Next

7 In the *Standard ID* entry field, select the simplified MatrixCheck with the CombiCheck standard solution or enter a designation for another standard solution used.

If the CombiCheck is selected, no more entries are required (continue with step 10).

8 Enter the concentration of the used standard solution in the *Standard concentration* entry field.

Specifying the series of measurements:

- 9 Enter the volumes of sample and standard of the individual test sample solutions in the columns, *Sample [ml]* and *Standard [ml]*. The nominal value is calculated after each entry.
  - You can delete a measurement from the series of measurements with [Delete].

Note that all nominal values have to be within the measuring range of the test.

**10** Using *[Next]*, accept all entries on the page and switch to the next page. The entries are checked by the photometer.

The photometer is ready to carry out the series of measurements.

# Carrying out the series of measurements:

According to the program, the samples are measured top down. You can, however, select the samples yourself and thus change the order with  $< > < \forall >$ .

**11** Use *[Measurement]* to proceed to the measurement of the (first) sample.

The measurement display appears.

**12** Insert the cell with the respective sample.

The sample is measured.

MatrixCh	eck (Sp	pike)			16.01.15 9:52
Method Sample c	oncent	ration	1: C3/2 45 mg/l		
Sample [ml]	Stanc [ml]	lard	Target value [mg/l]	non [mg	ninal I/I]
10	0.5		62	58	
10 10	<b>0.5</b> 1		62 77	58	
	0.5 1 1.5			58	
10	1 1.5	Measur	77 91	58	Complete

MatrixCheck		16.01.15	9:52
Method Sample concentration Sample Standard	1: C3/25 45 mg/ICOD 10 ml 0.5 ml		
To start measurer insert cell or press	nent, s <start enter=""></start>		
16 mm			
Back			

MatrixCh	leck			6.01.15 9:52
Method Sample o	oncentratio	1: C3/3 n 45 mg		
Sample [ml]	Standard [ml]	Target value [mg/l]	nomin [mg/l]	
10	0.5	62	58	94 % 🗸
10	1	77		
10	1.5	91		
Bac	k Mea	sureme		Complete

After the measurement, the recovery rate is displayed in the right table column.

The assessment of the recovery rate is displayed next to the recovery rate ( $\checkmark$  or  $\checkmark$ ).

The criteria for the assessment are determined in the menu, *AQA3/MatrixCheck setup / Maximum difference*.

- **13** If necessary, repeat the steps 11 and 12 for the remaining samples.
- 14 Use [Complete] to complete the MatrixCheck.

The Save dialog box pops up.

 15 If necessary, change the storage location with [Location]: Internal DataB folder.
 Exchange folder in the instrument or USB memory: USB memory device connected at

the USB-A connection.

- **16** If necessary, change the file name.
- 17 Save the file with <START.ENTER>.

The display returns to the measured value display of the original sample without spiking / dilution.

The [MC] status indicator is displayed. A MatrixCheck was carried out for this measured value.

**Test record** The result of the MatrixCheck is displayed in a test record. You can print this record and save it as a file.

COD

10 - 150 mg/l

Unit

16.01.15 9:52

[MC]

45

Method list Citation form

To save the file on the photometer, select *Internal DataB folder* as the location. To save the file to an external USB memory device connected to the USB-A connection, select *USB memory* as the location (see section 4.12.1).

Sample printout:

Concentration

18: C3/25 (445 nm)

Setup

16 m m

```
photoLab 7600 UV-VIS 09130512 1.30-WTW-1.60 Administrator
MatrixCheckOK
Protocol ID7
Method 1: C3/25 COD
Sample concentration45 mg/lCOD
Standard IDCOD 1500
Standard concentration400 mg/l COD
Sample
       StandardTarget valueActual value
ml
        ml
                mg/lmg/l
10
        0.5
                625894% OK
10
        1
                 777192% OK
10
        1.5
                 917785% OK
```

### 4.17 User management

The functions of the user management are only available for users of the user group, *Administrator*.

An administrator can

- activate or deactivate the user management for the meter
- create, change or delete individual user accounts.

### 4.17.1 User levels and user rights

The photoLab<sup>®</sup> 7600 UV-VIS allows the management of up to 100 users. Every user is member of a user group with defined user rights.

**User groups** There are three hierarchical user groups:

- Administrator (top level)
- User (user account registered by the administrator)
- Guest (user without user account)

Administrators and users log in to the photometer with their user name and password. Guests can optionally enter a name for their login. Thus, documented measured values can later be assigned to the user.

User rights in detail	Action	Administrator	User	Guest
	Select methods	✓	1	✓
	Carry out measurements	1	1	✓
	Store measurement data	1	1	<ul> <li>✓</li> </ul>
	Check photometer (AQA1)	1	1	$\Diamond$
	Check total system (AQA2)	1	1	$\Diamond$
	AQA1 measured value labeling	1	1	<ul> <li>✓</li> </ul>
	AQA2 measured value labeling	1	1	$\Diamond$
	Edit user-defined methods	✓	1	Ø
	Exchanging methods / profiles	✓	Ø	Ø
	Change AQA settings	✓	Ø	Ø
	Clear the memory	✓	Ø	Ø
	Set the date and time	<ul> <li>✓</li> </ul>	0	$\Diamond$
	Administrate users	<ul> <li>✓</li> </ul>	0	$\Diamond$
	Reset photometer settings	✓	Ø	$\Diamond$
	Carry out software update		0	Ø



<HOME>

[Conoral actual

You can also switch off the user management and reactivate it as necessary. To do so, you need administrator rights. If the user management is switched off, the user name and password do not have to be entered. Each user has full rights.

### 4.17.2 Activating or deactivating the User management function

Each user can activate the user management function. If the function is deactivated, each user has administrator rights.

Only members of the user group, administrator can deactivate the user management function.

If the function is active, each user has to log in to the photometer. After the login, the user has certain rights depending on the user group.

### Activating the user management function

User m	nanagement	16.01.15 9:5
-		
	User management not active	
	User management not active Activate user management?	
	Activate user management?	

Deactivating the user management function

<HOME> [General setup] |- User management |- [Setup] |- Deactivate user management **1** Select and confirm *Yes*.

The user management function is active.

Activating the user management creates an administrator user account. The user name is "Administrator". The preset password is "admin". Change this password as soon as possible.

The user management function is inactive.

Each user has administrator rights.



If the user management is deactivated by a user of the *Administrator* user group, all user accounts that were set up are lost. The password for the administrator is reset to "admin".

#### 4.17.3 Creating, changing or deleting a user account

When the user management function is active, a user with administrator rights can administrate user accounts.

Creating a user<br/>accountDuring the creation of a user account, the Name, whether or not the user<br/>belongs to a User group and the Password are defined.

	<home> [General setup]  - User management  - [Add]</home>			
	User management Name User group	16.01.15 9:52		The input field for the new user name pops up.
	Administrator Administrator Admini2 Enter user name A_		1	Enter the user name ( <b><a9></a9></b> ) and confirm.
				The selection field for the user group ( <i>Administrator / User</i> ) pops up.
			2	Select and confirm the user group.
	Setup Add Delete	Change		The input field for the password pops up.
			3	Enter the password ( <b><a9< b="">&gt;) and confirm.</a9<></b>
				The user account is created and appears in the list of user accounts.
Editing a user account	When a user account is ch changed.	nanged, th	ne Us	<i>er group</i> and <i>Password</i> can be
	<home> [General setup] - User management</home>			

User manag	ement		16.01.15 9:52
Name	Usergroup		
Administrato	r Administrator		
Admin2	Administrator		
User	group		
User			
Admi	nistrator		
		r <u></u>	
Setup	Add	Delete	Change

- 1 Select a user account.
- 2 Press [Change] to edit the user account.

The selection field for the user group (*Administrator / User*) pops up.

**3** If necessary, select and confirm another user group.

The input field for the password pops up.

4 If necessary, enter (**<A...9>**) and confirm another password.

The user account is changed and appears in the list of user accounts.

### Deleting a user account



- 1 Select a user account.
- 2 Delete the user account with [Delete].

A security prompt appears: *Confirm deletion ?* 

**3** Confirm the security prompt. The user account is deleted.

### 4.17.4 Login with active user management

To be able to always assign measurement data to a user, the administrator can activate the user management function. After doing so, the photometer can only be operated after login with a user name. Depending on the authorization class (administrator, user, guest), important settings are released for changing or locked.



The user management function is not active in the delivery condition of the photoLab<sup>®</sup> 7600 UV-VIS. Every user can carry out all functions.

Activating the user management creates an administrator user account. The user name is "Administrator". The preset password is "admin". Change this password as soon as possible.

Make sure to use the correct spelling (upper and lower case) of user name and password for the login.

After logging in to the *Administrator* group with a user name, you can create further users or administrators or switch off the user management function.

The *Login* window with the *Enter user name* prompt appears after the meter has been switched on and after a user has logged off. In the following example, a user will log in with the user name, "Administrator".

Login	16.01.15 9:52
Enter user name	
Administrator	

The photometer is switched on. The *Login* dialog is displayed.

1 Enter the user name (**<A...9>**) and confirm.

The input field for the password pops up.

If the user name is not known (or incorrectly spelled) it is possible to log in without a password as a guest with restricted rights (see section 4.17.1).

Login				1	6.01.15 9:52
		ssword			
ac	dmin				
Home (Adn	ninistrat	or)			05/04/15 12:56
	Conc	entration		OptRF	
25.1% 0.600		rbance ransmission		Spectr	'um
λ.	Spec Mult	ial i wavelengths	T	Kineti	cs

2 Enter the password (**<A...9>**) and confirm.

If the password is written correctly (note upper and lower case), the main menu *Home* opens up. The user name that was entered is displayed.


### 4.17.5 Changing the password

AQA Info

General setup

The administrator sets up user accounts and assigns a password to each user account.

As soon as any user has successfully logged in with the password, they can change the password for their user accounts themselves.



- 1 Enter and confirm the old password.
- 2 Enter and confirm the new password.

The password is changed.

### 4.18 Reset

You can reset (initialize) the measurement settings or all settings.



The *Reset* function is available for users of the user group Administrator only.

You have the following options of resetting the photometer settings:

<ul> <li>Reset configuration</li> </ul>	All settings except for the measurement data memory, user- defined methods and measured blank values are deleted.
<ul> <li>Delivery condition</li> </ul>	All settings (including measurement data memory and user-defined methods) are deleted and the photometer is reset to the delivery condition.
<ul> <li>Internal DataB folder</li> </ul>	The measurement data memory is erased. All other settings are retained.
	Save your measurement data (e.g. to a USB memory device) before erasing the internal data memory of the photometer.

### 

The menu where to select the reset type (*Delivery condition* / *Reset configuration*) is displayed.

1 Select and confirm the reset type. The reset is carried out.

### 4.19 Photometer information ([Info])

The following photometer information is displayed:

- Photometer designation
- Version number of the meter software/method data
- Hardware version
- Series number of the meter
- Registered user
- Hardware status (for service purposes)
- Memory status

### <HOME> [Info]

Info	16.01.15 9:52
Model designation:	photoLab <sup>®</sup> 7600 UV-VIS
Serial number:	15150001
Software/methods version:	2.18-WTW-2.20
IQ-LabLink version:	1.20
Build:	25.02.15 11:57
Hardware version:	0-1-0-22-50
Hardware status:	FF 0000000
Lamp counter	12
System test	11
Filter test	11
Lamp test	11
Wavelength calibration	
Free internal memory space	5046 KB
Registered user	1

The meter information and result of the self-test are displayed and can be printed.

### 4.20 Lamp counter

The photometer counts the operating hours of the lamp. The information on the operating hours of the lamp is given in the *Info* menu.

The number quoted there corresponds to the number of flashes.

### 4.21 Software and methods update

The software and method update is used to continuously update your photometer.



If the user administration is activated, only users of the user group *Administrator* may carry out any software and method updates.

The update comprises

- the newest firmware (meter software)
- new or changed method data



User-defined data (such as settings, user-defined methods or measured data) are not changed by a software and methods update.

The current software version is available on the Internet at http://www.WTW.com.

The software can be transmitted to the photometer as follows:

- by means of a USB memory device as a temporary storage (section 4.21.1).
- by means of a USB connection between a PC and the photometer (section 4.21.2).

### 4.21.1 Software- and method update using a USB memory device

Store the new software required for the update on the USB memory and connect it to the photometer.

### Execution

- 1 Connect the USB memory device to the PC.
- 2 Unpack the contents of the downloaded exe or zip file with the entire folder structure to the main directory (top level) of the USB memory.



Make sure the folder structure of the files is retained during the unpacking process.

If you use a program such as WinZip for unpacking, the option "Nutze Ordnernamen" or "Use Folder Names" must be set. Details are given in the documentation of your unpack program.

The USB memory must have the folder "Update" on the top level. The Update folder comprises several subfolders.

The following steps are carried out at the photometer.

- **3** Connect the USB memory device to the photometer.
- 4 Switch on the photometer if necessary.

# <HOME>

#### [General setup]

- Software/methods update

Software/methods update	16.01.15	9:52
Software/methods update		
Select source of update data:		
USB memory device		
PC		
Cancel		
	1	

5 Using <▲><▼>, select USB memory device as the source and press <START.ENTER>.

The update process takes approx. five minutes. Subsequently, the photometer switches itself off and then on again.



If the update cannot be carried out, an error message appears on the display.

Check whether the "Update" folder with its subfolders is stored on the USB memory device (top level).

If on the photometer there is not enough free memory capacity for the update, you can create memory capacity by erasing measurement data. Save your data to a USB memory device before erasing them on the photometer.

### 4.21.2 Software- and method update via PC

The software- and method update via PC is done with the PC software photoLab<sup>®</sup> Data *spectral* (for details, see operating manual of the PC

software photoLab<sup>®</sup> Data spectral).


# 5 Maintenance and cleaning

5.1 Exchanging the buffer battery

### CAUTION

There is a risk of explosion if unsuitable batteries are used. Only use leakproof alkaline manganese batteries.

If you leave the photometer switched on during the exchange or insert the new batteries within a minute after taking out the old ones, the date and time are retained in the photometer.



- **1** Turn the photometer upside down and place it on a soft surface.
- **2** Open the lid of the battery compartment (1).
- 3 Remove the old batteries from the battery compartment.
- Insert the four new batteries in the battery compartment. Make sure that the poles of the batteries are in the correct position.
   The ± signs on the batteries must correspond to the ± signs in the battery compartment.
- 5 Close the lid of the battery compartment.

**Battery service life** The power consumption of the clock is very low. The lifetime of high quality batteries is at least 5 years.

Disposal of<br/>batteriesDispose of the batteries at a suitable facility according to local legal require-<br/>ments. It is illegal to dispose of the batteries with household refuse.

Within the European Union, the batteries are removed at a specialized treatment center at the instrument's end of life. The instruments are taken to one of those specialized treatment centers via the recycling system set up for this purpose.

### 5.2 Cleaning

Especially after a cell has broken or after a reagent accident, the photometer should immediately be cleaned (see also section 6.1 ACTIONS IN THE CASE OF A BROKEN CELL).

### 5.2.1 Cleaning the enclosure

### NOTE

The housing components are made out of synthetic materials (ABS, PMMA and PC). Thus, avoid contact with acetone, ethyl alcohol and similar detergents that contain solvents. Remove any splashes immediately.

Clean the photometer enclosure as follows:

- If the housing surface is dirty, wipe it with a soft cloth and mild soapy water.
- Remove any chemicals splashes as soon as possible.
- For disinfection, you can use isopropanol for cleaning for a short time.

### 5.2.2 Cleaning the cell shaft

### NOTE

The housing components are made out of synthetic materials (ABS, PMMA and PC). Thus, avoid contact with acetone, ethyl alcohol and similar detergents that contain solvents. Remove any splashes immediately.



If a cell has broken, the cell shaft has to be cleaned immediately. To do so, proceed as described in section 6.1.

Normally, it is not required to clean the cell shaft routinely. Remove dust and slight contamination with a moist, lint free cloth. Use isopropanol <u>briefly</u> to remove persistent coatings (e.g. reagent remains). Especially clean the bottom parts of the lateral surfaces of the rectangular cell shaft where the light barriers for the automatic cell recognition are located.

### 5.2.3 Cleaning the detector lens


Normally, it is not required to clean the detector lens routinely. Cleaning the detector lens can be necessary in the following cases:

- If the lens is visibly smudged, e.g. after a cell has broken or after a reagent accident (see also section 6.1 ACTIONS IN THE CASE OF A BROKEN CELL).
- If, due to environmental influences or reagent contamination, the photometer displays the message, *Wavelength calibration* during the self-test after being switched on (see section 6.2)



If the lens is often smudged (error, *Wavelength calibration* during the self-test), check whether the correct operating conditions are observed. Follow the details in section 3.2 for this purpose.

Proceed as follows to clean the detector lens:



The detector lens is on the front left side of the rectangular cell shaft (pos. 1).

- **1** Switch off the photometer.
- 2 Cut off one end of a customary cotton swab (approx. 2 cm).
- 3 Grasp the cut-off end with the tip of a pair of tweezers or small pliers. Clean the lens with the dry head of the cotton swab. To do so, move the head from the center of the lens outward in circles. If there are persistent coatings, moisten the cotton swab with a little deionized water or isopropanol.



After recommissioning, carry out the photometer monitoring for all measurements (see section 4.16.2).

# 6 What to do if...

6.1 Actions in the case of a broken cell



### WARNING

Cells can contain dangerous substances. If the contents are released, follow the safety instructions of the package insert. If necessary, take corresponding protective measures (protective goggles, protective gloves etc.).



### CAUTION

Do not turn the photometer upside down or laterally to remove the liquid!

When doing so, the liquid could come into contact with electronic components and damage the photometer.

The photometer has a drain device through which the contents of a broken cell can drain off without causing any damage.

# Proceeding after a cell has broken

- 1 Switch off the photometer and disconnect it from the power supply.
- 2 Let the liquid drain off into a suitable container and dispose of it properly according to the instructions of the reagent package.
- **3** Carefully remove all broken glass, e.g. with tweezers.
- 4 Carefully clean the cell shaft using a moist, lint-free cloth. If there are persistent coatings, use isopropanol <u>for a short time</u>. Especially clean the bottom parts of the lateral surfaces of the rectangular cell shaft where the light barriers for the automatic cell recognition are located.
- 5 Let the cell shaft dry.



After recommissioning, carry out the photometer monitoring for all measurements (see section 4.16.2).

If, after recommissioning, an error occurs during the wavelength calibration,

**Cleaning the** 

detector lens

the detector lens is probably smudged. In such a case, clean the lens as follows:

The detector lens is on the front left side of the rectangular cell shaft (pos. 1).

- Switch off the photometer.
- 2 Cut off one end of a customary cotton swab (approx. 2 cm).
- 3 Grasp the cut-off end with the tip of a pair of tweezers or small pliers. Clean the lens with the dry head of the cotton swab. To do so, move the head from the center of the lens outward in circles. If there are persistent coatings, moisten the cotton swab with a little deionized water or isopropanol.

### 6.2 Error causes and remedies

Meter does not react	Cause	Remedy
to keystroke	<ul> <li>Operating condition undefined or EMC load unallowed</li> </ul>	<ul> <li>Processor reset:</li> <li>Press the <b><on off=""></on></b> and</li> <li><b><esc></esc></b> key simultaneously.</li> </ul>
Audio signal on	Cause	Remedy
keystroke	<ul> <li>The key does not have any func- tion in the current operating state</li> </ul>	<ul> <li>Press the appropriate key</li> </ul>
Measuring range	Cause	Remedy
undercut or exceeded	<ul> <li>Measuring range or method not suitable</li> </ul>	<ul> <li>Select method with suitable mea- suring range</li> </ul>
		<ul> <li>Dilute the sample</li> </ul>



In *Concentration* mode you can display the current absorbance value as an additional information (*[Setup]/Display absorbance*, see also section 4.5.6).

1

Self-test does not start.	Cause	Remedy
The instrument displays Please remove cell	<ul> <li>A cell is inserted in one of the cell shafts</li> </ul>	<ul> <li>Remove the cell</li> <li>Then press the</li> <li><start-enter> key</start-enter></li> </ul>
	<ul> <li>A foreign object is inserted in one of the cell shafts</li> </ul>	<ul> <li>Remove foreign object</li> <li>Then press the</li> <li><start·enter> key</start·enter></li> </ul>
	<ul> <li>Occasionally, the instrument carries out an automatic readjustment for the rectangular cell recognition. The informative message <i>Please remove cell</i> is displayed even when no cell is inserted.</li> </ul>	<ul> <li>Press the <b>START·ENTER</b></li> <li>key.</li> </ul>
	<ul> <li>The cell shaft is contaminated</li> </ul>	<ul> <li>Clean the cell shaft (see section 5.2.2 and section 6.1)</li> <li>Restart the instrument</li> <li>If necessary, confirm the <i>Please remove cell</i> message with</li> </ul>
	<ul> <li>Instrument defective</li> </ul>	<pre> START·ENTER&gt;.  Please contact the service department. </pre>

Obviously incorrect measured values	Cause	Remedy
measured values	<ul> <li>Cell contaminated</li> </ul>	- Clean the cell
	<ul> <li>Dilution set incorrectly</li> </ul>	<ul> <li>Set the dilution</li> </ul>
	<ul> <li>Selected method not suitable</li> </ul>	<ul> <li>Select different method</li> </ul>
	- Zero measurement incorrect	<ul> <li>Perform zero measurement</li> </ul>
	<ul> <li>Blank value incorrect</li> </ul>	- Remeasure the blank value

Fluctuating measured values	Cause	Remedy
measured values	<ul> <li>Cell shaft cover open</li> </ul>	<ul> <li>Close the cell shaft cover</li> </ul>

Self test failed.	Cause	Remedy
	- System test: Instrument defective	<ul> <li>Please contact the service department.</li> </ul>
	- Filter test: Instrument defective	<ul> <li>Please contact the service department.</li> </ul>
	- Wavelength calibration:	
	<ul> <li>Foreign particle in the cell shaft</li> </ul>	<ul> <li>Remove foreign object</li> </ul>
	<ul> <li>Lens smudged</li> </ul>	<ul> <li>Clean the lens (see section 5.2.3 or section 6.1).</li> <li>If this happens repeatedly, check the operating conditions (see section 3.2)</li> </ul>
	<ul> <li>Instrument defective</li> </ul>	<ul> <li>Please contact the service department.</li> </ul>
Connected printer	Cause	Remedy
does not print	<ul> <li>Printer not suitable</li> </ul>	<ul> <li>Connect a printer that can inter- pret the printer control language PCL-3</li> </ul>
Data transmission	Cause	Remedy
to USB memory device does not work	<ul> <li>Connected USB memory device was not recognized</li> <li>The USB memory device has been formatted to a file system which is not supported, e. g. NTFS</li> </ul>	<ul> <li>Use other USB memory device</li> <li>Format the USB memory device to the FAT 32 file system</li> </ul>

# 7 Technical data

### 7.1 Measurement characteristics

Measuring principle Spectrophotometer with reference beam tech
--

Light source	Lamp type	Xenon flashlamp
	Average lifetime	$5 \times 10^8$ flashes, corresponding to at least 13000 h in permanent operation
Monochromator	Туре	Grating monochromator with step motor
	Wavelength range	190 - 1100 nm
	Max. scan speed	approx. 1000 nm/min
	Wavelengths calibration	Automatic
	Accuracy	± 1 nm
	Reproducibility	$\pm$ 0.5 $\pm$ nm (checkable, e.g. with holmium oxide filter)
	Resolution	1 nm
	Spectral band width	4 nm
Photometric mea-	Light sensor	Photo diode
surement	Measuring range	A = -3.300 to A = +3.300
	Linearity	< 1 % for A ≤ 2.000 in the range 340 900 nm
	Accuracy*	± 0.003 A for A < 0.600
		$\pm$ 0.5 % of the reading for 0.600 ${\leq}A{\leq}$ 2.000
	Reproducibility *	± 0.002 at A = 1.000
	Resolution	ΔA = 0.001
	Scattered light	< 0.05 % transmission at 340 and 408 nm

\* in the range 200 nm ... 1000 nm

Usable cells	Round cells	<ul> <li>Outer diameter: 16 mm</li> <li>Inner diameter: 13.6 mm</li> <li>Flat cell bottom</li> </ul>
	Rectangular cells*	<ul><li>Path length: 10 mm, 20 mm and 50 mm</li><li>Maximum width: 12.6 mm</li></ul>
	Minimum filling level	20 mm
	Minimum filling volume	Round cell 16 mm: 4 ml
		Rectangular cell, 10 mm: 2 ml
		Rectangular cell, 20 mm: 4 ml
		Rectangular cell, 50 mm: 10 ml
	Cell recognition	Automatic for most types

Suitable are rectangular glass cells with opaque lateral surfaces (see section 8.1). Plastic cells with clear or serrated lateral surfaces are not reliably recognized by the automatic cell recognition.

**Measuring modes** 

#### Concentration

- Measurement with permanently programmed methods, adapted to the WTW test set program
- Automatic method selection if test sets with barcodes are used
- Program support for the creation of additional user-defined methods (max. 100)
- Citation forms and units method dependent
- Display of the absorbance value can be added
- Method data update possible via Internet
- Absorbance / % Transmission
  - Measurement against own reference absorbance value possible
- Multi wavelengths
  - Freely definable calculations from absorbance values at up to 10 different wavelengths
  - Calculations can be stored as methods (max. 499)

Measuring modes	<ul> <li>Spectrum</li> </ul>
-----------------	------------------------------

- Absorbance or % transmission mode
- Limits freely selectable within the wavelength range
- Increment: 1 nm
- Recording duration for the complete wavelength range: < 7 min
- Settings can be stored as profiles (max. 20)
- Evaluation functions: Cursor scanning, zoom, min./max. recognition, peak area determination, derivation, smoothing, multiplication by constants, addition of constants, addition and subtraction of spectra, formation of the quotient of two spectra
- Kinetics
  - Absorbance or % transmission mode
  - Minimal adjustable scan interval: 1 s (if the absorbance of the test sample is high, the scan interval is extended due to the longer duration of the individual measurements)
  - Settings can be stored as profiles (max. 20)
  - Evaluation functions: Cursor scanning, zoom, min./max. determination, slope calculation (for an interval or total), enzymatic activity
- OptRF measurement
  - Measurement with permanently programmed methods for which no test sets are required
  - Citation forms method dependent

### 7.2 Measured value documentation and quality assurance

Memory for mea- sured values	Memory capacity	<ul> <li>5000 single measured values from the measuring modes, concentration, absor- bance / % transmission and multi wave- lengths</li> </ul>
		<ul> <li>40 MByte internal memory, sufficient for approx. 500 spectra and 400 kinetic curves (sample values based on the following assumptions: All spectra over a wave- length range of 600 nm and all kinetic curves with 150 single values each)</li> </ul>
	Output options	USB memory device, printer, PC
	File formats	ASCII, *.csv

Monitoring func-	AQA1	Check of the photometer
tions	AQA2	Check of the total system
	AQA3	Check of the sample matrix
User management	Can be switched off	yes
User management	Can be switched off User accounts	yes 3 hierarchical levels (administrator, user, guest)

### 7.3 General meter data

Dimensions	404 x 197 x 314 mm (width x height x depth)		
Weight	approx. 4.5 kg (without plug-in power supply)		
Housing type of protection	IP 30		
Electrical protec- tive class	III		
Test certificates	CE		
Allowed environmental	Temperature	Operation: +10 °C to + 35 °C (41 °F to 95 °F)	
conditions		Storage: -25 °C to + 65 °C (-13 °F to 268 °F)	
	Humidity	Yearly mean: $\leq 75 \%$ 30 days/year:95 %	
		Other days: 85%	
	Climatic class	2	
Power	Power pack	Type: EDACPOWER EA1036R Input: 100 - 240 V ~ / 50 - 60 Hz / 1 A	
supply			
		Output: 12 V = / 3 A	
		Length of the connection cable to the photometer: 2.0 m (In compliance with Eco-design directive 2009/125/EG, EuP step 2)	

Guidelines and norms used	EMC	EC directive 2004/108/EC EN 61326-1	
		<ul> <li>Interference emission: Class B</li> </ul>	
		<ul> <li>Interference immunity: IEC 61000-4-3</li> </ul>	
		- Tolerance extension: 0.008 E	
		FCC Class A	
	Meter safety	EC directive 2006/95/EC	
		EN 61010-1	
	Climatic class	VDI/VDE 3540	
	IP protection class	EN 60529	
Communication	Ethernet	RJ45	
interfaces	USB	<ul> <li>1 x USB-A (for printer, USB memory devices, keyboard or bar code reader)</li> </ul>	
		- 1 x USB-B (for PC)	
Other features	<ul> <li>Drain for spilled cell contents</li> </ul>		

• Photometer software update and method data update possible via Internet

Available languages

- German
- English
- Français
- Español
- Italiano
- Bulgarian/Български
- Česko
- Simplified Chinese/ 中文
- Traditional Chinese/ 繁體中文
- Greek/Ελληνικά
- Indonesian/Indonesia
- Japanese/日本語
- Magyar
- Malay/Melayu
- Norsk
- Polski
- Portuguése
- Russian/Русский
- Slovenščina
- Thai/ ภาษาไทย
- Turkish/Turkce
- Dansk
- Română
- Nederlands

# 8 Accessories and options

### 8.1 Accessories

Cells for the WTW test set program	Description	Model	Order no.
	25 empty round cells (16 mm)	RK 14/25	250 621
	Rectangular cell, 10 mm	REK 10	250 605
	Rectangular cell, 20 mm	REK 20	250 600
	Rectangular cell, 50 mm	REK 50	250 614
	Rectangular cell, quartz, 10 mm	REK 10 Quarz	250 606
	Rectangular cell, quartz, 20 mm	REK 20 Quarz	250 601
	Rectangular cell, quartz, 50 mm	REK 50 Quarz	250 615
	100 PMMA disposable rectangular cells, 10 mm	REC 10 PMMA	250 607
	Half micro cells, optical special glass, 50 mm	REC 50/2	250 616
	Stand aid for 10 mm plastic cells	PL6-10 SIC	250 213
		1	1
Case and cable for mobile use	Description	Model	Order no.
mobile use	Case for photoLab <sup>®</sup> 7xxx	FC spectral 6000	250 212
PC software	Description	Model	Order no.

PC software	Description	Model	Order no.
	PC software package for convenient data management	photoLab <sup>®</sup> Data spectral	902 761

Application packages	Description	Model	Order no.
	Brewing application package (Ger- man/English)	PL6-BREW	250 214

#### 8.2 **Test equipment**

Test equipment	Description	Model	Order no.
	Test equipment for AQA1	PhotoCheck 14693	250 490
	Test equipment for pipette volume	PipeCheck 14692	250 498



Standard solutions for the WTW test set program are listed in the WTW catalog or on the Internet.

#### 8.3 **Optional equipment**

The following optional extensions are available in specialist shops:

- USB barcode reader (hand scanner)
- USB PC keyboard
- 12 V power supply unit (e.g. portable power source or car battery) for mobile operation of the photometer
- 12 V vehicle connection cable for the mobile operation of the photometer by means of a 12 V power source.

Suitable are vehicle connection cables with the following characteristics:

Voltage	12 V
Amperage	8 A
Barrel connector	2.5 x 5.5mm
Positive pole at the inner contact	

#### **Connection cable:** 8.4

You can connect a PC to the photoLab<sup>®</sup> 7600 UV-VIS in one of the following PC ways:

Description	Order no.
<ul> <li>Cable with USB-B and USB-A plug</li> </ul>	Specialist shops

USB printer	You can connect a USB printer to the photoLab <sup>®</sup> 7600 UV-VIS:
-------------	---

Description	Order no.
<ul> <li>Cable with USB-B and USB-A plug</li> </ul>	Specialist shops

# 9 Appendix

g	.1 Glossary
Absorbance	Logarithmic dimension for the absorption of the sample; negative decadic logarithm of the transmission.
Analysis instructions	The exact proceeding to carry out the detection procedure is described in the analysis instructions.
AQA	Analytical Quality Assurance.
AQA labeling	In the documentation, measured values are given an AQA labeling (AQA1 or AQA2), depending on whether or not the measurement was carried out with AQA and with which AQA level.
AQA1	1st step of the analytical quality assurance: Monitoring of the instrument.
AQA2	2nd step of the analytical quality assurance: Monitoring of the total system.
AQA3	3rd step of the analytical quality assurance: Check of whether the photometric determination is disturbed by other sample ingredients (sample matrix). The MatrixCheck can be carried out by spiking or diluting:
AutoSelector	Plastic cylinder with bar code. It is inserted in the round cell shaft and transmits the code for a reagent test set to the photometer.
Barcode	Optical code (black and white bars) of the method that can be read by light barriers in the photometer.
Baseline	Reference value for the spectrum of reference absorbances or reference transmissions.
Cell	Vessel to take a liquid sample for measurement in a photometer. The cell material (mostly glass) must have certain optical features to be suitable for photometry.
Citation forms	Different display formats that can be derived from each other of the measured value for a concentration.
	The method for the determination of phosphate, e.g. supplies a mea- sured value for phosphorous P. This measured value can alterna- tively be given in the citation forms PO4, PO4-P or P2O5.
CombiCheck	Multiparameter standards used to check the total system for a method.
Concentration	Mass or amount of a dissolved substance per volume, e.g. in g/L or mol/L.

Correlation coefficient	Specifies the extent of the linear relationship of value pairs when determining the zero point and slope for a user-defined method.
Detection procedure	The detection procedure designates the general principle of how a sample is brought into a form suitable for measurement.
	Different methods can be based on the same detection procedure.
Kinetics	Measurement over a period of time.
MatrixCheck	see AQA3.
Measured value	The measured value is the special value of a measured parameter to be determined. It is given as a combination of the numerical value and unit (e. g. 3 m; 0.5 s; 5.2 A; 373.15 K).
Method	A method comprises a chemical detection procedure and special method data (calibration line) required to evaluate the measurement results.
	How to carry out the method up to measuring with the photometer is described in the analysis instructions.
	The photoLab $^{\ensuremath{\mathbb{R}}}$ 7600 UV-VIS contains a database with methods. Furthermore, user-defined methods can be entered in the database as well.
PhotoCheck standard	Stable color solution with defined absorbance values for the check of the photometer.
Reagent blank value	The evaluation of the photometric measurement always refers to the comparison value of a test sample without the substance to be determined (reagent blank value). Thus the influence of the basic absorbance of the reagents on photometric measurement is compensated for.
Recovery	The recovery rate is the found measured value divided by the default value (percentage). Example: Default value 20 mg/l; Found 19.7 mg/l => recovery 0.985 or recovery rate 98.5%.
Reference absorbance	With the reference absorbance, the basic absorbance stored in the photometer can be replaced by a measurement of your own.
Reset	Restoring the original condition of all settings of a measuring system.
Sample blank value	Absorbance value of the original test sample, handled according to prescription but without color reagent.
Spectrum	Distribution of the intensity, transmission or absorbance depending on the wavelength.
Standard	Sample with a defined concentration of the analyte to be determined.

Test sample	Designation of the test sample ready to be measured. Normally, a test sample is made by processing the original sample. The test sample and original sample are identical if the test sample was not processed.
Test set (test)	A test set contains all reagents that are required for the photometric determination of the sample according to the analysis instructions.
Transmission	The part of the light that goes through the sample.
Turbidity	Light attenuation caused by diffuse scattering at undissolved sub- stances.
Zero adjustment	Adjusting a photometer with a water-filled cell.

### 9.2 List of trademarks

Trademark	Owner
Certipur <sup>®</sup>	Merck KGaA
Microsoft <sup>®</sup>	Microsoft Corporation
Spectroquant <sup>®</sup>	Merck KGaA
Windows <sup>®</sup>	Microsoft Corporation

## 9.3 Index

# Α

Absorbance / % Transmission, measuring	71
Accessories	
Analysis timer	
Analytical quality assurance (AQA)	
AQA	

## В

21
49
47

# С

Cell breakage	
Cleaning	
Commissioning	
Concentration measurement	
Connections	7
Connectors	7
Copying files	
Copyright	2

## D

Dataset	
Date/time	
Disinfection	

# G

183
l

## I

Initial commissioning	15
Initialize	
Instrument settings	
IQ LabLink	69

# Κ

Keypad8
Kinetics105

### Μ

Measure diluted sample	45, 92
Measurement data memory	
Measurement dataset	
Meter information	

Method	
Methods update	
Multi-wavelengths methods	75

# 0

Operating elements	7
Operating principles	
Overview	7

# Ρ

Print	
Printer	
Profile (kinetic)	
Profile (spectrum)	
(	-

# R

Reagent blank value	
Measuring the	50
Reference absorbance	
Reset	

# S

Safety	11
Sample blank value	47
Scope of delivery	14
Self-test	
Software update	
Software version number	
Standard adjustment	52
Storing in memory	120
Switching on	23
System info	
System management	

## т

Technical data	
Timer	
Turbidity correction	58

# U

Update	
USB memory device	
User calibration	
User-defined methods	
Concentration	
Multi wavelengths	75
-	

# Ζ

Zero adjustment		32	2
-----------------	--	----	---

# Xylem |'zīləm|

1) The tissue in plants that brings water upward from the roots; 2) a leading global water technology company.

We're 12,500 people unified in a common purpose: creating innovative solutions to meet our world's water needs. Developing new technologies that will improve the way water is used, conserved, and re-used in the future is central to our work. We move, treat, analyze, and return water to the environment, and we help people use water efficiently, in their homes, buildings, factories and farms. In more than 150 countries, we have strong, long-standing relationships with customers who know us for our powerful combination of leading product brands and applications expertise, backed by a legacy of innovation.



a xylem brand

Distributed by:



Tallaght Business Park Whitestown, Dublin 24, Ireland D24 RFK3

Tel: (01) 4523432 Fax: (01) 4523967 Quatro House, Frimley Road, Camberley, United Kingdom GU16 7ER

Tel: 08452 30 40 30 Fax: 08452 30 50 30 E-mail: info@labunlimited.com E-mail: info@labunlimited.co.uk Web: www.labunlimited.com Web: www.labunlimited.co.uk