

Wide neck bottles series 310 HDPE - square, natural colour and white

The wide neck opening means that the bottles are easy to fill and clean.

Square cross-section.

Recessed labelling field.

Three different types of closures available.

Ideal for liquids, pastes and granulates.

Application:

- packaging - sending samples - transport - storage

Volume Discounts Apply



UN approved

Available ex stock with UN approval for solids. The approvals are only valid with the closures offered here.

Wide neck bottles series 310 HDPE, natural colour and white

Nominal capacity ml	A mm	B mm	C mm	D mm	Mould no.	No. per case* pcs.	UN-approval	Order no. natural col.	Order no. white
50	36 x 36	21,4	28	88,5	310-115	700	-	9.072 869	9.072 879
100	45 x 45	24	32	101	310-88	582	UN 3H1/X0.3/S	9.072 870	9.072 880
250	59 x 59	35	45	129,5	310-86	273	UN 3H1/X0.8/S	9.072 871	9.072 881
350	60 x 60	35	45	164	310-46	190	-	9.072 877	-
500	75 x 75	42,5	54	162	310-89	124	UN 3H1/X1.5/S	9.072 872	9.072 882
750	83 x 83	42,5	54	188	310-48	90	UN 3H1/X2.5/S	9.072 878	9.072 883
1000	85 x 85	42,5	54	214	310-85	81	UN 3H1/X2.5/S	9.072 873	9.072 884
1500	108 x 108	68,5	80	209	310-93	42	UN 3H1/X4/S	9.072 874	9.072 885
2500	122 x 122	68,5	80	260	310-51	27	UN 3H1/X5/S	9.072 875	9.072 886
4000	142 x 142	68,5	80	294	310-52	15	UN 3H1/X5.5/S	9.072 876	9.072 887

* Packaging unit / minimum ordering quantity

A = Dimension cross section body B = Neck i.d. C =Cap size o.d. D = Height without closure

①	Screw closure PP, black, with PE-foam insert
②	Screw closure PP, blue, with tamper evident ring and conical seal
③	Screw closure PP, red, with tamper evident ring and PE-foam insert
④	Wrench: A tool suitable for all sizes of closures Order no. 805-85102

For bottle nom. cap. ml	C mm	① Order no.	② Order no.	③ Order no.
50	28	9.072 894	9.072 904	9.072 899
100	32	9.072 895	9.072 092	9.072 900
250/350	45	9.072 896	9.072 093	9.072 901
500/750/1000	54	9.072 897	9.072 094	9.072 902
1500/2500/4000	80	9.072 898	9.072 893	9.072 903



Volume Discounts Apply



Narrow neck bottles series 310 HDPE - square, natural colour

Square cross-section.
Recessed labelling field.
Four different types of closures available.
Ideal for liquids, leakproof.

Application:

- packaging - sending samples - transport - storage

Volume Discounts Apply

Narrow neck bottles series 310 HDPE - square, natural colour

Nominal capacity ml	A mm	B mm	C mm	D mm	Mould no.	No. per case* pcs.	Order no. natural col.
100	45 x 45	17	32	108	310-53 A	544	9.072 850
250	60 x 60	17	32	136	310-54 A	247	9.072 851
500	74 x 74	17	32	170	310-55 A	122	9.072 852
1000	85 x 85	17	32	221	310-56 A	72	9.072 853
2500	122 x 122	28	45	272	310-57 E	27	9.072 854

* Packaging unit / minimum ordering quantity

A = Dimension cross section body B = Neck i.d. C = Cap size o.d. D = Height without closure

UN Narrow neck bottles series 308/310 HDPE - natural colour and black

Narrow neck bottles with UN approval.
Particularly suitable for filling and transport of hazardous liquids.
Recessed labelling field with sizes 1000 and 2500 ml.
Three different types of closures available.

Application:

- packaging - sending samples - transport - storage



UN Narrow neck bottles series 308/310 HDPE - natural colour and black

Nominal capacity ml	A mm	B mm	C mm	D mm	Weight g	Mould no.	No. per case* pcs.	UN-approval	Order no. natural col.	Order no. black
500	77 x 77	28,7	45	176,5	52	308-257	120	UN 1H1/X1.9/250	9.073 373	9.073 376
1000	96 x 96	28,7	45	211,8	92	310-135	64	UN 3H1/X1.9/250	9.073 374	9.073 377
2500	127 x 127	28,7	45	282	167	310-136	26	UN 3H1/X1.9/250	9.073 375	9.073 378

* Packaging unit / minimum ordering quantity



Screw closures PP

- ① red, with pull-off ring and PE-foam insert
- ② blue, with pull-off ring and conical seal
- ③ white/blue, child-resistant, with tamper evident ring and conical seal
- ④ white/red, child-resistant, with tamper evident ring, conical seal and degassing membrane
- ⑤ **NEW:** Dedassing closure, red, with tamper evident ring, conical seal and degassing membrane (without picture)

Screw closures PP

Nominal capacity ml	C mm	① Order no.	② Order no.	③ Order no.	④ Order no.	⑤ Order no.
10 - 1000	32	9.072 863	9.072 041	9.072 866	9.072 867	9.072 905
2500	45	9.072 868	9.072 862	-	-	-

Volume Discounts Apply