

# LABORTOPS

valid until 30<sup>th</sup> of June 2017



**EXTRA**

**Safety Storage Cabinets**

## Safety Storage Cabinets S-CLASSIC-90 with Wing Doors



- Safe and approved storage of hazardous materials in workrooms
- Type proved according to EN 14470-1, fire resistance 90 minutes
- GS-approved, CE conformal
- Outer body made of sheet steel, colour light grey (RAL 7035) or warning yellow (RAL 1004), interior with chemical resistant surface
- Lockable doors with cylinder lock (integration into a locking system possible)
- Adjusting aids to compensate for uneven floor
- Ventilation connection (NW 75) in the cabinet ceiling
- Standard door open arrest system holds doors open, allowing for safe and convenient storage and removal

3 shelves, height adjustable (32 mm increments), load capacity 75 kg, bottom collecting sump with perforated insert as a first storage level.

### Specifications

Outer dimensions (W x D x H)

S90.196.120.WDAS:

1196 x 616 x 1968 mm

S90.196.090.WDAS:

896 x 616 x 1968 mm

S90.196.060.WDAS:

596 x 616 x 1968 mm

Inner dimensions (W x D x H)

S90.196.120.WDAS:

1050 x 524 x 1742 mm

S90.196.090.WDAS:

750 x 524 x 1742 mm

S90.196.060.WDAS:

450 x 524 x 1742 mm

Shelf dimensions (W x D x H)

S90.196.120.WDAS:

1040 x 432 x 30 mm

S90.196.090.WDAS:

740 x 432 x 30 mm

S90.196.060.WDAS:

440 x 432 x 30 mm

**A wide variety of cabinet versions is available. Please enquire.**

**-14%**

Description	PK	Cat. No.	Price	Now
			GBP	GBP
Safety Storage Cabinets S90.196.120.WDAS, light grey	1	6.260 928	<del>2,598.30</del>	2,234.54
Safety Storage Cabinets S90.196.120.WDAS, warning yellow	1	6.270 910	<del>2,598.30</del>	2,234.54
Safety Storage Cabinets S90.196.090.WDAS, light grey	1	6.286 539	<del>2,494.93</del>	2,145.64
Safety Storage Cabinets S90.196.090.WDAS, warning yellow	1	6.286 540	<del>2,494.93</del>	2,145.64
Safety Storage Cabinets S90.196.060.WDAS, light grey	1	6.269 868	<del>2,016.43</del>	1,734.13
Safety Storage Cabinets S90.196.060.WDAS, warning yellow	1	6.286 542	<del>2,016.43</del>	1,734.13
Shelf for S90.196.120.WDAS	1	9.840 532	<del>22.58</del>	28.02
Shelf for S90.196.090.WDAS	1	6.286 541	<del>22.58</del>	28.02
Shelf for S90.196.060.WDAS	1	9.840 514	<del>26.49</del>	22.78
Recirculating air filter system (option)	1	9.840 650	<del>1,431.42</del>	1,231.02
Ventilation system (option)	1	9.840 616	<del>549.72</del>	472.76

Safety Storage Cabinets K-PHOENIX-90



Combined and unrestricted storage of flammable and aggressive hazardous materials in working areas, in accordance with EN 14470-1

- Type-tested in accordance with EN 14470-1
- Fire resistance of 90 minutes (type 90)
- GS-tested
- Two in one: 2 compartments (vertically separated) allow the storage of flammable liquids and aggressive chemicals in one cabinet
- Robust construction and longevity: body made of powder - coated sheet steel, safety elements assembled outside the cabinet for increased protection against corrosion
- User-friendly and safe (folding door compartment): door open arrest system (standard), allow for convenient and safe removal or storage of containers
- No unauthorised use: doors lockable with profile cylinder (integration in an existing locking system possible)
- Easy alignment: adjusting aids to compensate for uneven floor
- Ventilation: integrated air ducts ready for connection (DN 75) to a technical exhaust system

Model K90.196.120.MF.FWAS in Light grey (RAL 7035), interior equipment with 6x drawer, 1x bottom collecting sump (sheet steel powder-coated), 6x pull-out shelf (melamine/PP) extraction unit optional

Extraction unit, safe ventilation of safety storage cabinets, plug-in - fast and easy to install



-14%

Description	Dimensions (W x D x H) mm	PK	Cat. No.	Price GBP	Now GBP
Safety storage cabinet, left: 3x shelf, 1 perforated insert, 1 bottom collecting sump, right: 4x pull-out shelf	1197 x 617 x 1969	1	6.259 719	<del>2,967.56</del>	2,552.10
Safety storage cabinet, left: 6x drawer, 1x bottom collecting sump, right: 6x pull-out shelf	1197 x 617 x 1969	1	6.284 342	<del>3,692.31</del>	3,175.39
Sump inliner for drawer, PP grey	-	1	9.840 626	<del>59.99</del>	51.59
Anti-slip rubber mat for drawer, black	-	1	6.237 260	<del>11.88</del>	10.22
Shelf for 6.259 719	-	1	9.840 514	<del>26.49</del>	22.78
Ventilation system	-	1	9.840 616	<del>549.72</del>	472.76

Disposal Bin, PE



- Not approved for waste materials containing solvents
- Lid opened by foot pedal
- Lid self-closing

-14%

Capacity litres	Dimensions (dia. x H) mm	PK	Cat. No.	Price GBP	Now GBP
23	420 x 400	1	9.840 640	<del>89.38</del>	76.87
38	420 x 470	1	9.840 641	<del>96.68</del>	83.14
53	560 x 530	1	9.840 642	<del>99.58</del>	85.64

Disposal Bin, stainless steel



- Approved for waste materials containing solvents
- Lid opened by foot pedal
- Lid is self-closing
- Base edge has ventilation openings; prevents heating up of the base and spontaneous combustion of contents

-9%

Capacity litres	Dimensions (dia. x H) mm	PK	Cat. No.	Price GBP	Now GBP
20	270 x 500	1	9.840 643	<del>369.92</del>	336.63
50	390 x 600	1	9.840 644	<del>486.11</del>	442.36

## Fire Resistant Gas Cylinder Cabinets G90 Series



- Colour: grey RAL 7035
- Exhaust air connection (NW75) on top of the cabinet
- Adjusting aids to compensate for uneven floor
- Large interior height
- Lead-through points for pipes, sensors etc. in the ceiling of the cabinet
- Complete with mounting rails, rolling ramp, bottle retainer with ratchet straps

Approved and safe

- Labelling in accordance with EN 14470-2, fire resistance 90 minutes
- GS approved, CE compliant



Lateral cylinder retainer for 50 l gas cylinders, sheet steel powder-coated, height adjustable without tools



Radial ventilator 230/400 V, IP 55, EX-Zone 2  
Technical Specifications on request.



**-14%**

Description	Dimensions (W x D x H) mm	Weight kg	PK	Cat. No.	Price GBP	Now GBP
For one 50 litre bottle	600 x 615 x 2050	365	1	6.265 442	<del>2,446.84</del>	<b>2,104.28</b>
For two 50 litre bottles	600 x 615 x 2050	365	1	6.261 842	<del>2,833.42</del>	<b>2,436.74</b>
For three 50 litre bottles	900 x 615 x 2050	490	1	9.840 609	<del>2,787.84</del>	<b>2,397.54</b>
For four 50 litre bottles	1200 x 615 x 2050	610	1	9.840 608	<del>3,035.34</del>	<b>2,610.39</b>
For four 50 litre bottles	1400 x 615 x 2050	690	1	6.267 279	<del>3,447.83</del>	<b>2,965.13</b>
Lateral cylinder retainer			1	6.286 543	<del>115.67</del>	<b>99.48</b>
Radial ventilator			1	6.253 657	<del>340.54</del>	<b>292.86</b>

Price plus disposal in the lab and installation.

\* Upon conclusion of an asecos service and maintenance agreement (BasicPlus tariff) with a fixed term of 5 years.



# LABORTOPS



## WORLD NOVELTY:

### V-LINE

### THE NEW FORMULA FOR THE STORAGE OF HAZARDOUS MATERIALS

Vertical drawer technology for  
simplified operation.

COMING SOON!



[www.asecos.global](http://www.asecos.global)