



WTW Safety Data Sheet
According to EC Directives 91/155/EEC

1. Identification of the substance/preparation and of the company/undertaking

Product name:

pH 4.006 Standard Buffer

Ord. No.: 109 110

Model: PL 4

Manufacturer:

WTW Wissenschaftlich-Technische Werkstätten GmbH, Dr.-Karl-Slevogt-Straße 1, D-82362 Weilheim
Germany, Tel.: (0881) 183-0, Fax.: (0881) 62539

Emergency telephone:
(001)-703-527-3887

2. Composition/information on ingredients

Chemical characterization: Aqueous saline solution

Ingredients: Potassium hydrogen phthalate

Hazard symbols: No hazardous product

R-phrases: /

Content: 1 %

CAS-No.: 877-24-7

EINECS-No.: 212-889-4

EG-No.: 212-889-4

3. Hazards identification

No hazardous product.

4. First aid measures

After inhalation: Fresh air

After skin contact: Wash off with plenty of water.
Remove contaminated clothing.

After eye contact: Rinse out with plenty of water
with the eyelid held wide open.

After ingestion: Make injured person drink plenty of water,
induce vomiting, summon doctor.

WTW Safety Data Sheet

Identification of the substance/preparation

pH 4.006 Standard Buffer

5. Fire-fighting measures

Suitable extinguishing media: To suit environment.
Unsuitable extinguishing media: No information available.
Special risks: None
Other information: Non combustible

6. Accidental release measures

Personal-related precautionary measures: None
Environmental-protection measures: None
Procedures for cleaning/absorption: Take up with liquid-absorbent material.
Forward for disposal.
Clean up affected area with water.

7. Handling and storage

Handling: No further requirements.
Storage: Store tightly closed.

8. Exposure control/personal protection

Specific control parameters: None
Personal protective equipment
Respiratory protection: Required when dusts are generated.
Eye protection: Protective glasses
Hand protection: Gloves
Industrial hygiene: Wash hands and face after working with solution.
Change contaminated clothing.

WTW Safety Data Sheet

Identification of the substance/preparation

pH 4.006 Standard Buffer

9. Physical and chemical properties

<i>Form:</i> Liquid	<i>Colour:</i> Colourless	<i>Odour:</i> Odourless
<i>pH-value</i>	4,006 at 25°C	
<i>Melting temperature:</i>	°C	
<i>Boiling temperature:</i>	about 100°C	
<i>Ignition temperature:</i>	°C	
<i>Flash point:</i>	°C	
<i>Explosion limite lower:</i>	/	
<i> upper:</i>	/	
<i>Density:</i>	/	
<i>Solubility in water:</i>	miscible	
<i>Thermal decomposition:</i>	/	

10. Stability and reactivity

<i>Conditions to be avoided:</i>	No information available.
<i>Substances to be avoided:</i>	No information available.
<i>Hazardous decomposition products:</i>	No information available.
<i>Further information:</i>	No information available.

11. Toxicological information

<i>Acute toxicity:</i>	Quantitative data on the toxicity are not available.
<i>Further toxicological information:</i>	The product should be handled with the care usual when dealing with chemicals.

12. Ecological information

<i>Ecotoxic effects:</i>	No ecological problems are to be expected when the product is handled and used with due care and attention.
--------------------------	---

WTW Safety Data Sheet

Identification of the substance/preparation

pH 4.006 Standard Buffer

13. Disposal considerations

Products: There are no uniform regulations for the disposal of chemicals or residues. The disposal is regulated through corresponding regional laws and regulations. We recommend that you contact either the authorities in charge or approved waste disposal companies which will advise you on how to dispose the special waste.

Packaging: Handle contaminated packaging in the same way as the substance itself. If not officially specified differently, noncontaminated packaging may be treated like household waste or recycled.

14. Transport information

Not subject to the transport regulations.

15. Regulatory information

Labelling according to EC directives

Symbol: No hazardous product
R-phrases: /
S-phrases: /

German regulations

Water pollution class; nonpolluting solution

Regulation on combustible liquids; VbF: Non combustible

Other regulations: None

16. Other information

Date of issue: 22.1.2004
Supersedes edition of 4.12.2002
Reason for alteration: General update

The information contained herein is based on the present state of our knowledge. It characterizes the product with regard to the appropriate safety precautions. It does not represent a guarantee of the properties of the product.