

WTW Safety Data Sheet
According to EC Directives 91/155/EEC

1. Identification of the substance/preparation and of the company/undertaking

Product name:

Kal.Kit Turb 430 T

Ord. No.: 600 561
Model: Kal.Kit Turb 430 T

Manufacturer:

WTW Wissenschaftlich-Technische Werkstätten GmbH, Dr.-Karl-Slevogt-Straße 1, D-82362 Weilheim, Germany, Tel.: (0881) 183-0, Fax.: (0881) 62539

Emergency telephone:
(001)-703-527-3887

2. Composition/information on ingredients

Chemical characterization: Aqueous solution

Hazardous ingredients: Sodium azide and Divinyl benzene

Sodium azide	Divinyl benzene
<i>Hazard symbols:</i> very toxic (T+)	<i>Hazard symbol:</i> irritant, dangerous for the environment (Xi, N)
<i>R-phrases:</i> 28-32-50/53	<i>R-phrases:</i> 36/37/38-51/53
<i>Content:</i> < 0,01 %	<i>Content:</i> < 1 %
<i>CAS-No.:</i> 26628-22-8	<i>CAS-No.:</i> 1321-74-0
<i>EC-No.:</i> 247-852-1	<i>EC-No.:</i> 215-325-5

3. Hazards identification

No hazardous product

4. First aid measures

After inhalation: Fresh air

After skin contact: Wash off with plenty of water Remove contaminated clothing. Get medical attention immediately

After eye contact: Rinse out with plenty of water. Held the eyelid wide open. Summon eye specialist.

After ingestion: Consult doctor if feeling unwell.

WTW Safety Data Sheet

Identification of the substance/preparation

Kal.Kit Turb 430 T

5. Fire-fighting measures

Suitable extinguishing media: Water, powder, foam and carbondioxide

Unsuitable extinguishing media: No information available

Special risks: None

Other information: Non combustible

6. Accidental release measures

Personal-related precautionary measures: Avoid generation of aerosols.
Avoid substance contact.

Environmental-protection measures: Do not allow to enter sewerage system.

Procedures for cleaning/absorption: Take up with liquid-absorbent material.
Forward for disposal.
Clean up effected area with water.

7. Handling and storage

Handling: Do not open containers for no apparent reason.
If containers are opened, the value of the standard is not longer valid.
Do not shake. Shaking induces air bubbles and increases turbidity.

Storage: Store away from extremes in temperature. Prevent from freezing.
Freezing causes flocculation.

8. Exposure control/personal Protection

Specific control parameters: MAK Germany (max. workplace conc.): 0,2 mg/m³ for sodium azide.

Personal protective equipment

Respiratory protection: Required when aerosols are generated.

Eye protection: Protective glasses

Hand protection: Gloves

Industrial hygiene: Wash hands and face after working with solution.
Change contaminated clothing.
Do not eat, drink and smoke at workplace.

WTW Safety Data Sheet

Identification of the substance/preparation

Kal.Kit Turb 430 T

9. Physical and chemical properties

<i>Form:</i> liquid	<i>Colour:</i> colourless	<i>Odour:</i> odourless
<i>pH-value</i>	6.7 / 25°C	
<i>Melting temperature:</i>	0°C	
<i>Boiling temperature:</i>	100°C	
<i>Ignition temperature:</i>	/ °C	
<i>Flash point:</i>	/ °C	
<i>Explosion limite lower:</i>	/	
<i> upper:</i>	/	
<i>Density:</i>	1,00 g/cm ³	
<i>Solubility in water:</i>	miscible	
<i>Thermal decomposition:</i>	/	

10. Stability and reactivity

<i>Conditions to be avoided:</i>	No extremes in temperature
<i>Substances to be avoided:</i>	No acids, No ions of heavy metals (Generation of precipitations is possible. These are explosive in dry state.)
<i>Hazardous decomposition products:</i>	None
<i>Further information:</i>	None

11. Toxicological information

<i>Acute toxicity:</i>	LD 50 (oral, rat): 27 mg/ kg for sodium azide.
<i>Further toxicological information:</i>	No teratogenic effect in animal experiments. After inhalation of aerosols and substance contact: Irritation effects in the respiratory tract and of eyes and skin. After ingestion: Nausea, vomiting, headache, restlessness, and diarrhea.

12. Ecological information

<i>Ecotoxic effects:</i>	Do not allow to enter sewerage system. No ecological problems are to be expected when the product is handled and used with due care and attention.
--------------------------	---

WTW Safety Data Sheet

Identification of the substance/preparation

Kal.Kit Turb 430 T

13. Disposal considerations

Products: There are no uniform regulations for the disposal of chemicals or residues. The disposal is regulated through corresponding regional laws and regulations. We recommend that you contact either the authorities in charge or approved waste disposal companies which will advise you on how to dispose the special waste.

Packaging: Handle contaminated packaging in the same way as the substance itself. If not officially specified differently, noncontaminated packaging may be treated like household waste or recycled.

14. Transport information

Not subject to transport regulations.

15. Regulatory information

Labelling according to EC directives

Symbol: No hazardous product
R-phrases: /
S-phrases: /

German regulations

Water pollution class; nonpolluting substance

Regulation on combustible liquids; VbF: Non combustible

Other regulations: None

16. Other information

Date of issue: 5.12.2005
Supersedes edition of: First edition
Reason for alteration: First edition

The information contained herein is based on the present state of our knowledge. It characterizes the product with regard to the appropriate safety precautions. It does not represent a guarantee of the properties of the product.