

## DrySyn Classic

### Single Position Heating Blocks

*25 mL to 1000 mL*

The original member of the DrySyn family, the DrySyn Classic allows the user to heat any size standard round-bottom flask from 25 mL to 1000 mL.

The base unit of the DrySyn Classic accepts a standard 1000 mL round bottom flask and allows the inserts to be used - one at a time - to adapt it to all standard flask sizes from 500 mL to 25 mL.

All single-position base units include two heat-resistant handles for safe handling when hot and adjustable feet to allow complete compatibility with any hotplate stirrer.



**Ask us about other single position heating options!**



## DrySyn Classic

### Single Position Heating Blocks

25 mL to 1000 mL

- ✓ Removes risk of oil fires
- ✓ No risks associated with oil spillages
- ✓ No toxic fume emission
- ✓ Reduction in electrical energy consumption
- ✓ No costly maintenance or disposal of oil



#### Popular Purchase Options\*

##### ADS1: DrySyn Classic Kit

Includes 1000 mL base and 4 x adapter inserts for 50, 100, 250 & 500 mL round bottom flasks.


##### ADS1-25: DrySyn Classic 25 mL Adapter Insert

##### ADS1-A: DrySyn Classic Starter Kit

Includes 1000 mL base and 1 x adapter insert (specify size).

### DrySyn Classic Heat Up Times

100 mL Flasks & Inserts, 50 mL Water



Ask us about other single position heating options!