

Chest freezer LT C150

- Low cost of ownership
- Limited maintenance
- On-site service

- Easily installable & movable
- Efficient cyclopentane insulation
- ISO 9001, 14001 & 18001 OHSAS Certified

Specifications

Technical

| | |
|--------------------------|----------------|
| Gross volume | 140 litres |
| Net Volume | 133 liters |
| Temp. range | -10°C to -45°C |
| Temp. probes | 1/PTC |
| Compressor type | Secop |
| Baskets (incl./optional) | 1/2 |
| Refrigerant | R 290 |
| Noise level | 55 dB |
| Power consumption 24/h | 2,44 kWh |
| Power consumption | 891 kWh / year |
| Voltage | 220V-240V |
| Frequency | 50Hz |

Features

| | |
|--|--------|
| Audible & visual alarm | YES |
| Castors | 4 |
| Integrated lock | YES |
| Inner lids | NO |
| Defrost | Manual |
| Battery Backup | NO |
| USB | NO |
| Data connection | NO |
| Digital controller with display | XR30CX |
| Data logging | NO |
| High/Low temperature alarm | YES |
| Porthole for external probe & CO2 backup | Ø12 mm |
| Prepared for external logging equipment | YES |

Material

| | |
|-------|-----------------------------|
| Inner | Precoated white steel |
| Outer | Galvanized pre-coated steel |

Outer dimensions

| | |
|-------------------|--------|
| Height | 850 mm |
| Height w. castors | 890 mm |
| Width | 720 mm |
| Depth | 695 mm |
| Net weight | 55 kg |

Inner dimensions

| | |
|-------------------|------------|
| Height | 635/385 mm |
| Width | 350/210 mm |
| Depth | 440 mm |
| Insulation | 80 mm |
| Insulation type | PUF |
| Internal humidity | <5 % |

LT C150



Product no.

N111002

Storage capacity

| Box height | No. of boxes | No. of racks | Total no. of boxes | Product no. |
|------------|--------------|--------------|--------------------|-------------|
| 50 mm | 6 | 3 | 18 | N311200 |
| | 11 | 6 | 66 | N311220 |
| | | | = 84 | |
| 75 mm | 4 | 3 | 12 | N311300 |
| | 8 | 6 | 48 | N311320 |
| | | | = 60 | |
| 100 mm | 3 | 3 | 9 | N311400 |
| | 5 | 6 | 30 | N311420 |
| | | | = 39 | |