Comparison of Gel/PCR Extraction Kits- IBI and Brand Q

	Purity & Stability			Each Sample Tested 3 X		
<u>IBI Brand</u>	DH5a/TA			BL21/ pET20b		
Volume (ml) DNA Conc. (ng/ λ) Total DNA (μg) A260/ A280	35.50 69.80 2.48 1.85	36.50 66.60 2.43 1.97	35.00 71.30 2.50 1.86	43.00 25.20 1.08 1.80	44.50 24.50 1.10 1.70	41.00 36.20 1.48 1.88
Brand Q	DH5a/TA			BL21/ pET20b		
Volume (ml) DNA Conc. (ng/ λ) Total DNA (μg) A260/ A280	37.00 48.80 1.81 1.63	37.00 62.20 2.30 1.63	37.00 50.40 1.86 1.66	46.00 13.90 0.64 1.74	46.00 17.60 0.81 1.79	46.00 18.20 0.84 1.83





This graph is a helpful tool to determine what kit to acquire. Some kits can be used for both PCR Clean-up & Gel Extraction.





