DURAN® HPLC RESERVOIR BOTTLE GL 45

with conical base

- DURAN[®] reservoirs feature one central cavity for the inlet filter to allow the **most efficient delivery** of mobile phase to high-performance liquid chromatography (HPLC) analyzers.
- The conical bottom reservoirs allow a greater recovery of solvent, without having to tilt the reservoir and risk a solvent spill.
- On the two larger sizes, the central cavity can accommodate magnetic stir bars for improved mixing.
- Side vent / drain hole in the base for ultrasonic bath degassing, or water drainage after cleaning.
- Retrace Code and certificate for manufacturing lot traceability.
- Made from borosilicate glass conforming to USP Type I and ASTM E438, Type I, Class A requirements.
- **Note:** Bottle is provided without cap! Please order cap separately.



DURAN® HPLC reservoir bottle GL 45 with conical base					
Description	Volume [mL]	d [mm]	h [mm]	Pack Unit	Cat. No.
DURAN® HPLC reservoir bottle, clear, conical base, GL 45	1000	110	295	1	218415407
DURAN® HPLC reservoir bottle, clear, conical base, GL 45	2000	145	309	1	218416309
DURAN® HPLC reservoir bottle, clear, conical base, GL 45	5000	190	386	1	218417305
DURAN® HPLC reservoir bottle, clear, conical base, GL 45	10000	235	481	1	218418601

Accessories for the DURAN® HPLC reservoir bottle GL 45					
Description	Pack Unit	Cat. No.			
DURAN® 3-port GL 45 HPLC delivery cap 1/4"-28 thread for 3.2 mm OD tubing PTFE complete with three 1/4"-28 ETFE plugs and TFE/propylene O-ring	1	292312801			
DURAN® ETFE blanking plug 1/4"-28 thread for 3-port HPLC GL 45 delivery cap	1	292310003			