






**E 103/1****Insert**

For the optimum loading of test tubes up to 12 x 75 mm.

- Suitable for about 200 Test tubes (up to 12 x 75 mm)
- Divided into 6 surfaces
- The mesh width is 8 x 8 mm
- Including lid A 13

**Product details - E 103/1**

EAN: 4002514231051 / Article no. 69510302D / Mat. no. 06907630

Product category 	
Washer-disinfectors, laboratory	•
Machine/accessory categorisation 	
PG 8535	•
PG 8536	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category 	
Insert/module for test tubes	•
Product properties 	
Material	Stainless steel
Colour	Stainless steel
Area of application 	
Reprocessing laboratory glassware	•
Capacity 	
Test tubes [number]	200
Dimensions and weight 	
External dimensions, height in mm	122
External dimensions, width in mm	200
External dimensions, depth in mm	320

Net weight in kg

0,540

Subject to technical changes; no liability accepted for the accuracy of the information given!

Distributed by:



Tallaght Business Park
Whitestown, Dublin 24,
Ireland
D24 RFK3

Tel: (01) 4523432
Fax: (01) 4523967
E-mail: info@labunlimited.com
Web: www.labunlimited.com

Quatro House, Frimley Road,
Camberley,
United Kingdom
GU16 7ER

Tel: 08452 30 40 30
Fax: 08452 30 50 30
E-mail: info@labunlimited.co.uk
Web: www.labunlimited.co.uk