

Tallaght Business Park Whitestown, Dublin 24, Ireland **D24 RFK3**

Tel: (01) 4523432 Fax: (01) 4523967 Web: www.labunlimited.com

Quatro House, Frimley Road, Camberley, United Kingdom **GU16 7ER**

Tel: 08452 30 40 30 Fax: 08452 30 50 30 E-mail: info@labunlimited.com E-mail: info@labunlimited.co.uk Web: www.labunlimited.co.uk

Home Products Laboratory Technology Product selection accessories laboratory technology E 106/2



E 106/2

Insert

For the optimum loading of wide neck flasks, measuring cylinders, etc.

• equipped with 13 large spring hooks x 175 mm height

• The distance between the spring hooks is approx. 85 mm



Product details - E 106/2

EAN: 4002511115699 / Article no. 69510603D / Mat. no. 03808330

Product category		^
Washer-disinfectors, laboratory	•	
Large-capacity washer-disinfectors, laboratory		
Machine/accessory categorisation		^
G 7825	٥	
PG 8527	٠	
PG 8535	•	

PG 8536	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category	^
Insert/module for beakers	•
Insert/module for wide-necked glassware	•
Product properties	^
Material	Stainless steel
Colour	Stainless steel

Area of application		^
Reprocessing laboratory glassware	٠	
Dimensions and weight		^
External dimensions, height in mm	420	
External dimensions, width in mm	180	
External dimensions, depth in mm	186	
Net weight in kg	0.57	

Subject to technical changes; no liability accepted for the accuracy of the information given!

Distributed by:

Lab Unlimited

Tallaght Business Park Whitestown, Dublin 24, Ireland D24 RFK3 Quatro House, Frimley Road, Camberley, United Kingdom GU16 7ER

 Tel: (01) 4523432
 Tel: 08452 30 40 30

 Fax: (01) 4523967
 Fax: 08452 30 50 30

 E-mail: info@labunlimited.com
 E-mail: info@labunlimited.co.uk

 Web: www.labunlimited.com
 Web: www.labunlimited.co.uk