

E 118







Insert
for the optimum loading of half petri dishes.

- Suitable for 38 half petri dishes with Ø 100 mm
- equipped with 38 holders, height 70 mm, distance approx. 26 mm
- For use in upper basket



Product details - E 118

EAN: 4002511115767 / Article no. 69511801D / Mat. no. 03830270

Product category 	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine/accessory categorisation 	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category 	
Insert/module for petri dishes	•
Product properties 	
Material	Stainless steel
Colour	Stainless steel
Capacity 	
Petri dishes [number]	38
Dimensions and weight 	
External dimensions, height in mm	120
External dimensions, width in mm	460
External dimensions, depth in mm	445

Net weight in kg

1,600

Subject to technical changes; no liability accepted for the accuracy of the information given!

Distributed by:

Lab Unlimited
CARL STUART GROUP

Tallaght Business Park
Whitestown, Dublin 24,
Ireland
D24 RFK3

Tel: (01) 4523432
Fax: (01) 4523967

E-mail: info@labunlimited.com
Web: www.labunlimited.com

Quatro House, Frimley Road,
Camberley,
United Kingdom
GU16 7ER

Tel: 08452 30 40 30
Fax: 08452 30 50 30

E-mail: info@labunlimited.co.uk
Web: www.labunlimited.co.uk