



E 134








Insert
for optimum loading of 210 slides.

- suitable for 210 slides
- equipped with 26 compartments at 26 x 11 mm, wire thickness 3 mm



Product details - E 134

EAN: 4002511115897 / Article no. 69513401D / Mat. no. 03808600

Product category 	
Washer-disinfectors, laboratory	•
Machine/accessory categorisation 	
PG 8535	•
PG 8536	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category 	
Other inserts	•
Product properties 	
Material	Stainless steel
Colour	Stainless steel
Number of compartments	26
Area of application 	
Reprocessing laboratory glassware	•
Capacity 	
Slides [number]	210
Dimensions and weight 	
External dimensions, height in mm	73
External dimensions, width in mm	200

External dimensions, depth in mm

445

Net weight in kg

1,410

Subject to technical changes; no liability accepted for the accuracy of the information given!

Distributed by:

Lab Unlimited
CARL STUART GROUP

Tallaght Business Park
Whitestown, Dublin 24,
Ireland
D24 RFK3

Tel: (01) 4523432

Fax: (01) 4523967

E-mail: info@labunlimited.com

Web: www.labunlimited.com

Quatro House, Frimley Road,
Camberley,
United Kingdom
GU16 7ER

Tel: 08452 30 40 30

Fax: 08452 30 50 30

E-mail: info@labunlimited.co.uk

Web: www.labunlimited.co.uk