



## E 355

Mobile unit  
for the optimum loading of narrow necked  
glassware.

- equipped with 16 injector nozzles
- Half can be equipped with other inserts



## Product details - E 355

EAN: 4002511116191 / Article no. 69535501D / Mat. no. 03809560

Product category		^
Washer-disinfectors, laboratory		•
Machine/accessory categorisation		^
PG 8535		•
PG 8536		•
Product category		^
Mobile unit for narrow-necked glassware		•
Product properties		^
Material	Stainless steel	
Colour	Stainless steel	
Area of application		^
Reprocessing laboratory glassware		•
Capacity		^
Narrow-necked glassware [number]	16	
Dimensions and weight		^
External dimensions, height in mm	502	
External dimensions, width in mm	535	
External dimensions, depth in mm	517	
Net weight in kg	5,800	
Standard accessories		^
F 351 injector nozzles Ø 4 mm. length 160 mm [number]	7	

E 352 injector nozzles Ø 6 mm, length 220 mm [number]	9
E 353 clips for injector nozzles E 351 [number]	7
E 354 clips for injector nozzles E 352 [number]	9
Rinse nozzle for powder dispenser [number]	1

Subject to technical changes; no liability accepted for the accuracy of the information given!

Distributed by:

**Lab Unlimited**  
CARL STUART GROUP

Tallaght Business Park  
Whitestown, Dublin 24,  
Ireland  
D24 RFK3

Tel: (01) 4523432  
Fax: (01) 4523967  
E-mail: [info@labunlimited.com](mailto:info@labunlimited.com)  
Web: [www.labunlimited.com](http://www.labunlimited.com)

Quatro House, Frimley Road,  
Camberley,  
United Kingdom  
GU16 7ER

Tel: 08452 30 40 30  
Fax: 08452 30 50 30  
E-mail: [info@labunlimited.co.uk](mailto:info@labunlimited.co.uk)  
Web: [www.labunlimited.co.uk](http://www.labunlimited.co.uk)