



E 414








Mobile unit
for the optimum loading of narrow necked
glassware.

- equipped with 37 injector nozzles
- Equipped with connection for hot air drying
- Bracket for ML/2 for mobile unit recognition / spray arm control



Product details - E 414

EAN: 4002511116733 / Article no. 69541401D / Mat. no. 03810190

Product category 	
Washer-disinfectors, laboratory	•
Machine/accessory categorisation 	
PG 8535	•
PG 8536	•
Product category 	
Mobile unit with connection for hot air drying	•
Mobile unit for narrow-necked glassware	•
Product properties 	
Material	Stainless steel Plastic
Colour	Stainless steel Black
Area of application 	
Reprocessing laboratory glassware	•
Capacity 	
Narrow-necked glassware [number]	37
Dimensions and weight 	
External dimensions, height in mm	502
External dimensions, width in mm	535
External dimensions, depth in mm	517
Net weight in kg	6,100

Standard accessories	
ID 90, injector nozzles Ø 2.5 mm, length 90 mm [number]	4
ID 110, injector nozzles Ø 2.5 mm, length 110 mm [number]	3
ID 140, injector nozzles Ø 4 mm, length 140 mm [number]	5
ID 160, injector nozzles Ø 4 mm, length 160 mm [number]	5
ID 180, injector nozzles Ø 4 mm, length 180 mm [number]	5
ID 200, injector nozzles Ø 6 mm, length 200 mm [number]	5
ID 220, injector nozzles Ø 6 mm, length 220 mm [number]	5
ID 240, injector nozzles Ø 6 mm, length 240 mm [number]	5
Rinse nozzle for powder dispenser [number]	1

Subject to technical changes; no liability accepted for the accuracy of the information given!

Distributed by:

Lab Unlimited
CARL STUART GROUP

Tallaght Business Park
Whitestown, Dublin 24,
Ireland
D24 RFK3

Tel: (01) 4523432

Fax: (01) 4523967

E-mail: info@labunlimited.com

Web: www.labunlimited.com

Quatro House, Frimley Road,
Camberley,
United Kingdom
GU16 7ER

Tel: 08452 30 40 30

Fax: 08452 30 50 30

E-mail: info@labunlimited.co.uk

Web: www.labunlimited.co.uk