



ark Quatro House, Frimley Road 24, Camberley, United Kingdom GU16 7ER

 Tel: (01) 4523432
 Tel: 08452 30 40 30

 Fax: (01) 4523967
 Fax: 08452 30 50 30

 E-mail: info@labunlimited.com
 Heb: www.labunlimited.co.uk

 Web: www.labunlimited.com
 Web: www.labunlimited.co.uk



MINISTAR 20 control

/// Data Sheet

The new MINISTAR series by IKA: Developed using the latest cutting-edge technology, this high-tech overhead stirrer with its compact design is ideal for special applications.

Combining high performance with particular excellence, they require minimum space and come with a lifetime guarantee. See for yourself:

"The Fast One" in the high-tech mini class, perfect for low viscosities!

- Hardened glass enclosed, fast response display for maximum visibility and chemical resistance



- State-of-art vibration sensors detect deviations from permissible thresholds and automatically stop the process
- Clear display for all essential information at a glance
- Integrated timer / counter for the control of kinetic sensitive reactions and reminders
- Viscosities up to 10,000 mPas and volumes of up to 15 l
- Continuously adjustable speed between 0/30 2,000 rpm
- USB interface, e.g. for documenting parameters using labworldsoft® or installing firmware updates
- Intuitive and simple handling; touch-sensitive surface for long service life
- Temperature reading on display
- Chemical resistant housing
- Key lock function
- Microprocessor-controlled speed governor for constant rotational speed, even with changes in viscosity



Technical Data

15
60
46
Brushless DC
LCD
50
0/50
2000
10000
42
100
20
Turning knob
1
1
3
chuck
PT1000
yes
0.5 - 8
8.5
yes
extension arm
13
160
yes
electronic
0.2
trend
3
yes
LCD
0 - 6000
-10 - 350
0.1
±0.5 + tolerance PT1000 (DIN IEC 751 Class A)
≤ ± (0.15 + 0.002xITI)
alu-cast coating / thermoplastic polymer
150
70 x 193 x 154
1.56
5 - 40
80
IP 54
yes
100 - 240
50/60



DC Voltage [V=]	24
Current consumption [mA]	2900



Tallaght Business Park Whitestown, Dublin 24, Ireland **D24 RFK3**

Tel: (01) 4523432 Fax: (01) 4523967 Web: www.labunlimited.com Quatro House, Frimley Road, Camberley, United Kingdom **GU16 7ER**

Tel: 08452 30 40 30 Fax: 08452 30 50 30 E-mail: info@labunlimited.com E-mail: info@labunlimited.co.uk Web: www.labunlimited.co.uk