



O 187









Upper basket
for the optimum loading of narrow necked
glassware.

- Equipped with 34 injector nozzles
- Vertical clearance 170 mm¹
- Powder dispensing not possible



Product details - O 187

EAN: 4002511808881 / Article no. 69118701D / Mat. no. 03881640

Product category 	
Washer-disinfectors, laboratory	•
Machine/accessory categorisation 	
PG 8535	•
PG 8536	•
Product category 	
Upper basket	•
Product properties 	
Material	Stainless steel
Colour	Stainless steel
Maximum vertical clearance in mm	170
Area of application 	
Reprocessing laboratory glassware	•
Capacity 	
Narrow-necked glassware [number]	34
Dimensions and weight 	
External dimensions, height in mm	250
External dimensions, width in mm	475
External dimensions, depth in mm	531
Net weight in kg	5,600
Standard accessories 	

Standard accessories

E 351 injector nozzles Ø 4 mm, length 160 mm [number]	34
E 353 clips for injector nozzles E 351 [number]	34

Subject to technical changes; no liability accepted for the accuracy of the information given!

¹ Depending on the laboratory glassware the max. vertical clearance can be lower.

Distributed by:

Lab Unlimited
CARL STUART GROUP

Tallaght Business Park
Whitestown, Dublin 24,
Ireland
D24 RFK3

Tel: (01) 4523432

Fax: (01) 4523967

E-mail: info@labunlimited.com

Web: www.labunlimited.com

Quatro House, Frimley Road,
Camberley,
United Kingdom
GU16 7ER

Tel: 08452 30 40 30

Fax: 08452 30 50 30

E-mail: info@labunlimited.co.uk

Web: www.labunlimited.co.uk