


Agilent Technologies

4500 & 7500a/i/s/c/ce/cs/T-mode Supplies

✓ Standard Configuration • Optional Configuration



CONES

Part #	Description	Cross Ref.	4500	7500a	7500i	7500s	7500c	7500ce	7500cs	T-Mode
HP1001-Ni	Ni Sampler Cone, Cu base	G1820-65238	✓	✓	✓	•	✓	✓	•	
HP1001T-Ni	Ni Sampler, T-Mode Only	G1820-65480	T-Mode	T-Mode	T-Mode					✓
HP1002CE-Ni	Ni Skimmer Cone for 7500ce	G3270-65024						✓		
HP1002C-Ni	Ni Skimmer, 7500c & 7500cs	G1833-65497					✓		•	
HP1002-Ni	Ni Skimmer Cone, 4500	G1820-65265		✓						
HP1002-Ni	Ni Skimmer Cone, 7500	G1820-65050			✓	✓	•			
HP1002T-Ni	Ni Skimmer, T-Mode Only	G1820-65481	T-Mode	T-Mode	T-Mode				✓	
HP1006-Pt	Pt Sampler Cone (10mm Pt insert)	G1820-65239	•	•	•	✓	•	•	✓	
HP1006-Pte	e = SPEC: Pt Sampler Cone* (10mm Pt/Ir insert)	G1820-65239 ¹	•	•	•	✓	•	•	✓	
HP1006A-Pt	Pt Sampler Cone (18mm Pt insert)	G1820-65360	•	•	•	•	•	•	•	•
HP1006A-Pte	e = SPEC: Pt Sampler Cone (18mm Pt/Ir insert)	G1820-65360	•	•	•	•	•	•	•	•
HP1008C-Pt	Pt Skimmer, 7500c	G1833-65092					•			
HP1008CS-Pt	Pt Skimmer Cone, 7500ce & 7500cs	G1833-65132						•	✓	
HP1008CS-Pte	e = SPEC: Pt Skimmer Cone	G1833-65132						•	✓	
HP1008CS-Pt/Ni	Pt skimmer Cone, Nickel Base, 7500ce/7500cs	G3270-60106						•	✓	
HP1008-Pt	Pt Skimmer Cone	G1820-65237	•	•	•	✓				
HP5003e	e = SPEC: Long Life Platinum Shield Plate, 4500	G1820-65357	✓							
HP5004e	e = SPEC: Long Life Platinum Shield Plate, 7500	G1833-65419		✓	✓	✓	✓	✓	✓	✓
HP5001	O-ring for sampler cone	G1820-65025	✓	✓	✓	✓	✓	✓	✓	✓

¹ The 10mm Pt tip of a HP1006-Pte cone is manufactured in a 25mm Ni core for greater durability

Distributed by:



Carl Stuart Limited

ADVANCED APPLIED TECHNOLOGIES

Contact Us:

Irl Ph: 01 4523432

UK Ph: 08452 30 40 30

Web: www.carlstuart.com

Email: info@carlstuart.com