

**U 184/1**

Lower basket  
for the optimum loading of centrifuge tubes,  
phials and test tubes.

- Equipped with 96 injector nozzles
- Only suitable for use with upper baskets O 175 / O 187 / O 184
- Vertical clearance: Approx. 170 mm



## Product details - U 184/1

EAN: 4002513805659 / Article no. 69118403D / Mat. no. 06439180

|   |                            |   |
|---|----------------------------|---|
| Product category  |                            | ^ |
| Washer-disinfectors, laboratory                         |                            | • |
| Machine/accessory categorisation                        |                            | ^ |
| PG 8536   |                            | • |
| Product category  |                            | ^ |
| Lower basket  |                            | • |
| Product properties                                      |                            | ^ |
| Material  | Stainless steel<br>Plastic |   |
| Colour  | Stainless steel<br>Black   |   |
| Maximum vertical clearance in mm                        | 170                        |   |
| Area of application                                     |                            | ^ |
| Reprocessing pipettes/vials                             |                            | • |
| Dimensions and weight                                   |                            | ^ |
| External dimensions, height in mm                       | 250                        |   |
| External dimensions, width in mm                        | 515                        |   |
| External dimensions, depth in mm                        | 531                        |   |
| Net weight in kg  | 5,000                      |   |
| Standard accessories                                    |                            | ^ |
| ID 90, injector nozzles Ø 2.5 mm, length 90 mm [number] | 96                         |   |

Subject to technical changes; no liability accepted for the accuracy of the information given!

Distributed by:



Tallaght Business Park  
Whitestown, Dublin 24,  
Ireland  
D24 RFK3

Tel: (01) 4523432  
Fax: (01) 4523967  
E-mail: [info@labunlimited.com](mailto:info@labunlimited.com)  
Web: [www.labunlimited.com](http://www.labunlimited.com)

Quatro House, Frimley Road,  
Camberley,  
United Kingdom  
GU16 7ER

Tel: 08452 30 40 30  
Fax: 08452 30 50 30  
E-mail: [info@labunlimited.co.uk](mailto:info@labunlimited.co.uk)  
Web: [www.labunlimited.co.uk](http://www.labunlimited.co.uk)