

VE P 2000

Full demineralisation cartridge
For optimum water preparation.

- Max. through flow volume: 450 l/h
- State of depletion measurable with ConductivityMeter CM
- Max water inlet temperature 40 °C - pressure from 0.5-10 bar
- Resin content 12.5 litres



Product details - VE P 2000

EAN: 4002515405307 / Article no. 69431101D / Mat. no. 09773410

Product category	
Fresh-water dishwashers	•
Tank dishwashers	•
Throughfeed tank dishwashers	•
Washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine/accessory categorisation	
8066 g	•
PG 8067	•
PG 8055	•
PG 8055 U	•
PG 8056	•
PG 8056 U	•
PG 8057 TD	•
PG 8057 TD U	•
PG 8058	•
PG 8058 U	•
PG 8059	•
PG 8059 U	•
PG 8060	•
PG 8061	•
PG 8061 U	•
PG 8164	•

PG 8165	•
PG 8166	•
PG 8169	•
PG 8172	•
PG 8172 Eco	•
PG 8581	•
PG 8582	•
PG 8583	•
PG 8583 CD	•
PG 8591	•
PG 8592	•
PG 8593	•

Product category	^
------------------	---

Water pre-treatment	•
---------------------	---

Product properties	^
--------------------	---

Material	Stainless steel Plastic
Colour	Stainless steel Black

Dimensions and weight	^
-----------------------	---

External dimensions, height in mm	410
Diameter in mm	230
Net weight in kg	13,000

Subject to technical changes; no liability accepted for the accuracy of the information given!

Distributed by:

Lab Unlimited
CARL STUART GROUP

Tallaght Business Park
Whitestown, Dublin 24,
Ireland
D24 RFK3

Tel: (01) 4523432
Fax: (01) 4523967

E-mail: info@labunlimited.com
Web: www.labunlimited.com

Quatro House, Frimley Road,
Camberley,
United Kingdom
GU16 7ER

Tel: 08452 30 40 30
Fax: 08452 30 50 30

E-mail: info@labunlimited.co.uk
Web: www.labunlimited.co.uk